

Validation of a multi-gene qPCR-based pharmacogenomics panel across major ethnic groups in Singapore and Indonesia

Running Title

Multi-gene qPCR-based PGx panel for Asians

Authors

Anar Sanjaykumar KOTHARY^{1,*}, Caroline MAHENDRA^{1,*}, Mingchen TAN¹, Eunice Jia Min TAN², Jacelyn PHUA Hong Yi¹, GABRIELLA¹, TAN Xin Hui Jocelyn¹, Jessline Stephanie HARUMAN¹, Chun Kiat LEE², Benedict YAN^{2,3}, Astrid IRWANTO^{1,4,#}

Affiliations

¹Nalagenetics Pte Ltd, Singapore, Singapore

²Department of Laboratory Medicine, National University Health System, Singapore

³Stronghold Diagnostics Lab, Agency for Science, Technology and Research, Singapore

⁴Department of Pharmacy, Faculty of Science, National University of Singapore, Singapore

* equal contribution

corresponding author: astrid@nalagenetics.com

Acknowledgements

We would like to thank Mr. Xie Zhicheng and Mr. Heng Khai Koon at the Genome Institute of Singapore for their assistance with performing the Agena VeriDose® Core and *CYP2D6* CNV Panel studies. We also thank Dr. Tan Zhihao for his feedback on the data analysis approach and manuscript, and PT Nalagenetik Riset Indonesia for their assistance in study recruitment and experimentation performed in Indonesia.

1 **Abstract**

2
3 **Background** – With up to 70% of adverse drug reactions (ADRs) having high genetic associations, the clinical
4 utility of pharmacogenomics (PGx) has been gaining traction. Nala PGx Core™ is a multi-gene qPCR-based panel
5 that comprises 18 variants and 2 *CYP2D6* Copy Number markers across 4 pharmacogenes – *CYP2C9*, *CYP2C19*,
6 *CYP2D6* and *SLCO1B1*.

7
8 **Objectives** – In this study, we validated the performance of Nala PGx Core™ against benchmark methods, on the
9 Singaporean and Indonesian populations. Additionally, we examined the allele and diplotype frequencies across 5
10 major ethnic groups present in these populations namely, Indonesians, Chinese, Malays, Indians and Caucasians.

11
12 **Methods** – Human gDNA samples, extracted from the buccal swabs of 246 participants, were tested on Nala PGx
13 Core™ and two chosen benchmarks, Agena VeriDose® Core and *CYP2D6* Copy Number Variation (CNV) Panel,
14 and TaqMan® DME Genotyping Assays. Performance was evaluated based on assay robustness, precision and
15 accuracy at the genotype- and diplotype-level.

16
17 **Results** – Nala PGx Core™ demonstrated high genotype- and diplotype-level call rates of >97% and >95%
18 respectively in *CYP2D6*, and 100% for *CYP2C9*, *CYP2C19* and *SLCO1B1*. A precision rate of 100% was observed
19 on both intra- and inter-precision studies. Variant-level concordance to the benchmark methods was ≥96.9% across
20 all assays, which consequently resulted in a diplotype-level concordance of ≥94.7% across *CYP2C9*, *CYP2C19* and
21 *CYP2D6*. Overall, the allele frequencies of *CYP2D6**10 and *CYP2D6**36 were higher in our cohort as compared to
22 previous records. Notably, *CYP2D6* copy number variation (CNV) analysis demonstrated a *CYP2D6* *10/*36
23 frequency of 26.5% amongst the Indonesian cohort.

24
25 **Conclusion** –Nala PGx Core™ produced robust and accurate genotyping when compared to other established
26 benchmarks. Furthermore, the panel successfully characterized alleles of clinical relevance in the Singaporean and
27 Indonesian populations such as *CYP2D6**10 and *CYP2D6**36, suggesting its potential for adoption in clinical
28 workflows regionally.

29
30 **Keywords:** Pharmacogenomics; *CYP2C9*; *CYP2C19*; *CYP2D6*; *SLCO1B1*; qPCR; adverse drug reactions;
31 Singapore; Indonesia.

32
33 **Abstract count:** 290 words, 1764 characters

34 **Main text count:** 4337 words

35 Introduction

36 Pharmacogenomics (PGx) is a branch of medicine concerned with determining the response of an individual to
37 therapeutic drugs using genetic and/or genomic information¹. With up to 70% of adverse drug reactions (ADRs)
38 having high genetic associations², and the sunk cost of trial-and-error prescriptions amounting to USD 30 billion³,
39 the clinical utility⁴ of PGx appears to be gaining prominence.

40
41 Applications of PGx include the genotyping of Cytochrome P450 (P450) enzymes which are significantly involved
42 in drug metabolism and chemical toxicity⁵. *CYP2C9*, *CYP2C19* and *CYP2D6* display common polymorphisms that
43 produce altered enzyme activity, often resulting in deviant drug metabolism that leads to side effects from drugs
44 metabolized by them⁵.

45
46 Alterations in *CYP2C9* can disrupt the metabolism of warfarin, while *CYP2C19* is associated with the metabolism of
47 a number of drugs including clopidogrel, omeprazole and phenytoin^{6,7}. The highly polymorphic *CYP2D6* has more
48 than 100 variants and is responsible for the metabolism of approximately 25% of all clinically used drugs, including
49 codeine, fluoxetine and amitriptyline⁷. For this reason, pharmaceuticals are particularly cautious in the development
50 of drugs that are known to interact with these three P450s⁵. Additionally, pharmacogenes, such as *SLCO1B1*, are
51 responsible for encoding transporters of simvastatin and other statins into tissues also affect their response upon
52 intake in individuals. Variants of this gene are associated with increased risk of muscle pain (myalgia) and damage,
53 including myopathy and in severe cases rhabdomyolysis⁸. Previous studies suggest the role of specific *CYP2C19*
54 and *SLCO1B1* variants in the bioactivation of drugs such as clopidogrel and atorvastatin in the Singaporean
55 population^{9,10}. The *CYP2C19* variant, rs4986893, which is implicated in reduced clopidogrel metabolism, has been
56 observed to be 50 times more prevalent in the Singaporean Chinese population than in Europeans¹⁰. The *SLCO1B1*
57 variant, rs4149056, which is implicated in statin metabolism, has been observed to occur less frequently in the
58 Singaporean population⁹. Taken together, these cases highlight the need to harness PGx testing to account for
59 population differences in the personalization of prescriptions¹¹.

60
61 Current methods on PGx testing for the identification of genetic polymorphisms revolve around variant genotyping
62 that utilize technologies ranging from Next Generation Sequencing (NGS), Microarray, Mass Spectrometry to
63 Quantitative PCR (qPCR)^{1,6,12,13}. Main considerations in the adoption of a platform include its robustness,
64 effectiveness and cost-saving, especially in localities with limited resources. qPCR-based assays satisfy these criteria
65 and are a preferred platform for clinical laboratories¹⁴⁻¹⁹. However, some of the available commercial kits are
66 limited in their ability to test the 4 genes of interest mentioned above in a single panel²⁰⁻²⁵. Nala PGx Core™ is a
67 qPCR-based multi-gene panel consisting of 20 pharmacogenomic variants that are associated with altered *CYP2C9*,
68 *CYP2C19*, *CYP2D6* and *SLCO1B1* activity in Asian geographies. Our study aims to evaluate the performance of
69 Nala PGx Core™ in effectively identifying genetic variations that may contribute to PGx-related clinical outcomes
70 in Singapore and Indonesia, with a secondary objective of highlighting the allele and diplotype frequencies in the
71 geographies examined.

72 **Methods and Materials**

73 **Ethics Approval**

74 Institutional Review Board (IRB) approval was granted by the Parkway Independent Ethics Committee (Singapore)
75 under IRB Reference Number PIEC/2019/021 and Siloam Hospitals Ethics Review Committee (Jakarta, Indonesia)
76 under IRB Reference Number 001/EA/KEPKK/RSMRCCC/V/2019.

77 **Study Recruitment**

78 Patients were recruited on behalf of Nalagenetics Pte. Ltd. with written informed consent forms from recruitment
79 sites in Singapore and Indonesia. A total of 251 samples were evaluated from 5 major ethnic groups to ensure the
80 representation of the major ethnic groups residing in both countries – Chinese, Malays, Indians, Caucasians and
81 Indonesians. Participants identifying as one or more of the following ethnicities were categorized as Indonesians:
82 Ambon, Batak, Betawi, Jawa, Lampung, Manado, Minangkabau, Nusa Tenggara Timur, Palembang, Sulawesi,
83 Sunda, Timor Leste, Tolaki and Toraja.

84
85 Buccal samples were collected using OraCollect (Cat No. DNA OCR-100 from DNA Genotek) and genomic DNA
86 (gDNA) extracted using the Monarch Genomic DNA Purification Kit (Cat No. T3010 from NEB). The extraction
87 procedure followed manufacturer's instructions with additional dry-spin step at maximum speed for 1 minute after
88 the 2nd buffer washing step. The quality and concentration of gDNA extracts were quantified by NanoDrop 2000
89 Spectrophotometer (Singapore) and BioDrop- μ LITE (Indonesia). Samples that failed to meet the DNA quality
90 control criteria (n=5) of A260/280 >1.7 and DNA yield >500ng were excluded from the study. The remaining
91 extracted gDNA samples (n=246) were stored at -20°C for downstream application.

92 **Nala PGx Core™**

93 The Nala PGx Core™ kit from Nalagenetics Pte. Ltd. consists of 20 qPCR-based variant assays across four genes –
94 *CYP2C9*, *CYP2C19*, *CYP2D6* and *SLCO1B1*. Whilst assays for *CYP2C9*, *CYP2C19* and *CYP2D6* have been
95 designed to enable the detection of specific star alleles, the *SLCO1B1* assay has been designed to detect the variant
96 rs4149056, which is present in three reduced function haplotypes namely, *SLCO1B1**5, *SLCO1B1**15 and
97 *SLCO1B1**17. The *SLCO1B1* assay is thus, unable to differentiate between each of the three aforementioned
98 haplotypes. The variants covered by the kit are outlined in Table 1.

99
100 Assays were set up on a 96-well plate. Human gDNA was added at a concentration of 2 ng/ μ L as template for the
101 qPCR reaction, which was then performed on the Bio-Rad CFX96 IVD Touch™ Real-Time PCR Detection System
102 per the product insert. Run analysis was performed using the application CFX Manager 3.1 or CFX Maestro, and
103 exported as raw .csv files. Exported files were uploaded into the companion software, Nala Clinical Decision
104 Support™ (Nala CDS™) for further analysis of variant genotyping, diplotype determination and phenotype

105 translation. Genotyping using Nala PGx Core™ was performed at the Molecular Diagnosis Centre, National
106 University Health System, Singapore (NUHS MDC) and PT Nalagenetik Riset Indonesia.

107 **Agena VeriDose® Core and *CYP2D6* Copy Number Variation (CNV) Panel**

108 The VeriDose® Core and *CYP2D6* Copy Number Variation (CNV) Panel from Agena Bioscience® consists of 68
109 variant assays in 20 genes and 5 *CYP2D6* CNV assays, accompanied by a reporting software that automatically
110 analyzes each variation²⁶. Genotyping using Agena VeriDose Core and *CYP2D6* CNV Panel was performed at the
111 Genome Institute of Singapore. Variants evaluated using this platform are listed in Table 1. The Agena VeriDose®
112 Panel has been utilized by the United States Centers for Disease Control and Prevention (CDC) as part of their
113 Genetic Testing Reference Material (GeT-RM) Coordination Program²⁷.

114 **TaqMan® Drug Metabolism Enzyme (DME) Genotyping Assay**

115 TaqMan® DME Genotyping Assays were utilized in the evaluation of *CYP2D6* rs769258 (TaqMan Assay ID
116 AH21B9N) and *CYP2D6* rs267608319 (TaqMan Assay ID C__27102444_F0). Assays were set up on a 384-well
117 plate²⁸ with a sample input of human gDNA at 2 ng/μL. The subsequent PCR reaction was performed on the
118 Applied Biosystems ViiA™ 7 Real-Time PCR System as per the recommended cycling conditions, at the Genome
119 Institute of Singapore. Post-PCR plate read was performed using the companion software, TaqMan® Genotyper™
120 Software for single nucleotide polymorphisms (SNP) genotyping. Similar to the Agena VeriDose® Panel,
121 TaqMan® DME Genotyping Assays were employed in the characterization of DNA samples as part of the CDC
122 GeT-RM program²⁷.

123 124 **Robustness**

125 Genotype- and diplotype-level call rates were defined as the percentage of samples that returned a genotype at the
126 variant-level or were assigned a distinct diplotype for the gene of interest, respectively. Failed tests were defined as
127 samples that did not return a genotype and/or diplotype call for the genes evaluated.

$$\text{Call Rates, \%} = \frac{\text{Total Sample Size} - \text{Failed Tests}}{\text{Total Sample Size}} \times 100\%$$

128

129 Precision

130 Three samples at 3 DNA concentrations were tested across 3 reagent lots on 2 machines. Each test condition was
131 repeated within the same plate for a triplicate. For variant assays that identified SNPs and indels, intra-precision was
132 performed within the same plate, run as triplicates across 47 tests. Inter-precision was assessed from 120 tests
133 performed across plate runs covering the 4 variables – samples, DNA concentration, reagent lots and machines.
134 Concordance rates across precision studies were calculated as the percentage of tests that returned a genotype call
135 concordant to the expected truth for each variant assay. Discordant genotype was defined as instances when the test
136 returned a genotype call that was different from the expected truth.

$$\text{Concordance Rate, \%} = \frac{\text{No. Of Tests Performed} - \text{No. Of Tests With Discordant Genotype}}{\text{No. Of Tests Performed}} \times 100\%$$

137
138 For *CYP2D6* CNV assays, copy number estimates for Intron 2 and Exon 9 of the three samples were derived based
139 on their cycle threshold (Ct) results across plate runs.

$$\text{Copy Number} = 2 \times 2^{-\Delta\Delta\text{Ct}}$$

$$\Delta\Delta\text{Ct} = (\text{Ct}_{\text{reference gene calibrator}} - \text{Ct}_{\text{CYP2D6 calibrator}}) - (\text{Ct}_{\text{reference gene}} - \text{Ct}_{\text{CYP2D6 sample}})$$

140
141
142 Testing of the three samples was repeated for a number of plate runs, n (Table 5), and calculated for the average
143 copy number of each sample and their coefficient of variation (CV). The CV for each plate run was calculated by
144 finding the standard deviation (σ_{plate}) between triplicates within the same plate run, and divided by the triplicate
145 mean (μ_{plate}). The average of the individual CVs was reported as the intra-precision CV. For inter-precision CV,
146 standard deviation population ($\sigma_{\text{plate means}}$) was divided by the mean population, i.e. average of means.

$$\text{Intra- CV, \%} = \frac{\sum_n \frac{\sigma_{\text{plate}}}{\mu_{\text{plate}}}}{n} \times 100\%$$

$$\text{Inter- CV, \%} = \frac{\sigma_{\text{plate means}}}{\sum_n \mu_{\text{plate}} / n} \times 100\%$$

148 Accuracy

149 Variant-level Concordance

150 The accuracy of Nala PGx Core™ in genotyping at a variant-level was evaluated by comparing calls produced by
151 Nala PGx Core™ assay against benchmark methods as listed in Table 1. Samples that successfully produced
152 genotype calls for all variants tested on Nala PGx Core™ and its benchmarks were considered for the evaluation (n
153 = 225). Samples that failed to produce a genotype call on one or more of the platforms were excluded from the
154 concordance calculation (n = 21, Supplementary Table 8). Discordant calls were defined as instances in which Nala
155 PGx Core™ provided a genotype call that was different from that of a call made by the corresponding benchmark.
156 Percentage concordance to the benchmark was calculated per variant as follows –

$$\text{Concordance To Benchmark Per Variant, \%} = \frac{\text{Total Sample Size} - \text{Discordant Calls By Nala PGx Core}^{\text{TM}}}{\text{Total Sample Size}} \times 100\%$$

157 *Diplotype-level Concordance*

158 The accuracy of Nala PGx CoreTM in assigning a diplotype call for *CYP2C9*, *CYP2C19*, and *CYP2D6*, was
159 evaluated by comparing calls against the Agena VeriDose[®] Core and *CYP2D6* CNV Panel. Samples that met the
160 following criteria were included in the sample size of each gene:

- 161 1. Successful genotype-level calls on the relevant platforms for all variants covered by the gene of interest
- 162 2. Successful assignment of a diplotype for the gene of interest on both Nala PGx CoreTM, and Agena
163 VeriDose[®] Core and *CYP2D6* CNV Panel

164 Discordant calls were defined as instances in which Nala PGx CoreTM assigned a diplotype that differed from the
165 call made by the Agena VeriDose[®] Core and *CYP2D6* CNV Panel.

$$\text{Concordance, \%} = \frac{\text{Total Sample Size} - \text{Discordant Calls By Nala PGx Core}^{\text{TM}}}{\text{Total Sample Size}} \times 100\%$$

166 **Frequencies By Ethnicity**

167 Ethnicities were obtained based on participant self-identification across both the population cohorts as part of the
168 recruitment questionnaire. Participants that did not report an ethnic group in the recruitment form were excluded
169 from the frequency analysis. The following were further excluded from the frequency analysis:

- 170 1. Samples in which participants did not report an ethnic group on the recruitment form (n=5)
- 171 2. Samples with one or more variant level failures across the 4 genes evaluated in Table 1 (n=18)
- 172 3. Samples with one or more diplotype-level failures (“No Call”) across *CYP2C9*, *CYP2C19*, and *CYP2D6* on
173 Nala PGx CoreTM and the Agena VeriDose[®] Core and *CYP2D6* CNV Panel (n=14)
- 174 4. Samples with discordant diplotype calls for the gene of interest (n=variable, Supplementary Table 6)

175
176 The remaining samples were included in the allele-level frequency analysis of *CYP2C9* (n=206), *CYP2C19* (n=201),
177 *CYP2D6* (n=195) and *SLCO1B1* (n=203), as well as in the diplotype-level frequency analysis of *CYP2C9* (n=206),
178 *CYP2C19* (n=201) and *CYP2D6* (n=195). Allele and diplotype frequency values were derived using the following
179 formulae^{29,30}, for both the overall study cohort as well as for each ethnic group.

$$\text{Frequency Of Allele 'X' In A Given Population} = \frac{\text{Total Copies Of Allele 'X'}}{\text{Total Copies Of All Alleles For The Gene Of Interest}}$$

180

$$\text{Frequency Of Diplotype 'X' In A Given Population} = \frac{\text{Total Instances Of Diplotype 'X'}}{\text{Total Number Of Individuals In The Population}}$$

181 **Results**

182 **Robustness**

183 We first sought to evaluate the observed genotype- and diplotype-level call rates of the platforms evaluated in this
184 study. 246 samples underwent variant genotyping and diplotype determination, across the four genes evaluated on
185 the genotyping platforms (Tables 2, 3). The genotype-level call rates for Nala PGx Core™ were at 100% for
186 *CYP2C9*, *CYP2C19* and *SLCO1B1*, and the diplotype-level call rates were at 100% for *CYP2C9* and *CYP2C19*. The
187 benchmark platform, Agena VeriDose® Core Panel, demonstrated call rates of >98% at the genotype-level and
188 >97% at the diplotype-level.

189
190 Most variants in *CYP2D6*, except for seven, achieved 100% call rates on Nala PGx Core™, while the corresponding
191 call rates of the benchmark platforms were observed to be between 95.9 – 99.2% on Agena VeriDose® Core and
192 *CYP2D6* CNV Panel, and 100% on TaqMan® DME Genotyping Assays. Out of the seven aforementioned variants,
193 Nala PGx Core™ demonstrated higher call rates than the benchmark for the genotyping of rs1065852, Intron 2 and
194 Exon 9 variants. For rs3892097, rs769258, rs5030865, and rs267608319, the accompanying benchmarks
195 demonstrated higher call rates. At the diplotype-level, Nala PGx Core™ demonstrated a *CYP2D6* call rate of 95.9%
196 as compared to the benchmark, which was observed to be at 90.7% (Table 3).

197 **Precision**

198 A precision study was conducted to assess the consistency of Nala PGx Core™ for samples tested under the same
199 conditions (intra-precision) and under different conditions (inter-precision). Both study resulted in 100%
200 concordance for all assays across replicates, demonstrating consistent genotyping results across a range of DNA
201 concentration, reagent lots and machine variations (Table 4). Precision of *CYP2D6* CNV assay was reported as the
202 average copy number obtained for Intron 2 and Exon 9 of three samples, and their CV calculated across the test
203 conditions (Table 5). The intra- CV ranged from 3-6% while inter- CV between 5-13%, demonstrating high
204 precision of the assays across variables, where acceptable ranges were intra- CV below 10% and inter- CV below
205 15%³¹.

206 **Accuracy**

208 *Variant-level Concordance*

209 To assess the accuracy of the panel, 18 SNPs and 2 *CYP2D6* Copy Number assays were genotyped on the panel,
210 Nala PGx Core™, against benchmark methods as listed in Table 1. The 225 sample cohort consisted of DNA
211 samples isolated from buccal swabs that had successfully produced genotype calls for all variants tested on Nala
212 PGx Core™ and its benchmarks. The sample size per genotype displayed in Table 6 refers to the genotype calls
213 made by the benchmark.

214

215 11 variants (*CYP2C9* rs1799853, rs1057910; *CYP2C19* rs12248560; *CYP2D6* rs5030655, rs3892097, rs35742686,
216 rs28371725, rs769258, rs5030656, rs59421388, rs267608319) were genotyped against Agena VeriDose® Core with
217 a resulting concordance rate of 100% (N = 225 samples). Discordance was observed for *CYP2C19* rs4244285 (n=7),
218 *CYP2C19* rs4986893 (n=3), *CYP2D6* rs1065852 (n=1), *CYP2D6* rs16947 (n=7), *CYP2D6* rs1135840 (n=5),
219 *CYP2D6* rs5030865 (n=1) and *SLCO1B1* rs4149056 (n=6), with a mismatch rate of 0.44% to 3.1% for the affected
220 assays. Overall, Nala PGx Core™ demonstrated >96% concordance to the benchmark, Agena VeriDose® Core, for
221 the 16 variants across 225 samples (Table 6).

222
223 Variants not present on Agena VeriDose® Core, *CYP2D6* rs769258 and *CYP2D6* rs267608319, were genotyped
224 using TaqMan® DME Genotyping Assays. Nala PGx Core™ demonstrated 100% concordance (N = 225) to the
225 benchmark for both SNPs.

226
227 For the *CYP2D6* Intron 2 and Exon 9 Copy Number assays, concordance was observed to be at 99.6% and 98.7%
228 respectively, against the Agena *CYP2D6* CNV Panel. Discordant calls were observed in samples with an Intron 2
229 copy number greater than 3 (n=1), and for samples with an Exon 9 copy number of one (n=1) and two (n=2). The
230 sample size per copy number displayed in Table 7 refers to the copy number calls made by the benchmark.

231
232 A summary of the concordance percentage as well as the distribution of genotypes of each variant are displayed in
233 Table 6. Similarly, the distribution of *CYP2D6* Intron 2 and Exon 9 copy numbers are displayed in Table 7.

234 235 *Diploidy-level Concordance*

236 Following successful genotyping at the variant level, we investigated the accuracy of Nala PGx Core™ in assigning
237 a diploidy call for *CYP2C9*, *CYP2C19* and *CYP2D6*, with reference to the Agena VeriDose® Core and *CYP2D6*
238 CNV Panel. Table 8 displays the percentage concordance after the further exclusion of samples that demonstrated
239 diploidy mismatches arising from technological differences, where technological differences refer to the varying
240 allele coverage of each platform (Supplementary Table 7). These differences were derived from the variant lists of
241 both Nala PGx Core™ (Table 1) and its benchmark, the Agena VeriDose® Core and *CYP2D6* CNV Panel²⁶.

242
243 Overall, a percentage agreement of 100% for *CYP2C9* (n=213), 96.4% for *CYP2C19* (n=223) and 94.7% for
244 *CYP2D6* (n=209) was observed between Nala PGx Core™ and the benchmark. Discordance was observed at n=1
245 for all diploidies listed in Table 8 except for the following with more than one discordant calls: *CYP2C19* *2/*2
246 (n=4), and *CYP2D6* *4/*36 (n=2).

247 248 **Frequencies By Ethnicity**

249 For samples that were concordant on Nala PGx Core™ and the benchmark platforms, we were able to observe the
250 allele frequencies amongst the populations residing in Singapore and Indonesia (Table 9). From the combination of
251 alleles present in individual's chromosome, we were able to observe both the diploidy and corresponding
252 phenotype frequencies amongst our study population (Table 10).

253
254 For *CYP2C9*, *3 allele was the most common amongst Chinese and Malay, and *2 allele amongst Caucasian, which
255 is in line with PharmGKB's reported distribution for the East Asian and European populations respectively. Our
256 study also reported *3 allele as the more common variant in Indian population than *2, as opposed to PharmGKB's
257 frequency. These allele frequencies translated to *1/*3 as a common diplotype observed in Chinese, Malay and
258 Indians, and *1/*2 in Caucasians.

259
260 For *CYP2C19*, we observed highest frequency of *CYP2C19**2 amongst Chinese, Malay and Indonesian which were
261 categorized as East Asian populations. This resulted into high frequency of *1/*2 heterozygous depicted as a
262 common diplotype amongst the population. The alleles *2 and *17 were observed as the common variants at equal
263 proportions of 0.229 in Indians and 0.138 in Caucasians. *CYP2C19**3 was a common minor allele least observed
264 amongst Indians and Caucasians, 0 and 0.017 respectively. As a result, *1/*2 and *1/*17 were common diplotypes
265 observed in Indian and Caucasian populations, and *2/*17 only seen in Indians.

266
267 Common polymorphisms of *CYP2D6* in our population were seen in *10 and *36 alleles, at almost two-fold higher
268 frequencies in Chinese, Malay and Indonesian than Indians. We noticed high frequencies of *36 in our East Asian
269 population. Additionally, two functional copies of *36 allele were present in 1.4% of Chinese population who
270 participated in our study (Figure 1). These alleles resulted in high frequency of *10/*36 as common diplotype
271 amongst Chinese, Malay and Indonesian population. Our Indian population observed equal distribution of alleles *2,
272 *10 and *36, although PharmGKB reported only *2 as the highest frequency in South Asians. The corresponding
273 common diplotypes observed in Indians were of equal proportions in *1/*2, *1/*36 and *1/*10 ranging from 0.12 to
274 0.16. The alleles *2 and *4 were most common amongst Caucasians as reported by PharmGKB in Europeans,
275 resulting in common diplotypes of *1/*2 and *1/*4 in the population.

276
277 For *SLCO1B1*, the frequencies of rs4149056 across all ethnicities were consistent with values reported in gnomAD,
278 with the variant being most common amongst Caucasians (0.167) and least amongst Indians (0.040). The frequency
279 amongst East Asians (0.125, gnomAD), as denoted by the Chinese and Indonesian ethnic groups in this study,
280 ranged between 0.074 and 0.125 respectively.

281

282 Discussion

283 In this study, we evaluated the performance of Nala PGx Core™, a qPCR-based panel that evaluates 18 variants and
284 2 *CYP2D6* Copy Number markers across 4 pharmacogenes with established relevance across major ethnic groups in
285 Singapore and Indonesia population. Nala PGx Core™ comes coupled with a reporting software that supports
286 variant detection, diplotype assignment, diplotype-to-phenotype translation and the generation of reports containing
287 clinical recommendations for each phenotype. The panel demonstrated high genotype-level call rates of >97% for
288 *CYP2D6*, and 100% for *CYP2C9*, *CYP2C19* AND *SLCO1B1*. Similarly, high diplotype-level call rates were
289 observed at >95% for *CYP2D6*, and 100% for *CYP2C9* and *CYP2C19*. A precision of 100% was observed under the

290 same conditions (intra) and across different conditions (inter). In comparison to other established platforms serving
291 as benchmarks during the study, Nala PGx Core™ had $\geq 96.9\%$ concordance rate for all variant level assays, which
292 consequently resulted in $\geq 94.7\%$ concordance at a diplotype level across *CYP2C9*, *CYP2C19* and *CYP2D6*.

293
294 Failures to produce a variant genotype call could be attributed to several reasons. Firstly, failures could potentially
295 stem from the quality of gDNA, despite the DNA quality checks (QC) performed prior to accepting a sample for
296 testing. Further QC that involves specific quantification of double-stranded non-fragmented DNA and traces of other
297 interfering materials like RNA, carryover carbohydrate, residual phenol, guanidine or other reagents could enhance
298 the call rate³². Regardless, the overall higher variant call rates on Nala PGx Core™ panel demonstrate high tolerance
299 of interfering substances, therefore alluding to the high robustness of the assay. Often, failures at variant genotyping
300 subsequently contribute to failures at determining diplotype, since an incomplete variant panel cannot translate into
301 a diplotype. Failures at diplotype calling could also arise from a combination of variants that do not map onto a
302 distinct diplotype, per the reference database, potentially indicating a novel combination.

303
304 We further sought to evaluate the allele frequency distribution in the study cohort across the 5 major ethnic groups
305 observed (Indonesian, Chinese, Malay, Indian and Caucasian). The data presented was limited strictly to the
306 geographical boundaries of Singapore and Indonesia, which could account for the difference in allele frequencies
307 observed in comparison to PharmGKB, which is representative of a more expansive and global cohort³³. Whilst
308 dissimilar to database figures, our study demonstrated the distributions for the following to be concordant with
309 previous studies^{7,34,35}, suggesting a niche in the PGx landscape of Singapore and Indonesia –

- 310 1. *CYP2C9**3 allele as the more common variant within Indians than *CYP2C9**2⁷
- 311 2. A *SLCO1B1* rs4149056 frequency of 12.5% amongst Indonesians³⁴
- 312 3. High frequencies of the *CYP2D6**10 amongst Indonesians³⁴
- 313 4. High frequencies of the *CYP2D6**36 allele as seen in the Indonesian, Chinese and Malay ethnicities³⁵
- 314 5. Two functional copies of *CYP2D6* within our Chinese population³⁵

315 Due to the lack of *CYP2D6* copy number references³⁴, it is our understanding that the frequency of *CYP2D6**36 in
316 Indonesia may not be well-represented. Our study revealed that the prevalence of *CYP2D6**36 to be approximately
317 seventeen times higher amongst Indonesians as compared to the corresponding East Asian allele frequency on
318 PharmGKB³⁶. Furthermore, our study provides insight on the frequency of the *CYP2D6* *10/*36 diplotype in the
319 archipelago, which may help inform the adoption of population-specific PGx workflows in the region. Taken
320 together, the data presents a case for extending tailored PGx testing across the 4 pharmacogenes studied, *CYP2C9*,
321 *CYP2C19*, *CYP2D6* and *SLCO1B1*, in South East Asia.

322 **Limitations**

323 The single largest reason for the diplotype-level discordance observed between Nala PGx Core™ and the Agena
324 VeriDose® Core and *CYP2D6* CNV Panel appears to be a difference in the technological offerings and hence, allele

325 coverage of each panel. From Supplementary Table 7, out of the 28 instances of discordance observed, 9 instances
326 were attributed to this cause (32.1%). This included two alleles each in *CYP2C9*, *CYP2C19* and *CYP2D6*.

327
328 The variants *CYP2C9*12*, *CYP2C9*13*, *CYP2C19*4A*, *CYP2C19*6* and *CYP2D6*13* were covered by Agena
329 VeriDose® Core but not by Nala PGx Core™. Due to lack of coverage for these alleles, the genotypes identified by
330 Nala PGx Core™ were inconsistent with the calls produced by Agena VeriDose® Core. It is important to note that
331 Nala PGx Core™ was curated specifically to include alleles that were relevant to Asian and European populations,
332 with a Minor Allele Frequency (MAF) >1% based on PharmGKB database for East Asians, South Asians and
333 Europeans³⁶. Similarly, *CYP2D6*35*, an allele present in 5.5% of Europeans and 1.2% of South Asians, was
334 identified by Nala PGx Core™ but not by Agena VeriDose® Core. The technological differences and subsequent
335 inconsistencies in reporting suggest the need for further review from a product development standpoint.

336
337 For the highly polymorphic *CYP2D6* gene, Nala PGx Core™ was designed to only cover the variant alleles listed in
338 Table 1, that could lead to overestimation of *1 allele in our population. This means star alleles that may be present
339 in the population were identified as normal genotype, which could result in overestimation of *1 allele in our
340 population as compared to the PharmGKB's reported frequency (Supplementary Table 4).

341
342 The study also relied on self-reported ethnicity, which may have contributed to racial bias in summarizing the
343 observed allele and diplotype frequencies per ethnic group. Future studies may consider screening the ethnicities of
344 participants as well as their next-of-kin, up to two generations prior¹⁰, to improve classification accuracy.

345
346 Finally, the study cohort consisted of 246 individuals, with each ethnic group consisting of approximately 24 to 75
347 individuals for the frequency analysis. Larger scale studies may be conducted to ensure that the trends observed as
348 part of this study may be generalized across the populations examined.

349 **Future Work**

350 Further investigation is required for the remaining samples that had discrepant genotypes across platforms. Nala
351 PGx Core™ was designed as a qPCR platform and genotyping was accompanied by Nala CDS™ that utilized
352 quantitative thresholds. There could be cases of samples displaying borderline or ambiguous qPCR amplification
353 curves and Ct-values, which could be misinterpreted by these thresholds. The benchmark platforms, Agena
354 VeriDose® Core and *CYP2D6* CNV Panel and TaqMan® DME Genotyping Assays, were also analyzed by their
355 accompanying software packages. Such ambiguous cases on either platform could lead to instances of discordance
356 and sequencing of the region of interest may be required to ascertain the truth. For gene variants caused by SNPs
357 and indels, the qPCR amplified product can be submitted for genome sequencing¹³. In the investigation of copy
358 number variants, a two-step approach can be adopted³⁷. Firstly, CNV calls observed in this study may be checked
359 against an alternative benchmark such as digital droplet PCR (ddPCR)³⁸. Secondly, where discordance in CNV
360 calling persists, methods such as XL-PCR can be adopted to identify the presence of *CYP2D6/CYP2D7* hybrid

361 states, with submitting samples for whole-genome sequencing³⁷ being an option. Investigation of copy number
362 variants may require amplification of region spanning across the gene of interest.

363 **Conflict of Interest**

364 A.S.K., M.T., J.P., J.T., C.M., G., J.S.H., and A.I. are employees of Nalagenetics Pte Ltd. A.I. has financial holdings
365 in Nalagenetics Pte Ltd. Nalagenetics Pte Ltd and PT Nalagenetik Riset Indonesia are the license holders of Nala
366 PGx Core™.

367 **Funding**

368 All expenses for this study were funded by Nalagenetics Pte Ltd and PT Nalagenetik Riset Indonesia.

370 **Data Availability Statement**

371 The data that supports the findings of this study are available as supplementary materials. Any additional data is
372 available from the corresponding author upon reasonable request.

373 **References**

- 374 1. Berm EJJ, Loeff M de, Wilffert B, et al. Economic Evaluations of Pharmacogenetic and Pharmacogenomic
375 Screening Tests: A Systematic Review. Second Update of the Literature. Bruns H, ed. *PLoS One*.
376 2016;11(1):e0146262. doi:10.1371/journal.pone.0146262
- 377 2. Chan SL, Ang X, Sani LL, et al. Prevalence and characteristics of adverse drug reactions at admission to
378 hospital: a prospective observational study. *Br J Clin Pharmacol*. 2016;82(6):1636-1646.
379 doi:10.1111/bcp.13081
- 380 3. Sultana J, Cutroneo P, Trifirò G. Clinical and economic burden of adverse drug reactions. *J Pharmacol*
381 *Pharmacother*. 2013;4(SUPPL.1). doi:10.4103/0976-500X.120957
- 382 4. Zhu Y, Lopes GS, Bielinski SJ, et al. Impact of Pharmacogenomic Information on Values of Care and
383 Quality of Life Associated with Codeine and Tramadol-Related Adverse Drug Events. *Mayo Clin Proc*
384 *Innov Qual Outcomes*. 2021;5(1):35-45. doi:10.1016/j.mayocpiqo.2020.08.009
- 385 5. Furge LL, Guengerich FP. Cytochrome P450 enzymes in drug metabolism and chemical toxicology: An
386 introduction. *Biochem Mol Biol Educ*. 2006;34(2):66-74. doi:10.1002/bmb.2006.49403402066
- 387 6. Malsagova KA, Butkova T V, Kopylov AT, et al. Pharmacogenetic Testing: A Tool for Personalized Drug
388 Therapy Optimization. *Pharmaceutics*. 2020;12(12):1-23. doi:10.3390/pharmaceutics12121240
- 389 7. Ling Goh L, Wei Lim C, Cheng Sim W, Xian Toh L, Pang Leong K. Analysis of Genetic Variation in
390 CYP450 Genes for Clinical Implementation. Published online 2017. doi:10.1371/journal.pone.0169233
- 391 8. Kisor DF, Hoefer C, Decker BS. Pharmacogenomics and precision medicine. In: *Clinical Pharmacy*
392 *Education, Practice and Research: Clinical Pharmacy, Drug Information, Pharmacovigilance,*
393 *Pharmacoeconomics and Clinical Research*. Elsevier; 2018:437-451. doi:10.1016/B978-0-12-814276-
394 9.00031-3
- 395 9. Ooi BNS, Raechell, Ying AF, et al. Robust Performance of Potentially Functional SNPs in Machine
396 Learning Models for the Prediction of Atorvastatin-Induced Myalgia. *Front Pharmacol*. 2021;12:605764.
397 doi:10.3389/fphar.2021.605764
- 398 10. Brunham LR, Chan SL, Li R, et al. Pharmacogenomic diversity in Singaporean populations and Europeans.
399 *Pharmacogenomics J*. 2014;14(6):555-563. doi:10.1038/tpj.2014.22
- 400 11. Bachtiar M, Ooi BNS, Wang J, et al. Towards precision medicine: interrogating the human genome to
401 identify drug pathways associated with potentially functional, population-differentiated polymorphisms.
402 *Pharmacogenomics J*. 2019;19(6):516-527. doi:10.1038/s41397-019-0096-y

- 403 12. Perkel J. SNP genotyping: Six technologies that keyed a revolution. *Nat Methods*. 2008;5(5):447-453.
404 doi:10.1038/nmeth0508-447
- 405 13. Yang W, Wu G, Broeckel U, et al. Comparison of genome sequencing and clinical genotyping for
406 pharmacogenes. *Clin Pharmacol Ther*. 2016;100(4):380-388. doi:10.1002/cpt.411
- 407 14. S.A. Deepak, K.R. Kottapalli, R. Rakwal, et al. Real-Time PCR: Revolutionizing Detection and Expression
408 Analysis of Genes. *Curr Genomics*. 2007;8(4):234-251. doi:10.2174/138920207781386960
- 409 15. Pratt VM, Zehnbauser B, Wilson JA, et al. Characterization of 107 genomic DNA reference materials for
410 CYP2D6, CYP2C19, CYP2C9, VKORC1, and UGT1A1: A GeT-RM and association for molecular
411 pathology collaborative project. *J Mol Diagnostics*. 2010;12(6):835-846. doi:10.2353/jmoldx.2010.100090
- 412 16. Pratt VM, Everts RE, Aggarwal P, et al. Characterization of 137 Genomic DNA Reference Materials for 28
413 Pharmacogenetic Genes: A GeT-RM Collaborative Project. *J Mol Diagnostics*. 2016;18(1):109-123.
414 doi:10.1016/j.jmoldx.2015.08.005
- 415 17. Syvänen AC. Accessing genetic variation: Genotyping single nucleotide polymorphisms. *Nat Rev Genet*.
416 2001;2(12):930-942. doi:10.1038/35103535
- 417 18. Suo W, Shi X, Xu S, Li X, Lin Y. Towards low cost, multiplex clinical genotyping: 4-fluorescent
418 Kompetitive Allele-Specific PCR and its application on pharmacogenetics. Kalendar R, ed. *PLoS One*.
419 2020;15(3):e0230445. doi:10.1371/journal.pone.0230445
- 420 19. Shen GQ, Abdullah KG, Wang QK. The TaqMan method for SNP genotyping. *Methods Mol Biol*.
421 2009;578:293-306. doi:10.1007/978-1-60327-411-1_19
- 422 20. Davis BH, DeFrank G, Limdi NA, Harada S. Validation of the Spartan RX CYP2C19 Genotyping Assay
423 Utilizing Blood Samples. *Clin Transl Sci*. 2020;13(2):260-264. doi:10.1111/cts.12714
- 424 21. Zhou Y, Armstead AR, Coshatt GM, Limdi NA, Harada S. Comparison of two point-of-care cyp2c19
425 genotyping assays for genotype-guided antiplatelet therapy. *Ann Clin Lab Sci*. 2017;47(6):738-743.
- 426 22. Bergmeijer TO, Vos GJ, Claassens DM, et al. Feasibility and implementation of CYP2C19 genotyping in
427 patients using antiplatelet therapy. *Pharmacogenomics*. 2018;19(7):621-628. doi:10.2217/pgs-2018-0013
- 428 23. Chae H, Kim M, Koh Y-S, et al. Feasibility of a Microarray-Based Point-of-Care CYP2C19 Genotyping
429 Test for Predicting Clopidogrel On-Treatment Platelet Reactivity. *Biomed Res Int*. 2013;2013:1-5.
430 doi:10.1155/2013/154073
- 431 24. Lee CC, McMillin GA, Babic N, Melis R, Yeo K-TJ. Evaluation of a CYP2C19 genotype panel on the
432 GenMark eSensor® platform and the comparison to the Autogenomics Infiniti™ and Luminex CYP2C19
433 panels. *Clin Chim Acta*. 2011;412(11-12):1133-1137. doi:10.1016/j.cca.2011.03.001
- 434 25. Vanwong N, Ngamsamut N, Hongkaew Y, et al. Detection of CYP2D6 polymorphism using Luminex
435 xTAG technology in autism spectrum disorder: CYP2D6 activity score and its association with risperidone
436 levels. *Drug Metab Pharmacokinet*. 2016;31(2):156-162. doi:10.1016/j.dmpk.2016.01.005
- 437 26. Lois A, Everts RE, Nakorchevsky A, et al. The VeriDose Core Panel: Strong Performance When Analyzing
438 Challenging Pharmacogenetic Samples. Accessed April 2, 2021. [https://agenabio.com/wp-](https://agenabio.com/wp-content/uploads/2020/06/Agena-Conference-Poster-Veridose-Core-V3-G020.pdf)
439 [content/uploads/2020/06/Agena-Conference-Poster-Veridose-Core-V3-G020.pdf](https://agenabio.com/wp-content/uploads/2020/06/Agena-Conference-Poster-Veridose-Core-V3-G020.pdf)
- 440 27. Gaedigk A, Turner A, Everts RE, et al. Characterization of Reference Materials for Genetic Testing of
441 CYP2D6 Alleles. *J Mol Diagnostics*. 2019;21(6):1034-1052. doi:10.1016/j.jmoldx.2019.06.007
- 442 28. Technologies L. *TaqMan® SNP Genotyping Assays*, *TaqMan® Predesigned SNP Genotyping Assays*,
443 *TaqMan® Custom SNP Genotyping Assays*, and *TaqMan® Drug Metabolism Enzyme Genotyping Assays*.;
444 2014.
- 445 29. Population Genetics and the Hardy-Weinberg Principle. Accessed May 2, 2021.
446 <https://www.cs.cmu.edu/~genetics/units/instructions/instructions-PGE.pdf>
- 447 30. Allele frequency. Accessed May 2, 2021. <https://www.nature.com/scitable/definition/allele-frequency-298/>
- 448 31. Calculating Inter- and Intra-Assay Coefficients of Variability. Accessed May 1, 2021.
449 <https://salimetrics.com/calculating-inter-and-intra-assay-coefficients-of-variability/>
- 450 32. Matlock B. Assessment of Nucleic Acid Purity. *Tech Bull NanoDrop Spectrophotometers*. Published online
451 2015:1-2.
- 452 33. Huddart R, Fohner AE, Whirl Carrillo M, et al. Standardized Biogeographic Grouping System for
453 Annotating Populations in Pharmacogenetic Research. *Clin Pharmacol Ther*. 2019;105(5):1256-1262.
454 doi:10.1002/cpt.1322
- 455 34. Runcharoen C, Fukunaga K, Sensorn I, et al. Prevalence of pharmacogenomic variants in 100
456 pharmacogenes among Southeast Asian populations under the collaboration of the Southeast Asian
457 Pharmacogenomics Research Network (SEAPharm). *Hum Genome Var*. 2021;8(1):7. doi:10.1038/s41439-
458 021-00135-z

- 459 35. Chan W, Li MS, Sundaram SK, Tomlinson B, Cheung PY, Tzang CH. CYP2D6 allele frequencies, copy
460 number variants, and tandems in the population of Hong Kong. *J Clin Lab Anal.* 2019;33(1).
461 doi:10.1002/jcla.22634
- 462 36. Whirl-Carrillo M, McDonagh EM, Hebert JM, et al. Pharmacogenomics Knowledge for Personalized
463 Medicine. *Clin Pharmacol Ther.* 2012;92(4):414-417. doi:10.1038/clpt.2012.96
- 464 37. Pös O, Radvanszky J, Styk J, et al. Copy Number Variation: Methods and Clinical Applications. *Appl Sci.*
465 2021;11(2):819. doi:10.3390/app11020819
- 466 38. Gross AM, Ajay SS, Rajan V, et al. Copy-number variants in clinical genome sequencing: deployment and
467 interpretation for rare and undiagnosed disease. *Genet Med.* 2019;21(5):1121-1130. doi:10.1038/s41436-
468 018-0295-y
469
- 470

471 **Tables**

472 **Table 1** Genes and variants evaluated

Gene	Star Allele	Variant	Nucleotide Changes	Effect On Protein	Allele Clinical Function Status	Genotyping Methods Utilized		
						Nala PGx Core™	Agena VeriDose® Core and CYP2D6 CNV Panel	TaqMan® DME Genotyping Assays
<i>CYP2C9</i>	*2	rs1799853	3608C>T	R144C	Decreased	Bi-allelic Assay	VeriDose Core	
	*3	rs1057910	42614A>C	I359L	None	Bi-allelic Assay	VeriDose Core	
<i>CYP2C19</i>	*2	rs4244285	19154G>A	Splicing Defect	None	Bi-allelic Assay	VeriDose Core	
	*3	rs4986893	17948G>A	W212X	None	Bi-allelic Assay	VeriDose Core	
	*17	rs12248560	-806C>T	5' Region	Increased	Bi-allelic Assay	VeriDose Core	
<i>CYP2D6</i>	*2	rs1135840	4181G>C	S486T	Normal	Bi-allelic Assay	VeriDose Core	
	*2	rs16947	2851C>T	R296C	Normal	Bi-allelic Assay	VeriDose Core	
	*3	rs35742686	2550delA	Frameshift	None	Bi-allelic Assay	VeriDose Core	
	*4	rs3892097, rs1065852	1847G>A, 100C>T	Splicing Defect, P34S	None	Bi-allelic Assay	VeriDose Core	
	*5	N/A	N/A	Gene Deletion	None	CNV Assay (Intron 2)	CYP2D6 CNV	
	*6	rs5030655	1708delT	Frameshift	None	Bi-allelic Assay	VeriDose Core	
	*8	rs5030865	1759G>T	G169X	None	Tri-allelic Assay	VeriDose Core	
	*9	rs5030656	2616delAAG	K281del	Decreased	Bi-allelic Assay	VeriDose Core	
	*10	rs3892097, rs1065852	1847G>A, 100C>T	Splicing Defect, P34S	Decreased	Bi-allelic Assay	VeriDose Core	
	*14	rs5030865	1759G>A	G169R	Decreased	Tri-allelic Assay	VeriDose Core	
	*21	rs72549352	2580_2581ins C	Frameshift	None	Bi-allelic Assay	VeriDose Core	
	*29	rs59421388	3184G>A	V338M	Decreased	Bi-allelic Assay	VeriDose Core	
	*31	rs267608319	4043G>A	R440H	None	Bi-allelic Assay	NA	AH21B9N
	*35	rs769258	31G>A	V11M	Normal	Bi-allelic Assay	NA	C__27102444_F0
	*36	N/A	Recombination at Exon 9	CYP2D6-2D7 Hybrid	None	CNV Assay (Exon 9)	CYP2D6 CNV	
*41	rs28371725	2989G>A	Splicing Defect	Decreased	Bi-allelic Assay	VeriDose Core		
<i>SLCO1B1</i>	NA [†]	rs4149056	g.52422T>C	V174A	Decreased	Bi-allelic Assay	VeriDose Core	

[†]Nala PGx Core™ detects the variant, rs4149056, which is associated with decreased enzymatic activity and is present in three known *SLCO1B1* haplotypes namely, *SLCO1B1**5, *SLCO1B1**15 and *SLCO1B1**17.

473
474
475

476 **Table 2** Observed genotype-level call rates per variant per gene per platform

Gene	Variant	Variant Call Rate, % (n=246)		
		Nala PGx Core™	Agena VeriDose® Core and CYP2D6 CNV Panel	TaqMan® DME Genotyping Assays
<i>CYP2C9</i>	rs1799853	100.0	99.2	NA
	rs1057910	100.0	99.6	NA
<i>CYP2C19</i>	rs4244285	100.0	99.6	NA
	rs4986893	100.0	99.6	NA
	rs12248560	100.0	98.8	NA
<i>CYP2D6</i>	rs1065852	98.4	95.9	NA
	rs5030655	100.0	99.2	NA
	rs3892097	98.8	99.2	NA
	rs35742686	100.0	98.8	NA
	rs16947	100.0	99.6	NA
	rs28371725	100.0	99.2	NA
	rs1135840	100.0	99.6	NA
	rs769258	98.8	NA	100.0
	rs5030865	97.2	99.2	NA
	rs5030656	100.0	99.6	NA
	rs59421388	100.0	99.2	NA
	rs267608319	99.6	NA	100.0
		CNV Assay (Intron 2)	99.6	99.2
	CNV Assay (Exon 9)	99.6	99.2	NA
<i>SLCO1B1</i>	rs4149056	100.0	99.2	NA

477
478 **Table 3** Observed diplotype-level call rates per gene per platform

Gene	Diplotype Call Rate, % (n=246)	
	Nala PGx Core™	Agena VeriDose® Core and CYP2D6 CNV Panel
<i>CYP2C9</i>	100.0	97.2
<i>CYP2C19</i>	100.0	98.8
<i>CYP2D6</i>	95.9	90.7

479
480 **Table 4** Intra-precision and inter-precision concordance rates

Gene	Variant	Concordance rate, %	
		Intra-precision (47 tests in triplicate)	Inter-precision (120 tests across 4 variables)
<i>CYP2C9</i>	rs1799853	100%	100%
	rs1057910	100%	100%
<i>CYP2C19</i>	rs4244285	100%	100%
	rs4986893	100%	100%
	rs12248560	100%	100%
<i>CYP2D6</i>	rs1065852	100%	100%
	rs5030655	100%	100%
	rs3892097	100%	100%
	rs35742686	100%	100%

	rs16947	100%	100%
	rs28371725	100%	100%
	rs1135840	100%	100%
	rs769258	100%	100%
	rs5030865	100%	100%
	rs5030656	100%	100%
	rs59421388	100%	100%
	rs267608319	100%	100%
<i>SLCO1B1</i>	rs4149056	100%	100%

481
482

Table 5 Intra-precision and inter-precision for *CYP2D6* Copy Number

Sample	Variant	Average Copy Number	Sample Size, n	Coefficient of Variation	
				Intra-precision	Inter-precision
221374	Intron 2	2	11	4%	6%
	Exon 9	1		5%	5%
222660	Intron 2	2	9	3%	5%
	Exon 9	2		4%	6%
222769	Intron 2	3	22	4%	10%
	Exon 9	2		6%	13%

483
484

Table 6 Genotype concordance for genes *CYP2C9*, *CYP2C19*, *CYP2D6* and *SLCO1B1*

Gene	Variant	Genotype	Sample Size Per Genotype [†]	Concordance to Benchmark, % [‡]	Discordant Calls By Nala PGx Core™ [§]
<i>CYP2C9</i>	rs1799853 N = 225	C/C	212	100.0	0
		C/T	13		0
		T/T	0		0
	rs1057910 N = 225	A/A	204	100.0	0
		A/C	21		0
		C/C	0		0
<i>CYP2C19</i>	rs4244285 N = 225	G/G	121	96.9	1
		A/G	89		2
		A/A	15		4
	rs4986893 N = 225	G/G	206	98.7	2
		A/G	19		1
		A/A	0		0
	rs12248560 N = 225	C/C	0	100.0	0
		C/T	0		0
		T/T	0		0
<i>CYP2D6</i>	rs1065852 N = 225	G/G	77	99.6	0
		A/G	76		1
		A/A	72		0
	rs5030655 N = 225	A/A	224	100.0	0
		delA/A	1		0
		delA/delA	0		0
	rs3892097 N = 225	C/C	205	100.0	0
		C/T	16		0
		T/T	4		0
	rs35742686 N = 225	T/T	225	100.0	0
		delT/T	0		0
		delT/delT	0		0
	rs16947 N = 225	G/G	145	96.9	2
		A/G	64		4
		A/A	16		1
	rs28371725 N = 225	C/C	204	100.0	0
		C/T	19		0
		T/T	2		0
	rs1135840	C/C	33	97.8	1

	N = 225	C/G	68	100.0	1
		G/G	124		3
	rs769258 N = 225	C/C	220	99.6	0
		C/T	5		0
		T/T	0		0
	rs5030865 N = 225	C/C	222	100.0	0
		C/T	3		1
		T/T	0		0
	rs5030656 N = 225	CT/CTT	224	100.0	0
		delCTT/CTT	1		0
		delCTT/delCTT	0		0
	rs59421388 N = 225	C/C	224	100.0	0
		C/T	1		0
		T/T	0		0
	rs267608319 N = 225	C/C	225	100.0	0
C/T		0	0		
T/T		0	0		
SLCO1B1 N = 225	T/T	182	97.3	2	
	C/T	42		4	
	C/C	1		0	

[†]“Sample Size Per Genotype” reflects the counts of each genotype as produced by the benchmark, and not by Nala PGx Core™.

[‡]Refer to Table 1 for an outline of the benchmark methods used per variant.

[§]Discordant calls refer to instances in which Nala PGx Core™ did not produce a genotype call concordant with the benchmark.

485

486

Table 7 Genotype concordance for *CYP2D6* Intron 2 And Exon 9 Copy Number Variations

Gene	Variant	Copy Number	Sample Size Per Copy Number [†]	Concordance to Benchmark, % [‡]	Discordant Calls By Nala PGx Core™ [§]
<i>CYP2D6</i>	Intron 2 N = 224	0	0	99.6	0
		1	12		0
		2	130		0
		>=3	82		1
	Exon 9 N = 224	0	0	98.7	0
		1	24		1
		2	193		2
		>=3	7	0	

[†]“Sample Size Per Copy Number” reflects the counts of each copy number call as produced by the benchmark, and not by Nala PGx Core™.

[‡]Refer to Table 1 for an outline of the benchmark methods used per variant.

[§]Discordant calls refer to instances in which Nala PGx Core™ did not produce a copy number call concordant with the benchmark.

487

488

489

Table 8 Diplotype concordance for *CYP2C9*, *CYP2C19* and *CYP2D6* between Nala PGx Core™, and Agena VeriDose® Core and *CYP2D6* CNV Panel

Genes	Concordance, % [†]	Discordant Diplotypes		
		Agena VeriDose® Core and <i>CYP2D6</i> CNV Panel	Nala PGx Core™	Instances
<i>CYP2C9</i>	100% (n = 213)	NA	NA	NA
<i>CYP2C19</i>	96.4% (n = 223)	*1/*1	*1/*3	1
		*1/*2	*1/*1	1
		*1/*2	*1/*3	1
		*1/*3	*2/*2	1
		*2/*2	*1/*1	2
		*2/*2	*1/*2	2
<i>CYP2D6</i>	94.7% (n = 209)	*1/*1	*2/*2	1
		*1/*10	*1/*10, CN >=3	1
		*1/*41	*39/*41	1
		*2/*2	*1/*1	1

		*2/*10	*1/*10	1
		*2/*36, CN >=3	*2/*10, CN >=3	1
		*4/*36	*4/*10	2
		*10/*36	*2/*36	1
		*13 *1/*10	*1/*10, CN >=3	1
		*13 *1/*41, CN = 1	*1/*41	1

[†]The concordance presented in this table excludes samples that have mismatches in diplotype calls arising from technological differences between platforms. Please refer to Supplementary Table 7 for a detailed breakdown of all diplotypes detected.

490
491

Table 9 Observed allele frequencies by ethnicity

Gene	Allele or Variant	Allele Frequencies (Per This Study)						Allele Frequencies (PharmGKB)		
		Indonesian	Chinese	Malay	Indian	Caucasian	Overall	East Asian	Central/South Asian	European
<i>CYP2C9</i> (n = 206)	*2	0.000	0.000	0.000	0.060	0.172	0.032	0.002	0.114	0.127
	*3	0.000	0.040	0.039	0.100	0.069	0.044	0.038	0.110	0.076
<i>CYP2C19</i> (n = 201)	*2	0.297	0.274	0.289	0.229	0.138	0.256	0.284	0.270	0.147
	*3	0.041	0.055	0.066	0.000	0.017	0.042	0.072	0.016	0.002
	*17	0.054	0.007	0.039	0.229	0.138	0.067	0.021	0.171	0.216
<i>CYP2D6</i> (n = 195)	*2	0.088	0.118	0.167	0.140	0.196	0.136	0.121	0.295	0.277
	*4	0.000	0.035	0.069	0.020	0.161	0.051	0.005	0.091	0.185
	*5	0.029	0.028	0.028	0.080	0.000	0.031	0.049	0.046	0.030
	*6	0.000	0.000	0.000	0.000	0.018	0.003	0.000	0.000	0.011
	*9	0.000	0.000	0.000	0.000	0.018	0.003	0.002	0.003	0.028
	*10	0.338	0.299	0.222	0.160	0.143	0.251	0.436	0.087	0.016
	*14	0.000	0.014	0.000	0.000	0.000	0.005	0.003	ND	0.000
	*29	0.000	0.000	0.000	0.000	0.018	0.003	0.000	0.003	0.001
*36	0.206	0.250	0.181	0.160	0.054	0.190	0.012	0.000	0.000	
*41	0.044	0.056	0.028	0.040	0.071	0.049	0.023	0.123	0.092	
<i>SLCO1B1</i> (n = 203)	rs4149056 [†]	0.125	0.074	0.026	0.040	0.167	0.084	0.125 [‡]	0.050 [‡]	0.159 [‡]

"ND" refers to instances in which no data is available for the given allele on PharmGKB.

[†]rs4149056 refers to the reduced function variant of *SLCO1B1* that is present in *SLCO1B1**5, *SLCO1B1**15 and *SLCO1B1**17.

[‡]Allele frequency values for rs4149056 have been obtained from gnomAD.

492
493

Table 10 Observed diplotype frequencies by ethnicity

Gene	Diplotype	Phenotype [‡]	Diplotype Frequencies (Per This Study) [†]											
			Indonesian		Chinese		Malay		Indian		Caucasian		Overall	
			Obs	Freq	Obs	Freq	Obs	Freq	Obs	Freq	Obs	Freq	Obs	Freq
<i>CYP2C9</i>	*1/*1	NM	39	1.000	69	0.920	35	0.921	17	0.680	16	0.552	176	0.854
	*1/*2	IM	0	0.000	0	0.000	0	0.000	3	0.120	9	0.310	12	0.058
	*1/*3	IM	0	0.000	6	0.080	3	0.079	5	0.200	3	0.103	17	0.083
	*2/*3	PM	0	0.000	0	0.000	0	0.000	0	0.000	1	0.034	1	0.005
<i>CYP2C19</i>	*1/*1	NM	10	0.270	31	0.425	14	0.368	8	0.333	14	0.483	77	0.383
	*1/*2	IM	19	0.514	28	0.384	14	0.368	5	0.208	7	0.241	73	0.363
	*1/*3	IM	3	0.081	6	0.082	3	0.079	0	0.000	1	0.034	13	0.065
	*1/*17	RM	3	0.081	1	0.014	1	0.026	5	0.208	5	0.172	15	0.075
	*2/*2	PM	1	0.027	5	0.068	3	0.079	1	0.042	0	0.000	10	0.050
	*2/*3	PM	0	0.000	2	0.027	2	0.053	0	0.000	0	0.000	4	0.020
	*2/*17	IM	1	0.027	0	0.000	0	0.000	4	0.167	1	0.034	6	0.030
	*17/*17	UM	0	0.000	0	0.000	1	0.026	1	0.042	1	0.034	3	0.015
<i>CYP2D6</i>	*1/*1	NM	6	0.176	7	0.097	5	0.139	4	0.160	2	0.071	24	0.123
	*1/*2	NM	0	0.000	2	0.028	2	0.056	3	0.120	3	0.107	10	0.051
	*1/*2, CN >=3	UM	0	0.000	0	0.000	0	0.000	1	0.040	1	0.036	2	0.010
	*1/*4	IM	0	0.000	2	0.028	1	0.028	0	0.000	4	0.143	7	0.036
	*1/*5	IM	0	0.000	2	0.028	1	0.028	1	0.040	0	0.000	4	0.021


*1/*6	IM	0	0.000	0	0.000	0	0.000	0	0.000	1	0.036	1	0.005
*1/*10	NM	4	0.118	3	0.042	2	0.056	2	0.080	2	0.071	13	0.067
*1/*14	NM	0	0.000	1	0.014	0	0.000	0	0.000	0	0.000	1	0.005
*1/*29	NM	0	0.000	0	0.000	0	0.000	0	0.000	1	0.036	1	0.005
*1/*36	IM	3	0.088	3	0.042	5	0.139	4	0.160	0	0.000	15	0.077
*1/*41	NM	1	0.029	2	0.028	1	0.028	0	0.000	2	0.071	6	0.031
*1/*41, CN >=3	UM	0	0.000	0	0.000	0	0.000	1	0.040	0	0.000	1	0.005
*2/*2	NM	0	0.000	1	0.014	1	0.028	1	0.040	1	0.036	4	0.021
*2/*4	IM	0	0.000	0	0.000	1	0.028	0	0.000	2	0.071	3	0.015
*2/*5	IM	0	0.000	0	0.000	1	0.028	1	0.040	0	0.000	2	0.010
*2/*10	NM	4	0.118	4	0.056	4	0.111	0	0.000	1	0.036	13	0.067
*2/*14	NM	0	0.000	1	0.014	0	0.000	0	0.000	0	0.000	1	0.005
*2/*36	IM	2	0.059	4	0.056	2	0.056	0	0.000	1	0.036	9	0.046
*2/*41	NM	0	0.000	4	0.056	0	0.000	0	0.000	1	0.036	5	0.026
*4/*4	PM	0	0.000	1	0.014	1	0.028	0	0.000	1	0.036	3	0.015
*4/*5	PM	0	0.000	0	0.000	0	0.000	1	0.040	0	0.000	1	0.005
*4/*9	IM	0	0.000	0	0.000	0	0.000	0	0.000	1	0.036	1	0.005
*4/*10	IM	0	0.000	1	0.014	1	0.028	0	0.000	0	0.000	2	0.010
*5/*10	IM	0	0.000	2	0.028	0	0.000	1	0.040	0	0.000	3	0.015
*5/*41	IM	2	0.059	0	0.000	0	0.000	0	0.000	0	0.000	2	0.010
*10/*10	IM	3	0.088	1	0.014	2	0.056	0	0.000	2	0.071	8	0.041
*10/*10, CN >=3	IM	0	0.000	1	0.014	0	0.000	0	0.000	0	0.000	1	0.005
*10/*36	IM	9	0.265	28	0.389	5	0.139	4	0.160	1	0.036	47	0.241
*10/*41	IM	0	0.000	1	0.014	0	0.000	1	0.040	0	0.000	2	0.010
*36/*41	IM	0	0.000	1	0.014	1	0.028	0	0.000	1	0.036	3	0.015

[†]"Obs" and "Freq" are abbreviations for "Observations" and "Frequency" respectively.

[‡]"NM," "IM," "PM" and "UM" are abbreviations for "Normal Metabolizer," "Intermediate Metabolizer," "Poor Metabolizer" and "Ultrarapid Metabolizer" respectively.


494 Figure Legends

495 **Figure 1. CYP2D6*36 Frequency By Ethnicity.** Distribution of individuals carrying one or two copies of the
496 CYP2D6*36 allele among the study cohort (n=205), grouped per ethnicity.



497
498

499 **Supplementary Figure 1. Observed *CYP2D6* Functional Copy Number Frequencies By Ethnicity.** Distribution
500 of individuals carrying one, two or three or more copies of the *CYP2D6* functional gene among the study cohort
501 (n=205), grouped per ethnicity.



502

Validation of a multi-gene qPCR-based pharmacogenomics panel across major ethnic groups in Singapore and Indonesia

A. Main Tables and Figures

Table 1 Genes and variants evaluated

Gene	Star Allele	Variant	Nucleotide Changes	Effect On Protein	Allele Clinical Function Status	Genotyping Methods Utilized		
						Nala PGx Core™	Agena VeriDose® Core and CYP2D6 CNV Panel	TaqMan® DME Genotyping Assays
<i>CYP2C9</i>	*2	rs1799853	3608C>T	R144C	Decreased	Bi-allelic Assay	VeriDose Core	
	*3	rs1057910	42614A>C	I359L	None	Bi-allelic Assay	VeriDose Core	
<i>CYP2C19</i>	*2	rs4244285	19154G>A	Splicing Defect	None	Bi-allelic Assay	VeriDose Core	
	*3	rs4986893	17948G>A	W212X	None	Bi-allelic Assay	VeriDose Core	
	*17	rs12248560	-806C>T	5' Region	Increased	Bi-allelic Assay	VeriDose Core	
<i>CYP2D6</i>	*2	rs1135840	4181G>C	S486T	Normal	Bi-allelic Assay	VeriDose Core	
	*2	rs16947	2851C>T	R296C	Normal	Bi-allelic Assay	VeriDose Core	
	*3	rs35742686	2550delA	Frameshift	None	Bi-allelic Assay	VeriDose Core	
	*4	rs3892097, rs1065852	1847G>A, 100C>T	Splicing Defect, P34S	None	Bi-allelic Assay	VeriDose Core	
	*5	N/A	N/A	Gene Deletion	None	CNV Assay (Intron 2)	CYP2D6 CNV	
	*6	rs5030655	1708delT	Frameshift	None	Bi-allelic Assay	VeriDose Core	
	*8	rs5030865	1759G>T	G169X	None	Tri-allelic Assay	VeriDose Core	
	*9	rs5030656	2616delAAG	K281del	Decreased	Bi-allelic Assay	VeriDose Core	
	*10	rs3892097, rs1065852	1847G>A, 100C>T	Splicing Defect, P34S	Decreased	Bi-allelic Assay	VeriDose Core	
	*14	rs5030865	1759G>A	G169R	Decreased	Tri-allelic Assay	VeriDose Core	
	*21	rs72549352	2580_2581insC	Frameshift	None	Bi-allelic Assay	VeriDose Core	
	*29	rs59421388	3184G>A	V338M	Decreased	Bi-allelic Assay	VeriDose Core	
	*31	rs267608319	4043G>A	R440H	None	Bi-allelic Assay	NA	AH21B9N
	*35	rs769258	31G>A	V11M	Normal	Bi-allelic Assay	NA	C__27102444_F0
	*36	N/A	Recombination at Exon 9	CYP2D6-2D7 Hybrid	None	CNV Assay (Exon 9)	CYP2D6 CNV	
*41	rs28371725	2989G>A	Splicing Defect	Decreased	Bi-allelic Assay	VeriDose Core		
<i>SLCO1B1</i>	NA [†]	rs4149056	g.52422T>C	V174A	Decreased	Bi-allelic Assay	VeriDose Core	

[†]Nala PGx Core™ detects the variant, rs4149056, which is associated with decreased enzymatic activity and is present in three known *SLCO1B1* haplotypes namely, *SLCO1B1**5, *SLCO1B1**15 and *SLCO1B1**17.

Table 2 Observed genotype-level call rates per variant per gene per platform

Gene	Variant	Variant Call Rate, % (n=246)		
		Nala PGx Core™	Agena VeriDose® Core and CYP2D6 CNV Panel	TaqMan® DME Genotyping Assays
<i>CYP2C9</i>	rs1799853	100.0	99.2	NA
	rs1057910	100.0	99.6	NA
<i>CYP2C19</i>	rs4244285	100.0	99.6	NA
	rs4986893	100.0	99.6	NA
<i>CYP2D6</i>	rs12248560	100.0	98.8	NA
	rs1065852	98.4	95.9	NA
	rs5030655	100.0	99.2	NA
	rs3892097	98.8	99.2	NA
	rs35742686	100.0	98.8	NA
	rs16947	100.0	99.6	NA
	rs28371725	100.0	99.2	NA
	rs1135840	100.0	99.6	NA
	rs769258	98.8	NA	100.0
	rs5030865	97.2	99.2	NA
	rs5030656	100.0	99.6	NA
	rs59421388	100.0	99.2	NA
	rs267608319	99.6	NA	100.0
	CNV Assay (Intron 2)	99.6	99.2	NA
	CNV Assay (Exon 9)	99.6	99.2	NA
<i>SLCO1B1</i>	rs4149056	100.0	99.2	NA

Table 3 Observed diplotype-level call rates per gene per platform

Gene	Diplotype Call Rate, % (n=246)	
	Nala PGx Core™	Agena VeriDose® Core and CYP2D6 CNV Panel
<i>CYP2C9</i>	100.0	97.2
<i>CYP2C19</i>	100.0	98.8
<i>CYP2D6</i>	95.9	90.7

Table 4 Intra-precision and inter-precision concordance rates

Gene	Variant	Concordance rate, %	
		Intra-precision (47 tests in triplicate)	Inter-precision (120 tests across 4 variables)
<i>CYP2C9</i>	rs1799853	100%	100%
	rs1057910	100%	100%
<i>CYP2C19</i>	rs4244285	100%	100%
	rs4986893	100%	100%
	rs12248560	100%	100%
<i>CYP2D6</i>	rs1065852	100%	100%
	rs5030655	100%	100%
	rs3892097	100%	100%
	rs35742686	100%	100%
	rs16947	100%	100%
	rs28371725	100%	100%
	rs1135840	100%	100%
	rs769258	100%	100%
	rs5030865	100%	100%
	rs5030656	100%	100%
	rs59421388	100%	100%
	rs267608319	100%	100%
<i>SLCO1B1</i>	rs4149056	100%	100%

Table 5 Intra-precision and inter-precision for *CYP2D6* Copy Number

Sample	Variant	Average Copy Number	Sample Size, n	Coefficient of Variation	
				Intra-precision	Inter-precision
221374	Intron 2	2	11	4%	6%
	Exon 9	1		5%	5%
222660	Intron 2	2	9	3%	5%
	Exon 9	2		4%	6%
222769	Intron 2	3	22	4%	10%
	Exon 9	2		6%	13%

Table 6 Genotype concordance for genes *CYP2C9*, *CYP2C19*, *CYP2D6* and *SLCO1B1*

Gene	Variant	Genotype	Sample Size Per Genotype [†]	Concordance to Benchmark, % [‡]	Discordant Calls By Nala PGx Core™ [§]
<i>CYP2C9</i>	rs1799853 N = 225	C/C	212	100.0	0
		C/T	13		0
		T/T	0		0
	rs1057910 N = 225	A/A	204	100.0	0
		A/C	21		0
		C/C	0		0
<i>CYP2C19</i>	rs4244285 N = 225	G/G	121	96.9	1
		A/G	89		2
		A/A	15		4
	rs4986893 N = 225	G/G	206	98.7	2
		A/G	19		1
		A/A	0		0
	rs12248560 N = 225	C/C	0	100.0	0
		C/T	0		0
		T/T	0		0
<i>CYP2D6</i>	rs1065852 N = 225	G/G	77	99.6	0
		A/G	76		1
		A/A	72		0
	rs5030655 N = 225	A/A	224	100.0	0
		delA/A	1		0
		delA/delA	0		0
	rs3892097 N = 225	C/C	205	100.0	0
		C/T	16		0
		T/T	4		0
	rs35742686 N = 225	T/T	225	100.0	0
		delT/T	0		0
		delT/delT	0		0
	rs16947 N = 225	G/G	145	96.9	2
		A/G	64		4
		A/A	16		1
	rs28371725 N = 225	C/C	204	100.0	0
		C/T	19		0
		T/T	2		0
	rs1135840 N = 225	C/C	33	97.8	1
		C/G	68		1
		G/G	124		3
	rs769258 N = 225	C/C	220	100.0	0
		C/T	5		0
		T/T	0		0
	rs5030865 N = 225	C/C	222	99.6	0
		C/T	3		1
		T/T	0		0
	rs5030656 N = 225	CTT/CTT	224	100.0	0
		delCTT/CTT	1		0
		delCTT/delCTT	0		0
	rs59421388 N = 225	C/C	224	100.0	0
		C/T	1		0
		T/T	0		0
rs267608319 N = 225	C/C	225	100.0	0	
	C/T	0		0	
	T/T	0		0	
<i>SLCO1B1</i>	rs4149056 N = 225	T/T	182	97.3	2
		C/T	42		4
		C/C	1		0

[†]“Sample Size Per Genotype” reflects the counts of each genotype as produced by the benchmark, and not by Nala PGx Core™.

[‡]Refer to Table 1 for an outline of the benchmark methods used per variant.

[§]Discordant calls refer to instances in which Nala PGx Core™ did not produce a genotype call concordant with the benchmark.

Table 7 Genotype concordance for *CYP2D6* Intron 2 And Exon 9 Copy Number Variations

Gene	Variant	Copy Number	Sample Size Per Copy Number [†]	Concordance to Benchmark, % [‡]	Discordant Calls By Nala PGx Core™ [§]
<i>CYP2D6</i>	Intron 2 N = 224	0	0	99.6	0
		1	12		0
		2	130		0
		>=3	82		1
	Exon 9 N = 224	0	0	98.7	0
		1	24		1
		2	193		2
		>=3	7		0

[†] Sample Size Per Copy Number[†] reflects the counts of each copy number call as produced by the benchmark, and not by Nala PGx Core™

[‡] Refer to Table 1 for an outline of the benchmark methods used per variant.

[§] Discordant calls refer to instances in which Nala PGx Core™ did not produce a copy number call concordant with the benchmark.

Table 8 Diplotype concordance for *CYP2C9*, *CYP2C19* and *CYP2D6* between Nala PGx Core™, and Agena VeriDose® Core and *CYP2D6* CNV Panel

Genes	Concordance, % [†]	Discordant Diplotypes		
		Agena VeriDose® Core and <i>CYP2D6</i> CNV Panel	Nala PGx Core™	Instances
<i>CYP2C9</i>	100% (n = 213)	NA	NA	NA
<i>CYP2C19</i>	96.4% (n = 223)	*1/*1	*1/*3	1
		*1/*2	*1/*1	1
		*1/*2	*1/*3	1
		*1/*3	*2/*2	1
		*2/*2	*1/*1	2
		*2/*2	*1/*2	2
<i>CYP2D6</i>	94.7% (n = 209)	*1/*1	*2/*2	1
		*1/*10	*1/*10, CN >=3	1
		*1/*41	*39/*41	1
		*2/*2	*1/*1	1
		*2/*10	*1/*10	1
		*2/*36, CN >=3	*2/*10, CN >=3	1
		*4/*36	*4/*10	2
		*10/*36	*2/*36	1
		*13 *1/*10	*1/*10, CN >=3	1
		*13 *1/*41, CN = 1	*1/*41	1

[†]The concordance presented in this table excludes samples that have mismatches in diplotype calls arising from technological differences between platforms. Please refer to Supplementary Table 7 for a detailed breakdown of all diplotypes detected.

Table 9 Observed allele frequencies by ethnicity

Gene	Allele or Variant	Allele Frequencies (Per This Study)						Allele Frequencies (PharmGKB)		
		Indonesian	Chinese	Malay	Indian	Caucasian	Overall	East Asian	Central/South Asian	European
<i>CYP2C9</i> (n = 206)	*2	0.000	0.000	0.000	0.060	0.172	0.032	0.002	0.114	0.127
	*3	0.000	0.040	0.039	0.100	0.069	0.044	0.038	0.110	0.076
<i>CYP2C19</i> (n = 201)	*2	0.297	0.274	0.289	0.229	0.138	0.256	0.284	0.270	0.147
	*3	0.041	0.055	0.066	0.000	0.017	0.042	0.072	0.016	0.002
	*17	0.054	0.007	0.039	0.229	0.138	0.067	0.021	0.171	0.216
<i>CYP2D6</i> (n = 195)	*2	0.088	0.118	0.167	0.140	0.196	0.136	0.121	0.295	0.277
	*4	0.000	0.035	0.069	0.020	0.161	0.051	0.005	0.091	0.185
	*5	0.029	0.028	0.028	0.080	0.000	0.031	0.049	0.046	0.030
	*6	0.000	0.000	0.000	0.000	0.018	0.003	0.000	0.000	0.011
	*9	0.000	0.000	0.000	0.000	0.018	0.003	0.002	0.003	0.028
	*10	0.338	0.299	0.222	0.160	0.143	0.251	0.436	0.087	0.016
	*14	0.000	0.014	0.000	0.000	0.000	0.005	0.003	ND	0.000
	*29	0.000	0.000	0.000	0.000	0.018	0.003	0.000	0.003	0.001
	*36	0.206	0.250	0.181	0.160	0.054	0.190	0.012	0.000	0.000
<i>SLCO1B1</i> (n = 203)	rs4149056 [†]	0.125	0.074	0.026	0.040	0.167	0.084	0.125 [‡]	0.050 [‡]	0.159 [‡]

"ND" refers to instances in which no data is available for the given allele on PharmGKB.

[†]rs4149056 refers to the reduced function variant of *SLCO1B1* that is present in *SLCO1B1**5, *SLCO1B1**15 and *SLCO1B1**17.

[‡]Allele frequency values for rs4149056 have been obtained from gnomAD.

Figure 1. *CYP2D636 Frequency By Ethnicity.** Distribution of individuals carrying one or two copies of the *CYP2D6**36 allele among the study cohort (n=205), grouped per ethnicity




Table 10 Observed diplotype frequencies by ethnicity

Gene	Diplotype	Phenotype [‡]	Diplotype Frequencies (Per This Study) [†]											
			Indonesian		Chinese		Malay		Indian		Caucasian		Overall	
			Obs	Freq	Obs	Freq	Obs	Freq	Obs	Freq	Obs	Freq	Obs	Freq
<i>CYP2C9</i>	*1/*1	NM	39	1.000	69	0.920	35	0.921	17	0.680	16	0.552	176	0.854
	*1/*2	IM	0	0.000	0	0.000	0	0.000	3	0.120	9	0.310	12	0.058
	*1/*3	IM	0	0.000	6	0.080	3	0.079	5	0.200	3	0.103	17	0.083
	*2/*3	PM	0	0.000	0	0.000	0	0.000	0	0.000	1	0.034	1	0.005
<i>CYP2C19</i>	*1/*1	NM	10	0.270	31	0.425	14	0.368	8	0.333	14	0.483	77	0.383
	*1/*2	IM	19	0.514	28	0.384	14	0.368	5	0.208	7	0.241	73	0.363
	*1/*3	IM	3	0.081	6	0.082	3	0.079	0	0.000	1	0.034	13	0.065
	*1/*17	RM	3	0.081	1	0.014	1	0.026	5	0.208	5	0.172	15	0.075
	*2/*2	PM	1	0.027	5	0.068	3	0.079	1	0.042	0	0.000	10	0.050
	*2/*3	PM	0	0.000	2	0.027	2	0.053	0	0.000	0	0.000	4	0.020
	*2/*17	IM	1	0.027	0	0.000	0	0.000	4	0.167	1	0.034	6	0.030
	*17/*17	UM	0	0.000	0	0.000	1	0.026	1	0.042	1	0.034	3	0.015
<i>CYP2D6</i>	*1/*1	NM	6	0.176	7	0.097	5	0.139	4	0.160	2	0.071	24	0.123
	*1/*2	NM	0	0.000	2	0.028	2	0.056	3	0.120	3	0.107	10	0.051
	*1/*2, CN >=3	UM	0	0.000	0	0.000	0	0.000	1	0.040	1	0.036	2	0.010
	*1/*4	IM	0	0.000	2	0.028	1	0.028	0	0.000	4	0.143	7	0.036
	*1/*5	IM	0	0.000	2	0.028	1	0.028	1	0.040	0	0.000	4	0.021
	*1/*6	IM	0	0.000	0	0.000	0	0.000	0	0.000	1	0.036	1	0.005
	*1/*10	NM	4	0.118	3	0.042	2	0.056	2	0.080	2	0.071	13	0.067
	*1/*14	NM	0	0.000	1	0.014	0	0.000	0	0.000	0	0.000	1	0.005
	*1/*29	NM	0	0.000	0	0.000	0	0.000	0	0.000	1	0.036	1	0.005
	*1/*36	IM	3	0.088	3	0.042	5	0.139	4	0.160	0	0.000	15	0.077
	*1/*41	NM	1	0.029	2	0.028	1	0.028	0	0.000	2	0.071	6	0.031
	*1/*41, CN >=3	UM	0	0.000	0	0.000	0	0.000	1	0.040	0	0.000	1	0.005
	*2/*2	NM	0	0.000	1	0.014	1	0.028	1	0.040	1	0.036	4	0.021
	*2/*4	IM	0	0.000	0	0.000	1	0.028	0	0.000	2	0.071	3	0.015
	*2/*5	IM	0	0.000	0	0.000	1	0.028	1	0.040	0	0.000	2	0.010
	*2/*10	NM	4	0.118	4	0.056	4	0.111	0	0.000	1	0.036	13	0.067
	*2/*14	NM	0	0.000	1	0.014	0	0.000	0	0.000	0	0.000	1	0.005
	*2/*36	IM	2	0.059	4	0.056	2	0.056	0	0.000	1	0.036	9	0.046
	*2/*41	NM	0	0.000	4	0.056	0	0.000	0	0.000	1	0.036	5	0.026
	*4/*4	PM	0	0.000	1	0.014	1	0.028	0	0.000	1	0.036	3	0.015
	*4/*5	PM	0	0.000	0	0.000	0	0.000	1	0.040	0	0.000	1	0.005
	*4/*9	IM	0	0.000	0	0.000	0	0.000	0	0.000	1	0.036	1	0.005
	*4/*10	IM	0	0.000	1	0.014	1	0.028	0	0.000	0	0.000	2	0.010
	*5/*10	IM	0	0.000	2	0.028	0	0.000	1	0.040	0	0.000	3	0.015
	*5/*41	IM	2	0.059	0	0.000	0	0.000	0	0.000	0	0.000	2	0.010
	*10/*10	IM	3	0.088	1	0.014	2	0.056	0	0.000	2	0.071	8	0.041
	*10/*10, CN >=3	IM	0	0.000	1	0.014	0	0.000	0	0.000	0	0.000	1	0.005
	*10/*36	IM	9	0.265	28	0.389	5	0.139	4	0.160	1	0.036	47	0.241
*10/*41	IM	0	0.000	1	0.014	0	0.000	1	0.040	0	0.000	2	0.010	
*36/*41	IM	0	0.000	1	0.014	1	0.028	0	0.000	1	0.036	3	0.015	

[†]"Obs" and "Freq" are abbreviations for "Observations" and "Frequency" respectively.

[‡]"NM", "IM", "PM" and "UM" are abbreviations for "Normal Metabolizer", "Intermediate Metabolizer", "Poor Metabolizer" and "Ultrarapid Metabolizer" respectively.

B. Supplementary Information

Supplementary Table 1. Observed Genotype-Level Call Rates Per Variant Per Gene Per Platform, With Counts Of Successful And No Calls

Supplementary Table 2. Observed Diplotype-Level Call Rates Per Gene Per Platform, With Counts Of Successful And No Calls

Supplementary Table 3. Precision Analysis Of Nala PGx Coretm

Supplementary Table 4. Observed allele frequencies by ethnicity

Supplementary Table 5. *CYP2D6* Intron 2 And Exon 9 Observations And Frequencies By Ethnicity

Supplementary Table 6. *CYP2D6**36 Observations And Frequencies By Ethnicity

Supplementary Table 7. Diplotype Calls Observed Across All Samples Tested On Agena Veridose® Core And *CYP2D6* CNV Panel And Nala PGx CoreTM

Supplementary Table 8. Raw And Reported Diplotype Concordance Values For Samples Tested On Agena Veridose® Core And *CYP2D6* CNV Panel And Nala PGx CoreTM

Supplementary Tables 9.1 to 9.18. Raw genotype calls for rs1799853, rs1057910, rs4244285, rs4986893, rs12248560, rs1065852, rs5030655, rs3892097, rs35742686, rs16947, rs28371725, rs1135840, rs769258, rs5030865, rs5030656, rs59421388 and rs267608319, on Nala PGx CoreTM and the corresponding benchmark methods, Agena VeriDose® Core and *CYP2D6* CNV Panel or TaqMan® DME Genotyping Assays

Supplementary Table 9.19. Raw variant calls for *CYP2D6* Intron 2 and Exon 9 on Agena VeriDose® Core and *CYP2D6* CNV Panel and Nala PGx CoreTM

Supplementary Figure 1. Observed *CYP2D6* Functional Copy Number Frequencies By Ethnicity. Distribution of individuals carrying one, two or three or more copies of the *CYP2D6* functional gene among the study cohort (n=205), grouped per ethnicity.