Title: The further development of cold plasma technology: The effectiveness of a
 contactless, indirect atmospheric cold plasma method for germ reduction on surfaces in

- 3 vitro and in vivo
- 4
- 5 Authors: Ulrich Schmelz<sup>1</sup>, Tom Schaal<sup>2</sup>, Gilbert Hämmerle<sup>3</sup>, Tim Tischendorf<sup>\*2</sup>
- <sup>6</sup> <sup>1</sup> University of Fulda, Fulda, Germany
- 7 <sup>2</sup> Faculty of Health and Healthcare Sciences, University of Applied Sciences Zwickau,
- 8 Zwickau, Germany
- <sup>9</sup> <sup>3</sup>DGKP, ZWM, AZWM, Head of Wound Outpatient Clinic LKH Bregenz
- <sup>10</sup> \*Corresponding author: Tim Tischendorf (TT), E-Mail: tim.tischendorf@fh-zwickau.de, Tel.:
- 11 +49 375 5363407
- 12

## 13 Authors' Contributions:

- 14 US: Writing original draft; Formal Analysis; Investigation; Project administration;
- 15 Conceptualization
- 16 TS: Writing original draft; Formal Analysis; Writing review & editing
- 17 GH: Writing original draft; Writing review & editing
- 18 TT: Writing original draft; Visualization; Data curation

NOTE: This preprint reports new research that has not been certified by peer review and should not be used to guide clinical practice.

#### Abstract 19

20 This study investigates the antimicrobial potential of an indirect cold plasma method for the treatment of wounds. Indirect plasma methods differ from direct methods in that the cold 21 plasma does not come into direct contact with the surface to be treated. The indirect plasma 22 method described here has been implemented in the PLASMOHEAL device. The device 23 generates an aerosol of liquid particles, which is conditioned with plasma reaction products 24 and passed over the areas to be treated without contact. In vitro tests show a significant 25 germ reduction of 3.4 to 4.5 log levels against various microorganisms. In vivo tests on 26 volunteers demonstrate a reduction in E. coli contamination of 4.06 to 5.15 log levels. These 27 28 results show that indirect plasma methods can achieve equivalent effects to direct methods. The highly effective, pain-free treatment at moderate costs make the indirect plasma 29 30 method a promising option in modern wound care.

31

#### 32 **Keywords**

Cold plasma; Wound treatment; Chronic wounds; Germ reduction 33

34

#### Introduction 35

36 What is a physical "plasma"?

"Plasma" in the sense of the physical definition, represents the state of a conductive gas. 37 In the plasma state, a gas that is practically non-conductive under normal conditions is 38 partially ionized by an ignition pulse. The gas ions are then electrically charged particles, 39 whereby the gas assumes the state of an ionic conductor (2nd order conductor). This state 40 is maintained by a further supply of electrical energy. 41

The ionization of gas molecules leads to a transfer of electrical energy to the gas molecules. 42

43 The ionized gas molecules can transfer this energy, for example to surfaces. As a result,

surfaces or even liquids or the surrounding air can be modified [1]. 44

Plasma can occur at atmospheric, elevated or reduced pressure. Plasma can also be 45 altered by the addition of heat energy [2]. Furthermore, plasmas are classified according to 46

the reacting gases; in the process currently under consideration, oxygen and water vapor
are brought to react as reactants in the plasma.

The process considered in this paper takes place at room temperature and atmospheric pressure, hence the term 'atmospheric cold plasma'. The gas used for the plasma reactions is ambient air, which contains approximately 21 percent oxygen and water vapour up to the (temperature-dependent) condensation limit. Accordingly, the ambient air supplies the gases (oxygen and water vapor), which are required for the plasma reaction [3]. Plasma reactions can be direct or indirect reactions (Figure 1).

55

56 <<Figure 1. Forms of plasma reactions: direct and indirect plasma reactions.>>

57

58 Until 2024, published studies have primarily documented direct plasma procedures or a combination of direct and indirect procedures in the field of wound sanitation [4]. The 59 60 antimicrobial effect of indirect cold plasma methods on surfaces was basically shown in 2024 [5]. In a direct plasma reaction, the direct contact of the plasma with a surface leads 61 to a modification of this surface. Energy is only required for the plasma reaction (Figure 1: 62 direct method: step 1 + step 2; indirect method: step 1). Direct methods are operating in 63 64 consideration of the surface. So direct plasma reactions are in contact to the surface needs to treat. 65

Additional steps are found with indirect plasma methods. In the case of indirect plasma reactions, additional steps are added that decouple the plasma reaction from the surface to be treated. Due to decoupling, the plasma reaction always takes place at a distance from the surface to be treated. This means that the energy required for the plasma reaction can be outsourced to an external device (Figure 1: step 1; indirect method).

In an indirect plasma reaction, the process takes place outside the relevant surface. Such an indirect plasma method is used, for example, in the "PLASMOHEAL" device (WK-MedTec GmbH, Bückeburg). The device produces the cold plasma reaction products, primarily hydroxyl radicals, which are generated by a targeted flow of ambient air via a

corresponding plasma source [3]. These are generated on the basis of the defined potential
difference (voltage; here: 1.45 kV as the effective value of the alternating voltage), the
frequency (38 kHz) and the frequency form (sinusoidal frequency). Ozone is only present
in very small traces.

Then distilled water is nebulized using an ultrasonic module. Ambient air which has been enriched by cold plasma products activate the nebulized water to a cold plasma aerosol (CAP). The plasma reaction products of the air lead to an increase in the electrophysical potential of the nebulized liquid without any material change [2] (Figure 2).

The aerosol is then passed over the areas to be disinfected (here: skin, in particular wounds) for three minutes at a distance of 7.5 cm.

85

86 <<Figure 2. Functional diagram of an indirect cold plasma method using the example of</p>

PLASMOHEAL. Note. Distilled water is used as the process water for aerosol formation.

4.5 liters of aerosol contain 1 mL of process water as an aerosol component.>>

89

90 Effect of direct and indirect plasma methods against microorganisms on surfaces:

Both direct and indirect plasma methods are surface-active methods. This focuses on the transient flora, which is preferably found on the surface of the skin, mucous membranes or wounds, while the resident flora (local flora, microbiome) is also localized in deeper layers. The plasma reaction products are effective on the surface of the skin and therefore not in the tissue. In the case of wounds, the surface is interrupted so that the method could also be effective in the contact area of the aerosol in superficial horny layers.

97 Cold plasma products (contain in aerosols or in direct reaction) support the disruption of 98 the redox processes of microorganisms. Specifically, a "short circuit" of the cell membrane 99 occurs, which leads to a collapse of the resting potential between the extracellular and 100 intracellular space, whereby all transmembrane transport processes in the microorganism 101 are damaged and the microorganism is destroyed [4]. As a result, there is generally no

longer any potential for infection. The shown short circuit can be described as a "selective
 short circuit", as the human cells are practically not affected (Figure 3).

104

105 <<Figure 3. Cold plasma eliminates microorganisms by a short circuit of the membrane 106 potential.>>

107

108 Absence of toxicological effects on human tissues and cells:

109 The effects occurring in microorganisms are not observed in higher organisms. Toxicologically relevant samples or by-products are not found in indirect plasma processes 110 111 [6]. To substantiate this, toxicity assessments can provide support. As an example, toxicity tests such as the Ames test (according to EN ISO 10993-3) and a cytotoxicity test 112 113 (according to EN ISO 10993-5) were performed for the CAP on which the device is based. Both assessments showed the absence of mutagenic effects in human cells. Furthermore, 114 115 no adverse effects were observed with multiple clinical applications of direct cold plasma systems in wound treatment. The presence of antioxidant enzymes (catalase, superoxide 116 dismutase, alpha-1-antitrypsin) in tissues and organisms that effectively neutralize plasma 117 reaction products before damage to the cell membrane of tissue cells is assumed [6]. It is 118 119 important to note that there are no recognizable acute or chronic toxic effects on higher organisms (eukaryotic tissues: plants, animals, humans). This lack of effect can be 120 attributed to the presence of beforehand mentioned antioxidant enzymes (catalase, 121 superoxide dismutase, alpha-1-antitrypsin) in these tissues and organisms [6]. 122

Hence, the indirect plasma process appears to be able to eliminate microorganisms without
 harming humans, thereby reducing the risk of infection, which has already been proven in
 comparison with direct plasma processes.

126

127 Proof of efficacy of direct plasma methods in the field of wound rehabilitation:

A germ-reducing effect of the direct plasma is generally recognized. In studies in vivo (on 18 patients), a reduction of germs in the wound was observed in swabs, regardless of the

localization and size of the wound, as well as previous illnesses and cofactors (e.g. tobacco
abuse, cardiovascular diseases) [7].

In this study, microbiocidal activity against both Gram-positive cocci bacteria and Gram-132 133 negative rod bacteria (both groups include the most common human infectious agents) is sufficiently demonstrated. At the same time, Daeschlein et al. have shown that there are 134 135 no differences in potency regarding inhibitor-resistant variants of certain microorganisms 136 (e.g. MRSA, MRGN). The microorganism is damaged "in toto" by the electrophysical 137 reaction, which means that a mutation of metabolic pathways has no influence on the microbiocidal potency [7]. Furthermore, Zimmermann et al. describe an in vitro virucidal 138 139 effect against adenoviruses in the range of approximately 3 log levels [8].

The evaluation of the results of Daeschlein et al. and Zimmermann et al. leads to the conclusion that a microbiocidal effect of > 3 log levels can be expected from the direct plasma methods.

The germ-reducing effect is also indirectly demonstrated by the observation that the pH value of the wound drops when treated with direct cold plasma methods. The typical bacterial wound infection pathogens result in an increase in the pH value of the wound to a range of 8 to 10. The subsequent killing of these alkalizing microorganisms as a result of the cold plasma method then brings about the pH reduction required for optimal wound healing [9].

In addition, a considerable acceleration of wound healing was observed with the use of
direct plasmas [10; 11; 12; 13, 14].

According to Strohhal et al. 2022, it is shown that wounds treated with plasma products have a higher probability of healing (approx. 60 per cent higher compared to supportive therapy alone). The individual localization and genesis of the wound is practically irrelevant. Furthermore, application studies with immunological skin diseases (using the example of rosacea) have shown that the method is also forward-looking in the treatment of such immunological diseases [15]. Finally, it should be noted that the claimed effectiveness of the direct methods can be attributed to the antimicrobial effect. The studies of the

effectiveness of the direct methods would therefore also be applicable to the indirect methods.

160

161 Derivation of the research question:

The antimicrobial potential of the CAP could be used for specialized disinfection, for example of wounds. The cold plasma disinfection methods for wounds used to date - with acceptable results - are almost exclusively direct-acting methods with a direct plasma jet or plasma field. If an equivalent microbiocidal efficacy of indirect plasma methods to direct plasma methods can be demonstrated, it can be concluded that indirect plasma methods also have an equivalent effect in the field of wound decontamination. This leads to the following research question: Does the indirect cold plasma process have

at least an equivalent antimicrobial effect compared to direct cold plasma processes?

170 To provide this proof, the microbiocidal efficacy of an indirect plasma method is tested in

vitro and in vivo in accordance with existing test standards.

172

### 173 Methods

174 Methodology of *in vitro* testing of the microbiocidal effect:

175 An in vitro test of the germ reduction performance against five normatively defined microorganisms (Staph. aureus ATCC 6538, Staph. epidermidis ATCC 14990, E. coli 176 NCTC 19538, Pseudomonas aeruginosa ATCC 10145, Candida albicans ATCC 10321, 177 BS) was carried out in accordance with DIN spec. 91315:2014. The identity of the test 178 organism was checked by using Gram staining in conjunction with transmitted light 179 microscopy and by determining the metabolic performance of the test organism in relation 180 to various substrates. Furthermore, the test organism was identified as indole-forming from 181 182 tryptophan (indole-positive). The universal medium Caso agar in 60 mm Petri dishes was used as a culture medium for cultivating the inoculates of the samples. After checking their 183 identity, the microorganisms were cultivated as a surface culture and washed with sterile 184 and pyrogen-free NaCl solution 0.9 per cent to obtain a microbial suspension with a turbidity 185

corresponding to McFarland Standard eight. The germ suspensions were then applied to 186 187 stainless steel test specimens (100\*10\*1mm) with a roughness depth of 100µm using a 188 sterile swab. The stainless steel test specimens were prepared by immersion for five seconds in the test germ suspension, draining and air drying for two hours. These prepared 189 test specimens were then exposed to the PLASMOHEAL process. The test specimen was 190 191 exposed to the aerosol of the device for three minutes. The microorganisms were then 192 washed off by resuspending them in 10 mL of sterile and pyrogen-free NaCl solution 0.9 193 per cent in a sterile test tube by shaking for 30 seconds using a VORTEX 3 shaker. From this suspension (undiluted, dilution factor therefore 10<sup>o</sup>) a geometric dilution series (in 194 195 powers of ten) up to 10<sup>-3</sup> is created. A volume of 0.1 mL is taken from each dilution stage and applied to the culture medium suitable for the respective test organism. The volume of 196 197 0.1mL is homogeneously distributed on the culture medium using a Drigalski spatula (Figure 4). 198

199

200 <<Figure 4. Overview geometric dilution row.>>

201

In parallel, an untreated test specimen is also analyzed for each test organism (Table 1).

Table 1. Microorganisms used for the in vitro test in accordance with DIN spec. 91315.

Microorganisms	Germ identification	Culture media with identification number	Incubation time and temperature
Staphylococcus aureus	ATCC 6538	Mannitol-salt-agar No.1050	48h / 36°C
Staphylococcus epidermidis	ATCC 14990	Caso-agar No. 4020	48h / 36°C
Escherichia coli	NCTC 19538	Colichromer-agar No. 4028	48h / 36°C
Pseudomonas aeruginosa	ATCC 10145	Cetrimid-agar No. 4025	48h / 36°C
Candida albicans	ATCC 10321	Sabouraud-agar 2% No. 4130 <sup>1</sup>	72h / 28°C

Note. <sup>1</sup>The microbiological culture media for this analysis were obtained from Dr. Möller & Schmelz, Göttingen, Germany.; ATCC = American Type Culture Collection. NCTC = National Collection of Type Cultures.

204

After the incubation time required for the respective test germ in the incubator, the colonies

of the culture medium of the dilution stage containing ten to 100 colonies are counted.

207 Considering the dilution level, the fact that 0.1mL of 10mL per dilution level was applied to 208 the culture medium and the initial volume of 10mL (in which the washout took place), the 209 bacterial count per test specimen is calculated. After forming the decadic logarithm of the 210 bacterial count, the log of the bacterial count after exposure is subtracted from the logarithm 211 of the bacterial count without exposure. This is used to calculate the log reduction factor, 212 which indicates by how many powers of ten an initial bacterial load is reduced by the 213 method under consideration.

214

215 Methodology of *in vivo* testing of the microbiocidal effect:

216 Based on the results of the in vitro test, the germ reduction performance was tested in the in vivo test in accordance with EN 1500 using the test germ Escherichia coli K12. EN 1500 217 218 is used for microbiocidal testing of hand sanitizers and therefore allows in vivo testing on test subjects. The EN 1500 standard on which the test is based describes the wetting of 219 220 the fingertips as the location of the wash-off, because these can be mobilised mechanically by rubbing with the opposing thumb. The hands were immersed in the test germ suspension 221 for five seconds in the area of the fingers, then drained and loaded with the test germ for 222 three minutes. Before and after application of the PLASMOHEAL method, the bacterial 223 224 count of the skin of the hands was methodically determined in accordance with EN 1500, with both hands contaminated by the bacterial suspension in the area of the fingers. The 225 "pre-value" without the use of the device was determined from the left hand, then the right 226 hand was exposed to the aerosol in the area of the contaminated fingers for three minutes. 227 The "after value" was then determined from the exposed right hand. The dilution series and 228 cultivation were carried out analogue to the method described above for the in vitro test of 229 the microbiocidal effect. After cultivation and evaluation, the log reduction factor of the cold 230 231 plasma process when applied to human skin was determined. For evaluation, the test result of the PLASMOHEAL method was compared with the log reduction factor of standard 232 alcohol as a reference method following the analogy of EN 1500. 233

EN 1500 specifies E. coli K12 as the test germ, as this microorganism is apathogenic (it is 234 235 also a natural commensal of the human colon, among other things) and furthermore its tenacity (resistance to inactivation) includes sufficient obligate pathogenic germs, E. coli is 236 a sufficient disinfection surrogate and at the same time there is no infectiological risk. E. 237 238 coli on the hand is therefore used as a surrogate for a contaminated wound. Therefore, the performance of the disinfection efficacy test is not a clinical test in the sense of requiring 239 240 approval by an ethics committee. Rather, this test is a standard procedure, i.e. "state of the 241 art or science", which is required according to EN 1500 in order to determine the aforementioned disinfection efficacy in a technically reproducible manner. The test is 242 243 carried out in this form many times a day in Germany.

The test subjects have declared to the person carrying out the test that the test is voluntary, 244 245 furthermore there are no direct dependencies, and the test subjects also declared that they were not, did not want to be, or could not be pregnant. We confirm that all methods of this 246 247 study were conducted in accordance with the relevant guidelines and regulations. The participants have given their informed consent to participate in the study in advance. 248 Written informed consent to participate was obtained from the participants. The consents 249 were documented in a standardized laboratory list and signed by the participants. No 250 251 minors were included. A written exemption from ethics approval was obtained from the institutional review board, as the test was carried out in accordance with the requirements 252 of the EN 1500 test standard and its effectiveness was confirmed in laboratory tests. 253

254

### 255 Results

The aim of this study is to investigate the germ-reducing effect of the PLASMOHEAL process both in vitro and in vivo. By analyzing in detail the results of the in vitro tests on the germ-reducing effect of the PLASMOHEAL procedure (Table 2) and the in vivo tests (Table 3), the dynamic effect of this innovative wound treatment method is comprehensively examined. The germ-reducing properties of the procedure are presented in detail at various levels of testing.

#### Table 2. Results of the in vitro test of the germ-reducing effect of the PLASMOHEAL 262

263 process.

No.	V1	V2	CfU/ Agar plate	CfU/ Approach		CfU/ Proband	Log CfU/ Proband	Log Reductions faktor
				phylococcu	is aureus 3			
1.1.A-3	10,0 mL		87	8,700	10,000	87,000,000	7.94	
1.1.B-3	10,0 mL	0,1 mL	28	2,800	1	2,800	3.45	4.49
1.2.A-3	10,0 mL	0,1 mL	56	5,600	10,000	56,000,000	7.75	
1.2.B-3	10,0 mL	0,1 mL	26	2,600	1	2,600	3.41	4.33
1.3.A-3	10,0 mL	0,1 mL	78	7,800	10,000	78,000,000	7.89	
1.3.B-3	10,0 mL	0,1 mL	22	2,200	1	2,200	3.34	4.55
			Staph	ylococcus e	epidermidi	s 3min:		
1.1.A-3	10,0 mL	0,1 mL	56	5,600	10,000	56,000,000	7.75	
1.1.B-3	10,0 mL	0,1 mL	31	3,100	1	3100	3.49	4.26
1.2.A-3	10,0 mL	0,1 mL	44	4,400	10,000	44,000,000	7.64	
1.2.B-3	10,0 mL	0,1 mL	25	2,500	1	2,500	3.40	4.25
1.3.A-3	10,0 mL	0,1 mL	76	7,600	10,000	76,000,000	7.88	
1.3.B-3	10,0 mL	0,1 mL	37	3,700	1	3,700	3.57	4.31
				Escherichia	a coli 3min	:		
1.1.A-3	10,0 mL	0,1 mL	92	9,200	10,000	92,000,000	7.96	
1.1.B-3	10,0 mL	0,1 mL	54	5,400	1	5,400	3.73	4.23
1.2.A-3	10,0 mL	0,1 mL	78	7,800	10,000	78,000,000	7.89	
1.2.B-3	10,0 mL	0,1 mL	34	3,400	1	3,400	3.53	4.36
1.3.A-3	10,0 mL	0,1 mL	65	6,500	10,000	65,000,000	7.81	
1.3.B-3	10,0 mL	0,1 mL	39	3,900	1	3,900	3.59	4.22
	· · · · ·		Pseu	idomonas a	eruginosa	3min:		
1.1.A-3	10,0 mL	0,1 mL	72	7,200	10,000	72,000,000	7.86	
1.1.B-3	10,0 mL	0,1 mL	26	2,600	1	2,600	3.41	4.44
1.2.A-3	10,0 mL	0,1 mL	61	6,100	10,000	61,000,000	7.79	
1.2.B-3	10,0 mL	0,1 mL	34	3,400	1	3,400	3.53	4.25
1.3.A-3	10,0 mL	0,1 mL	93	9,300	10,000	93,000,000	7.97	
1.3.B-3	10,0 mL	0,1 mL	45	4,500	1	4,500	3.65	4.32
				Candida alb	icans 3mii	n:		
1.1.A-3	10,0 mL	0,1 mL	78	7,800	10,000	7,800,000	6.89	
1.1.B-3	10,0 mL	0,1 mL	31	3,100	1	3,100	3.49	3.40
1.2.A-3	10,0 mL	0,1 mL	59	5,900	10,000	59,000,000	7.77	
1.2.B-3	10,0 mL	0,1 mL	27	2,700	1	2,700	3.43	4.34
1.3.A-3	10,0 mL	0,1 mL	86	8,600	10,000	86,000,000	7.93	
1.3.B-3	10,0 mL	0,1 mL	29	2,900	1	2,900	3.46	4.47

264

Note. No. = Identification of the test run, where the first two digits stand for the test person number and 'A' = test run before disinfection (initial bacterial count), 'B' = test run after disinfection (initial bacterial count); V1 = 265 Total volume Approach to leaching; V2 = Volume of the batch from V1, which was applied to the agar plate; 266 267 CfU = Colony forming Unit.

270 process.

<sup>268</sup> 

Table 3. Results of the in vivo test of the germ-reducing effect of the PLASMOHEAL 269

 No.	V1	V2	CfU/Agar plate	CfU/	Dilution	CfU/ Proband	Log CfU/ Proband	Log Reductions factor	Mean Values
NO.	10.0	0.1	plate	Approach	Dilution	Proband	Proband	Tactor	Log RF
01 before	mL	mL	46	4,600	10,000	46,000,000	7.66		
	10.0	0.1		.,	,				
01 after	mL	mL	3	300	1	300	2.48	5.19	5.19
	10.0	0.1							
02 before	mL	mL	18	1,800	10,000	18,000,000	7.26		
00 <i>f</i>	10.0	0.1	_			500			
02 after	mL	mL	5	500	1	500	2.70	4.56	4.56
03 before	10.0 mL	0.1 mL	16	1,600	100 000	160,000,000	8.20		
	10.0	0.1	10	1,000	100,000	100,000,000	0.20		
03 after	mL	mL	14	1,400	1	1400	3.15	5.06	5.06
	10.0	0.1		.,					
04 before	mL	mL	23	2,300	1,000	2,300,000	6.36		
	10.0	0.1							
04 after	mL	mL	2	200	1	200	2.30	4.06	4.06
	10.0	0.1							
05 before	mL	mL	84	8,400	10,000	84,000,000	7.92		
05 offer	10.0	0.1	5	500	1	500	2 70	E 02	E 22
05 after	mL 10.0	mL 0.1	5	500	I	500	2.70	5.23	5.23
06 before	mL	mL	15	1,500	10,000	15,000,000	7.18		
	10.0	0.1	10	1,000	10,000	10,000,000	1.10		
06 after	mL	mL	3	300	1	300	2.48	4.70	4.70
	10.0	0.1							
07 before	mL	mL	4	400	10,000	400,000	6.60		
	10.0	0.1							
07 after	mL	mL	1	100	1	100	2.00	4.60	4.60
001.0	10.0	0.1			10.000	~~ ~~ ~~ ~~			
08 before	mL	mL	92	9,200	10,000	92,000,000	7.96		
08 after	10.0 mL	0.1 mL	17	1,700	1	1700	3.23	4.73	4.73
	10.0	0.1	17	1,700	•	1700	5.25	4.75	4.75
09 before	mL	mL	39	3,900	10,000	39,000,000	7.59		
	10.0	0.1		-,	-,	, ,			
09 after	mL	mL	5	500	1	500	2.70	4.89	4.89
	10.0	0.1							
10 before	mL	mL	24	2,400	10,000	24,000,000	7.38		
10 11	10.0	0.1		400		100			
10 after	mL	mL	1	100	1	100	2.00	5.38	5.38
11 before	10.0 mL	0.1 mL	17	1,700	10,000	17,000,000	7.23		
	10.0	0.1	17	1,700	10,000	17,000,000	1.23		
11 after	mL	mL	5	500	1	500	2.70	4.53	4.53
	10.0	0.1	-				-		
12 before	mL	mL	11	1,100	10,000	11,000,000	7.04		
	10.0	0.1							
12 after	mL	mL	1	100	1	100	2.00	5.04	5.04
	10.0	0.1							
13 before	mL	mL	26	2,600	10,000	26,000,000	7.41		
12	10.0	0.1	n	200	4	200	0.00	E 44	E 44
13 after	mL 10.0	mL 0.1	2	200	1	200	2.30	5.11	5.11
14 before	mL	mL	14	1,400	10,000	14,000,000	7.15		
14 Deluie	111L	111L	14	1,+00	10,000	14,000,000	1.10		

	10.0	0.1							
14 after	mL	mL	4	400	1	400	2.60	4.54	4.54
	10.0	0.1							
15 before	mL	mL	14	1,400	10,000	14,000,000	7.15		
	10.0	0.1							
15 after	mL	mL	1	100	1	100	2.00	5.15	5.15
	10.0	0.1							
16 before	mL	mL	39	3,900	10,000	39,000,000	7.59		
	10.0	0.1							
16 after	mL	mL	9	900	1	900	2.95	4.64	4.64
	10.0	0.1							
17 before	mL	mL	21	2,100	10,000	21,000,000	7.32		
	10.0	0.1							
17 after	mL	mL	3	300	1	300	2.48	4.85	4.85
	10.0	0.1							
18 before	mL	mL	14	1,400	10,000	14,000,000	7.15		
	10.0	0.1							
18 after	mL	mL	4	400	1	400	2.60	4.54	4.54
							Mean		
							value		
							log RF	4.82	
							Median		
							log RF	4.79	
Compariso	on contr	ol with s	standard al	cohol acco	ording to E	N 1500 (60%)	n-propanol,	40% water)	
	10.0	0.1							
19 before	mL	mL	9	900	10,000	9,000,000	6.95		
	10.0	0.1							
19 after	mL	mL	1	100	1	100	2.00		
	10.0	0.1							
19 after	mL	mL	3	300	1	300	2.48		
	10.0	0.1							
19 after	mL	mL	1	100	1	100	2.00		
	10.0	0.1							
19 after	mL	mL	4	400	1	400	2.60		
	10.0	0.1							
19 after	mL	mL	2	200	1	200	2.30	2.28	4.68

Note. No. = Identification of the test run; V1 = Total volume Approach to leaching; V2 = Volume of the batch from V1, which was applied to the agar plate; CfU = Colony forming Unit; RF = Reduction factor.

271

#### Discussion 272

The results obtained show a microbiocidal reduction performance of the PLASMOHEAL 273 process tested in vitro against the test germs Staph. aureus, Staph. epidermidis, E. coli, 274 Pseudomonas aeruginosa and Candida albicans of 3.4 to 4.5 log levels. When considering 275 the reduction factors of all five microorganisms tested, a mean value of 4.28 and a median 276 of 4.32 log levels were calculated. 277

The tested method thus achieves a germ reduction of > 4 powers of ten, which corresponds 278

to a reduction in the germ count density of > 99.99 percent. 279

The in vivo test with E. coli in accordance with EN 1500 was carried out on the real skin of the test subjects. Due to the similarity in bacterial colonization between skin and wounds, the skin can be regarded as a surrogate for wound testing [16]. A reduction of the initial E. coli contamination in the range between 4.06 and 5.15 log levels was determined. The mean value of log reduction factor across the results of the 18 test subjects was 4.82 log levels, while the median of this group of values was determined to be 4.79 log levels.

The targeted reduction of microorganisms is an important method of disinfection. Disinfection describes an antiseptic process that leads to a state of asepsis and is characterized by the fact that the risk of infection is largely eliminated. This is achieved by reducing the density of microorganisms, which is so low after disinfection that an infection is unlikely [17].

The German Association for Applied Hygiene (VAH) initially recommends the highest possible germ reduction. At least 3 log levels are considered necessary as a quantitative value for germ reducing processes. Depending on the measure and the purpose of use, 4 to 5 log levels are required [18].

The PLASMOHEAL process shows a reduction of 4.28 log levels against the test germs 295 Staph. aureus, Staph. epidermidis, E. coli, Pseudomonas aeruginosa and Candida albicans 296 297 when tested according to DIN spec. 91315. The in vitro test showed a mean value of 4.28 log levels of germ count reduction. This demonstrates that the disinfection effect of the 298 PLASMOHEAL process is safe under test conditions. Furthermore, the in vivo test based 299 on EN 1500 shows an average reduction performance of 4.82 log levels. This meets both 300 the VAH recommendation of 3 to 5 log levels and the standard alcohol also analyzed as a 301 reference in accordance with EN 1500, which showed a reduction of 4.68 log levels in the 302 303 current study.

304

The initially formulated research question can be answered to the effect that the sufficient germ reduction performance of indirect plasma methods is achieved and that indirect plasma methods in the field of wound decontamination can be regarded as at least

equivalent to the effect of direct plasma methods. The claimed effect of direct cold plasma
methods is the antimicrobial effect. The use of PLASMOHEAL as a proxy for an indirect
method shows the equivalence in the efficacy of indirect methods in relation to the assumed
efficacy claim of direct methods. Consequently, the studies of direct cold plasma methods
can also be transferred to indirect cold plasma methods as part of a metrological
traceability.

The energy in the indirect method is only applied in step 1 (Figure 1). As a result, the treatment is painless and non-irritating [9]. By adding the aerosol, a large surface area (up to 10 \* 10 cm) can be treated simultaneously within three minutes using the indirect method compared to the direct method.

318

The overall view of the results and the additional aspects of wound treatment with indirect plasma methods shows that the tested indirect plasma method can be used in a forwardlooking manner in wound treatment, particularly in the case of microbiological colonization. The purpose of the analysis was to determine the microbiocidal effect and its quantification, therefore the analysis focused on these parameters. The consideration of other and secondary parameters is another aspect for future analysis.

325

### 326 **References**

[1] Ehlbeck J, Schnabel U, Polák M, Winter J, Von Woedtke T, Brandenburg R, et al. Low
 temperature atmospheric pressure plasma sources for microbial decontamination. J Phys
 D Appl Phys. 2011;44(1):013002. https://doi.org/10.1088/0022-3727/44/1/013002.

[2] Lerouge S, Wertheimer MR, Yahia L. Plasma Sterilisation: A Review of Parameters,
 Mechanisms, and Limitations. Plasmas Polym. 2001;6:175-88.
 https://doi.org/10.1023/A:1013196629791.

[3] Krömker, W, Schorling, T. Vorrichtung zur Ioniserung von Umgebungsluft [Device for
 ionising ambient air] (EP 4 175 082 A1). European Patent Office 2023.

[4] Braný D, Dvorská D, Halašová E, Škovierová H. Cold Atmospheric Plasma: A Powerful
Tool for Modern Medicine. Int J Mol Sci. 2020;21(8):2932.
https://doi.org/10.3390/ijms21082932.

[5] Tischendorf T, Schaal T, Schmelz U. Study on hand disinfection in inpatient geriatric
care on the superiority of cold plasma aerosol versus alcohol-based disinfection methods
in a parallel group design. Sci Rep. 2024 Sep 17;14(1):21703. doi: 10.1038/s41598-02472524-7. PMID: 39289454: PMCID: PMC11408535.

[6] Masur K, Schmidt J, Stürmer E, von Woedtke T. Kaltes Plasma zur Heilung chronischer
Wunden [Cold plasmas for the healing of chronic wounds]. WundManagement.
2018;12:253-9.

[7] Daeschlein G, Napp M, Lutze S, Arnold A, Von Podewils S, Guembel D, et al. Skin and
wound decontamination of multidrug-resistant bacteria by cold atmospheric plasma
coagulation. J Ger Dermatol Soc. 2015;13(2):143-9. https://doi.org/10.1111/ddg.12559.

[8] Zimmermann J, Dumler K, Shimizu T, Morfill GE, Wolf A, Boxhammer V, et al. Effects
of cold atmospheric plasmas on adenoviruses in solution. J Phys D Appl Phys.
2011;44(50):505201. https://doi.org/10.1088/0022-3727/44/50/505201.

[9] Hämmerle G, Ascher S, Gebhardt L. Positive effects of cold atmospheric plasma on pH
in wounds: a pilot study. J Wound Care. 2023;32(9):530-6.
https://doi.org/10.12968/jowc.2023.32.9.530.


[10] Rached NA, Kley S, Storck M, Meyer T, Stücker M. Cold Plasma Therapy in Chronic
Wounds-A Multicenter, Randomised Controlled Clinical Trial (Plasma on Chronic Wounds
for Epidermal Regeneration Study): Preliminary Results. J Clin Med. 2023;12(15):5121.
https://doi.org/10.3390/jcm12155121.

[11] Strohal R, Dietrich S, Mittlböck M, Hämmerle G. Chronic wounds treated with cold
atmospheric plasmajet versus best practice wound dressings: a multicenter, randomized,
non-inferiority trial. Sci Rep. 2022;12(1). https://doi.org/10.1038/s41598-022-07333-x.

[12] Stratmann B, Costea TC, Nolte C, Hiller J, Schmidt J, Reindel J, et al. Effect of Cold
 Atmospheric Plasma Therapy vs Standard Therapy Placebo on Wound Healing in Patients

With Diabetic Foot Ulcers: A Randomized Clinical Trial. JAMA Netw Open. 2020;3(7).
 https://doi.org/10.1001/jamanetworkopen.2020.10411.

- [13] Isbary G, Heinlin J, Shimizu T, Zimmermann J, Morfill G, Schmidt H, et al. Successful
  and safe use of 2 min cold atmospheric argon plasma in chronic wounds: results of a
  randomised controlled trial. Br J Dermatol. 2012;167(2):404-10.
  https://doi.org/10.1111/j.1365-2133.2012.10923.x.
- [14] Schachl J, Königshofer M, Stoiber M, Socha M, Grasl C, Abart T, Michel-Behnke I,
- 370 Wiedemann D, Riebandt J, Zimpfer D, Schlöglhofer T. Cold atmospheric plasma therapy
- as a novel treatment for Berlin Heart EXCOR pediatric cannula infections. Artif Organs.
- 372 2024. doi: 10.1111/aor.14869.
- [15] Hofmeyer S, Weber F, Gerds S, Emmert S, Thiem A. A Prospective Randomised
  Controlled Pilot Study to Assess the Response and Tolerability of Cold Atmospheric
  Plasma for Rosacea. Skin Pharmacol Physiol. 2023;36(4):205-13.
  https://doi.org/10.1159/000533190.
- [16] Tomic-Canic M, Burgess JL, O'Neill KE, et al. Skin Microbiota and its Interplay with
  Wound Healing. Am J Clin Dermatol. 2020;21(Suppl 1):36-43.
  https://doi.org/10.1007/s40257-020-00536-w.
- [17] Dunkelberg H, Wedekind S. Eine neue Methode zur Wirksamkeitsprüfung von
   Sterilgutverpackungen in der Praxis [A new method of testing the effectiveness of sterile
   packaging in general practice]. Biomed Tech (Berl). 2002;47(11):290-3.
- [18] Disinfectant Commission of the VAH (ed.). Anforderungen und Methoden zur VAH Zertifizierung chemischer Desinfektionsverfahren [Requirements and methods for VAH
   certification of chemical disinfection procedures]. Chapters 10-11. Available from:
   https://vah-online.de/files/download/ebooks/eBook\_VAH\_Methoden\_Anforderungen.pdf.
- 387 Accessed 2024 Aug 17.




## **Comparison of Direct vs. Indirect Cold Plasma Methods**

### **Direct Cold Plasma Method**



### Indirect Cold Plasma Method



medRxiv preprint doi: https://doi.org/10.1101/2024.10.12.24315382; this version posted October 15, 2024. The copyright holder for this preprint (which was not certified by peer review) is the author/funder, who has granted medRxiv a license to display the preprint in

perpetuity. It is made available under a CC-BY 4.0 International license .

