

Supplementary Table 1: Internal reliability of survey question measured with Cronbach's alpha.

Survey	Question No.	alpha	Question text	
Observer	1	0.3799	<i>Did the study participant read the instructions?</i>	
	2	0.9684	<i>Did the study participant appear confident in their ability to follow the instructions?</i>	
	3	0.5534	<i>Did the study participant properly wash their hands before and after sample collection?</i>	
	4	0.4434	<i>Did the study participant appear to properly follow instructions for sample collection set up?</i>	
	5	0.8518	<i>Did the study participant appear to properly follow instructions for adequate sample collection?</i>	
	8	0.8256	<i>Did the study participant securely fasten the collection tube?</i>	
	9	0.6961	<i>Did the study participant clean down the outside of the sample tube following collection?</i>	
	10	0.8208	<i>Did the study participant properly store their sample in the biohazard bag after collection?</i>	
	11	0.4703	<i>Did the study participant appear to struggle with any particular step? If so, explain which.</i>	
	Participant	1	0.5617	<i>Did you read all the instructions carefully prior to collecting the sample?</i>
		2	0.7480	<i>Did you feel that you needed help while collecting your sample?</i>
3		0.5382	<i>Do you feel confident that you collected the sample properly?</i>	
4		0.4737	<i>Did sample collection cause you physical discomfort?</i>	
5		0.2517	<i>Was collecting the sample difficult in general?</i>	
6		0.2229	<i>Was it difficult to put the appropriate amount of saliva into the tube?</i>	
7		0.2251	<i>Did you get any saliva on the outside of the collection tube?</i>	
8		0.7217	<i>Did you know what to do if saliva came into contact with the outside of the collection tube?</i>	
9		0.9411	<i>Did you understand that if you did not follow the procedure exactly, you might get a false result?</i>	
10		0.9121	<i>Did the instructions clearly explain how to collect the sample?</i>	
Laboratory	1	0	<i>The sample was of sufficient volume (200-500 ul)</i>	
	2	1	<i>The sample was easy to pipette</i>	
	3	1	<i>The sample was normal, true saliva</i>	
	4	1	<i>The sample was free from food particles</i>	
	5	0.916	<i>The sample was not unusually discolored</i>	