

Supplementary Table I

A.

PCR Master Mix:

Master mix components	Volume per reaction (μl)/well
RT-PCR enzyme mix (125x)	0.08
SYBR™ master mix (2x)	5
4μM Forward Primer	0.5
4μM Reverse Primer	0.5
DNase/RNase free water	1.92
Total	8

+ 2 μl (1:20 dilution) of RNA template, positive control plasmid, non-template control water.

Cyclin parameters:

Cycling Stage	Step	Temperature	Time
Holding	1	48°C	30:00
	2	95°C	10:00
Cycling	3	95°C	00:15
	4	60°C	01:00
	5	Go to Step 3, 44 more times	
Melt Curve, increment 0.5°C	6	95°C	00:15
		65°C	00:15

B.

Primer name	Forward sequence (5' to 3')	Reverse sequence (5' to 3')	Target gene	Gene product
N*	GCTGCAATCGTGCTACAAC	TGAACTGTTGCGACTACGTG	<i>N</i> gene	SARS-CoV-2 nucleocapsid protein
E*	TTCGGAAGAGACAGGTACGTT	CACACAATCGATGCGCAGTA	<i>E</i> gene	SARS-CoV-2 Small envelope protein
S	GCTGGTGCTGCAGCTTATTA	AGGGTCAAGTGCACAGTCTA	<i>S</i> gene	SARS-CoV-2 Spike surface glycoprotein
CDC_RNAseP (from CDC)	AGATTTGGACCTGCGAGCG	GAGCGGCTGTCTCCACAAGT	<i>RPP30</i> gene	Human Ribonuclease P (RNase P) protein