

## Towards adaptive deep brain stimulation: clinical and technical notes on a novel commercial device for chronic brain sensing

Yohann Thenaisie<sup>a,b,\*</sup>, Chiara Palmisano<sup>c,\*</sup>, Andrea Canessa<sup>c,d</sup>, Bart J. Keulen<sup>e,f</sup>, Philipp Capetian<sup>c</sup>, Mayte Castro Jiménez<sup>g</sup>, Julien F. Bally<sup>g</sup>, Elena Manferlotti<sup>c,h</sup>, Laura Beccaria<sup>c</sup>, Rodi Zutt<sup>e</sup>, Grégoire Courtine<sup>b,i,l</sup>, Jocelyne Bloch<sup>b,i</sup>, Niels A. van der Gaag<sup>m</sup>, Carel F. Hoffmann<sup>m</sup>, Eduardo Martin Moraud<sup>a,b,#</sup>, Ioannis U. Isaias<sup>c,#</sup> and M. Fiorella Contarino<sup>e,n,#,§</sup>

§ Corresponding Author M.F. Contarino, Els Borst-Eilersplein 275-2545AA, The Hague, The Netherlands. E-mail address: m.f.contarino@lumc.nl

**Supplementary Table 2: Patients with cardiac artefacts**

Patient	IPG implant site	Modes		Artefactual channel(s)	Monopolar contact impedance
NL1	Chest left	<i>Streaming and Setup</i>		0-2L 1-3R	0L: 967 Ω 1L: 1070 Ω (Stim) 2L: 1126 Ω
		Stimulation “on” (left or right)		—	1R: 971 Ω 2R: 962 Ω (Stim)
		Stimulation “off”		—	3R: 1006 Ω
		<i>Setup*</i>		All contact pairs	
NL2	Chest left	Session 1	<i>Streaming</i>	0-3L 1-3L 0-2R 0-3R	0R: 1224 Ω 1R: 1409 Ω (Stim or OFF/0mA) 2R: 1387 Ω (Stim or OFF/0mA) 3R: 1592 Ω
			Stimulation “on” (left or right)		
		Session 3	<i>Streaming</i>	0-3R	0R: 1643 Ω, 1R: 1189 Ω (Stim or OFF/0mA)
			Stimulation “on” (right)	—	2R: 1186 Ω (Stim or OFF/0mA)
		<i>Setup*</i>		All contact pairs	3R: 1685 Ω
PW1	Chest right	<i>Streaming</i>		1-3L 1-3R	1L: 1046 Ω 2L: 1128 Ω (Stim) 3L: 1445 Ω 1R: 829 Ω 2R: 996 (Stim) 3R: 815 Ω
CH6	Chest right	<i>Survey Indefinite Streaming</i>		1-3R	1R: 802Ω 3R: 785Ω

		<i>Setup*</i>	Five out of six contact pairs (including 1-3R)	
--	--	---------------	--	--

\* In the Setup mode, contact pairs were labelled as “artefactual” by the Percept system; in the other modalities, artefacts were detected offline by visual inspection.