

GECNAV pipeline

- 1) Go to <https://gdc.cancer.gov/about-data/publications/PanCan-CellOfOrigin>
- 2) Download the file Analyte level annotations - [merged_sample_quality_annotations.tsv](#)
- 3) Using information in the column **aliquot_barcode**, delete all samples named TCGA-XX-XXXX-YYX-XXX-XXXX-XX, where YY – any number, except 01, 03 and 09 (see https://docs.gdc.cancer.gov/Encyclopedia/pages/TCGA_Barcode/ and table <https://gdc.cancer.gov/resources-tcga-users/tcga-code-tables/sample-type-codes>), and also, using information in the column **Do_not_use**, delete all samples with False value, and save the resulting file as **merged_sample_quality_annotations_do_not_use.tsv**.
- 4) Download the file gzipped ISAR-corrected GISTIC2.0 all_thresholded.by_genes file - [ISAR_GISTIC.all_thresholded.by_genes.txt](#), extract to **ISAR_GISTIC.all_thresholded.by_genes.txt** and rename to **ISAR_GISTIC.all_thresholded.by_genes.tsv**
- 5) Delete from the file **ISAR_GISTIC.all_thresholded.by_genes.tsv** all samples named TCGA-XX-XXXX-YYX-XXX-XXXX-XX, where YY – any number, except 01, 03 и 09; and all samples present in the file **merged_sample_quality_annotations_do_not_use.tsv**, and save the resulting file as **ISAR_GISTIC.all_thresholded.by_genes_primary_whitelisted.tsv**
- 6) Download the file RNA batch corrected matrix - [EBPlusPlusAdjustPANCAN_IlluminaHiSeq_RNASeqV2.geneExp.tsv](#)
- 7) Delete from the file **EBPlusPlusAdjustPANCAN_IlluminaHiSeq_RNASeqV2-v2.geneExp.tsv** all samples named TCGA-XX-XXXX-YYX-XXX-XXXX-XX, where YY – any number, except 01, 03 и 09; and all samples present in the file **merged_sample_quality_annotations_do_not_use.tsv**, and save the resulting file as **EBPlusPlusAdjustPANCAN_IlluminaHiSeq_RNASeqV2-v2.geneExp_primary_whitelisted.tsv**
- 8) Download the file miRNA batch corrected matrix - [pancanMiRs_EBadjOnProtocolPlatformWithoutRepsWithUnCorrectMiRs_08_04_16.csv](#) and convert to **pancanMiRs_EBadjOnProtocolPlatformWithoutRepsWithUnCorrectMiRs_08_04_16.tsv**
- 9) Delete from the file **pancanMiRs_EBadjOnProtocolPlatformWithoutRepsWithUnCorrectMiRs_08_04_16.tsv** all samples named TCGA-XX-XXXX-YYX-XXX-XXXX-XX, where YY – any number, except 01, 03 и 09; and all samples present in the file **merged_sample_quality_annotations_do_not_use.tsv**, and save the resulting file as **pancanMiRs_EBadjOnProtocolPlatformWithoutRepsWithUnCorrectMiRs_08_04_16_primary_whitelisted.tsv**
- 10) Using the file **EBPlusPlusAdjustPANCAN_IlluminaHiSeq_RNASeqV2-v2.geneExp_primary_whitelisted.tsv**, determine the median expression level for each gene across patients. If the expression for a given gene in a given patient is below 0.05x median value, replace it with “-2”, if between 0.05x and 0.75x median value, replace it with “-1”, if between 1.25x and 1.75x median value, replace with “1”, if above 1.75x median value, replace with “2”, otherwise replace with “0”. Save the file as **EBPlusPlusAdjustPANCAN_IlluminaHiSeq_RNASeqV2-v2.geneExp_primary_whitelisted_median.tsv**
- 11) Using the file **pancanMiRs_EBadjOnProtocolPlatformWithoutRepsWithUnCorrectMiRs_08_04_16_primary_whitelisted.tsv**, determine the median expression level for each miRNA across patients. If the expression for a given miRNA in a given patient is below 0.05x median value, replace it with “-2”, if between 0.05x and 0.75x median value, replace it with “-1”, if between 1.25x and 1.75x median value, replace with “1”, if above 1.75x median value, replace with “2”, otherwise replace with “0”. Save the file as

pancanMiRs_EBadjOnProtocolPlatformWithoutRepsWithUnCorrectMiRs_08_04_16_primary_whitelisted_median.tsv

- 12) Process the file **ISAR_GISTIC.all_thresholded.by_genes_primary_whitelisted.tsv** according to the following table:

		Gene CNA status in a given patient				
		ISAR_GISTIC.all_thresholded.by_genes_primary_whitelisted.tsv				
		-2	-1	0	1	2
Gene expression status in the same patient EBPlusPlusAdjustPANCAN_IlluminaHiSeq_RNASeqV2-v2.geneExp_primary_whitelisted_median.tsv	2	0	0	0	2	2
	1	0	0	0	1	1
	0	0	0	0	0	0
	-1	-1	-1	0	0	0
	-2	-2	-2	0	0	0
If gene not found look in the miRNA file:						
miRNA expression status in the same patient pancanMiRs_EBadjOnProtocolPlatformWithoutRepsWithUnCorrectMiRs_08_04_16_primary_whitelisted_median.tsv	2	0	0	0	2	2
	1	0	0	0	1	1
	0	0	0	0	0	0
	-1	-1	-1	0	0	0
	-2	-2	-2	0	0	0
If still not found then:		-2	-1	0	1	2

Save the results as

ISAR_GISTIC.all_thresholded.by_genes_primary_whitelisted_RNAfiltered.tsv