

Supplemental Table S6: Estimation of SARS-CoV2 genome copy numbers for Ct test samples. Concentrations of the samples in SARS-CoV-2 genome copies per ml were calculated using internal data supplied by TIB MolBiol for the LightMix® Modular Sarbecovirus SARS-CoV-2 Kit, estimating 1000 viral genome copies per reaction at Ct30 and extrapolating therefrom. Numbers were corrected for the volume used in RT-qPCR (10µl) and the input/ eluate ratio of the RNA extraction process, and are given in logarithmic units base 10.

Ct test sample	Copies per ml (corrected for extraction and RT-qPCR reaction volume)	log10 [SARS-CoV-2 genome copies/ ml]
16	702 171 428	8,9
21	21 942 857	7,3
25	1 371 428	6,1
28	171 428	5,2
Reference: 30	42 857	4,6