

**Title: Tracking the temporal variation of COVID-19 surges through wastewater-based epidemiology during the peak of the pandemic: a six-month long study in Charlotte, North Carolina**

Authors: Visva Bharati Barua<sup>a\*</sup>, Md Ariful Islam Juel<sup>a\*</sup>, A. Denene Blackwood<sup>b</sup>, Thomas Clerkin<sup>b</sup>, Mark Ciesielski<sup>b</sup>, Adeola Julian Sorinolu<sup>a</sup>, David A. Holcomb<sup>c</sup>, Isaiah Young<sup>a</sup>, Gina Kimble<sup>d</sup>, Shannon Sypolt<sup>d</sup>, Lawrence S. Engel<sup>c</sup>, Rachel T. Noble<sup>b</sup>, Mariya Munir<sup>a</sup>

<sup>a</sup>Department of Civil and Environmental Engineering, University of North Carolina Charlotte, 9201 University City Boulevard, Charlotte, NC 28223, USA

<sup>b</sup>Institute of Marine Sciences, The University of North Carolina at Chapel Hill, Morehead City, NC 28557, USA

<sup>c</sup>Department of Epidemiology, Gillings School of Global Public Health, The University of North Carolina at Chapel Hill, Chapel Hill, NC 27599, USA


<sup>d</sup>Charlotte Water, 5100 Brookshire Blvd., Charlotte, NC 28216, USA

Corresponding author: [mmunir@uncc.edu](mailto:mmunir@uncc.edu) (Mariya Munir)


\* These authors contributed equally

**Supplementary file**

(a)



(b)



**Fig. S1.** Standard curve for (a) N1 and (b) N2 assay using workflow 1.

**Table S1.** Primers and probe sequences utilized in this WBE study.

Assay	Primer/Probe	Sequences	Reference
CDC N1	Forward	GACCCCAAATCAGCGAAAT	<sup>1</sup> CDC (2020)
	Reverse	TCTGGTACTGCCAGTTGAATCTG	
	Probe	FAM-ACCCCGCATTACGTTTGGTGGACC-BHQ1	
CDC N2	Forward	TTACAAACATTGGCCGCAA	<sup>1</sup> CDC (2020)
	Reverse	GCGCGACATTCCGAAGAA	
	Probe	FAM-ACAATTTGCCCCCAGCGCTTCAG-BHQ1	
BCoV	Forward	CTGGAAGTTGGTGGAGTT	<sup>2</sup> Decaro et al., 2008
	Reverse	ATTATCGGCCTAACATACATC	
	Probe	FAM-CCTTCATATCTATACACATCAAGTTGTT-BHQ1	

<b>Hep G</b>	Forward	CGG CCA AAA GGT GGT GGA TG	<sup>3</sup> Schlueter et al., 1996
	Reverse	CGA CGA GCC TGA CGT CGG G	
	Probe	HEX-AGG TCC CTC TGG CGC TTG TGG CGA G-BHQ1	

**Table S2.** Average recovery efficiency of BCoV and HepG for workflow 1

<b>Workflow 1</b>		
<b>Plant</b>	<b>BCoV recovery (%)</b>	<b>Hep G recovery (%)</b>
<b>A</b>	21.5±15	41.31±18.1
<b>B</b>	22±23	43.4±21
<b>C</b>	20.2±24	38±22.7
<b>D</b>	30.9±29	44±17.4

**Table S3.** Average recovery efficiency of BCoV and HepG for workflow 2

<b>Workflow 2</b>		
<b>Plant</b>	<b>BCoV recovery (%)</b>	<b>Hep G recovery (%)</b>
<b>A</b>	24.76±20	26.22±15.5
<b>B</b>	22.3±17	17.3±11.5
<b>C</b>	23.3±18	24±13.6
<b>D</b>	26±25	29.8±19.6

**Table S4.** Inhibition control data for N1 gene for workflow 1 (RT-qPCR)

Date	Plant	Undiluted				Diluted (1:5)			
		Cq				Cq			
		Rep#1	Rep#2	Rep#3	Avg. Cq	Rep#1	Rep#2	Rep#3	Avg. Cq
11/23/2020	A	33.67	31.89	32.34	32.63	35.54	ND	36.29	35.9
	B	34.24	33.37	33.92	33.77	ND	ND	ND	NA
	C	34.7	32.96	34.94	34.20	ND	ND	ND	NA
11/25/2020	A	33.24	33.27	33.33	33.28	ND	ND	ND	NA
	B	36.07	33.31	33.92	34.43	ND	ND	ND	NA
	C	32.00	32.54	32.01	32.18	36.85	ND	ND	NA
	D	ND	ND	ND	NA	ND	ND	ND	NA
11/30/2020	A	32.9	32.92	32.19	32.67	36.00	ND	ND	NA
	B	31.87	33.3	33.07	32.74	ND	ND	35.44	35.44
	C	32.33	32.1	32.24	32.22	36.24	ND	ND	36.24
Date	Plant	Undiluted				Diluted (1:3)			
		Cq				Cq			
		Rep#1	Rep#2	Rep#3	Avg. Cq	Rep#1	Rep#2	Rep#3	Avg. Cq
7/1/2020	A	33.93	32.92	-	33.42	33.87	35.55	35.87	35.10
	B	ND	ND	-	NA	ND	ND	ND	NA
	C	34.17	35.16	-	34.66	34.41	36.47	ND	35.44
	D	34.61	32.52	-	33.57	33.19	ND	34.74	33.97
7/8/2020	B	37.81	35.45	-	36.63	37.80	ND	ND	37.63
	C	36.22	34.36	-	35.29	ND	35.97	ND	35.97

	D	ND	ND	-	NA	ND	ND	ND	NA
7/22/2020	A	33.07	32.10	-	32.59	ND	33.01	34.02	34.81
	B	36.83	36.12	-	36.47	ND	ND	ND	NA
	C	37.21	32.47	-	34.84	ND	35.28	ND	35.28
	D	37.05	34.71	-	35.88	ND	ND	ND	NA
7/29/2020	A	ND	ND	-	NA	ND	ND	ND	NA
	B	33.73	ND	-	NA	ND	ND	ND	NA
	C	35.12	33.69	-	34.41	35.96	ND	ND	35.96
	D	ND	ND	-	NA	35.36	35.29	33.94	34.86

**Table S5.** Inhibition control data for N2 gene for workflow 2 (ddPCR)

Date	Plant	Undiluted				Diluted (1:3)			
		Cq				Cq			
		Rep#1	Rep#2	Rep#3	Avg. Cq	Rep#1	Rep#2	Rep#3	Avg. Cq
7/1/2020	A	35.23	34.22	37.43	35.63	35.29	34.76	ND	35.03
	B	ND	ND	ND	NA	ND	ND	ND	NA
	C	ND	37.15	36.06	36	34.61	ND	37.33	35.97
	D	ND	41.17	39.34	40.25	35.53	ND	33.85	34.69
7/8/2020	B	ND	ND	ND	NA	ND	ND	ND	NA
	C	ND	ND	35.39	35.39	ND	ND	ND	NA
	D	ND	ND	ND	NA	ND	ND	ND	NA
7/22/2020	A	33.60	33.79	34.33	35.91	ND	35.94	ND	35.94
	B	ND	ND	ND	NA	ND	ND	ND	NA
	C	33.97	36.93	34.29	35.06	ND	35.94	ND	35.94

	D	ND	ND	ND	NA	ND	36.74	ND	36.74
7/29/2020	A	36.75	ND	ND	36.75	ND	ND	ND	NA
	B	ND	ND	ND	NA	ND	ND	39.80	39.80
	C	36.46	34.16	36.52	35.71	ND	ND	36.43	36.43
	D	ND	ND	ND	NA	37.93	ND	ND	NA

Date	Plant	Undiluted				Diluted (1:2)			
		Cq		Avg. Cq		Cq		Avg. Cq	
7/15/2020	A	ND	ND	ND	NA	ND	ND	ND	NA
	C	ND	ND	ND	NA	ND	ND	ND	NA
8/19/2020	C	ND	ND	ND	NA	ND	ND	ND	NA
	D	ND	ND	ND	NA	ND	ND	ND	NA
8/26/2020	A	ND	ND	ND	NA	ND	ND	ND	NA
	B	ND	ND	ND	NA	ND	ND	ND	NA
	C	ND	ND	ND	NA	ND	ND	ND	NA
9/23/2020	A	ND	ND	ND	NA	ND	ND	ND	NA
	B	ND	ND	ND	NA	ND	ND	ND	NA
	C	ND	ND	ND	NA	ND	ND	ND	NA
	D	ND	ND	ND	NA	ND	ND	ND	NA
9/30/2020	C	ND	ND	ND	NA	ND	ND	ND	NA
10/7/2020	A	ND	ND	ND	NA	ND	ND	ND	NA
	B	ND	ND	ND	NA	ND	ND	ND	NA
	C	ND	ND	ND	NA	ND	ND	ND	NA
	D	ND	ND	ND	NA	ND	ND	ND	NA

11/23/2020	A	ND	ND	ND	NA	33.22	33.65	ND	33.44
	B	ND	ND	ND	NA	35.55	34.52	34.04	34.70
	C	ND	ND	ND	NA	33.51	33.41	ND	33.46
11/25/2020	A	ND	ND	ND	NA	34.26	33.78	33.9	33.98
	B	ND	ND	ND	NA	ND	35.16	35.67	35.42
	C	ND	ND	ND	NA	33.69	34.13	33.38	33.73
	D	ND	ND	ND	NA	ND	ND	ND	NA
11/30/2020	A	ND	ND	ND	NA	ND	33.66	33.5	33.58
	B	ND	ND	ND	NA	33.78	33.74	ND	33.76
	C	ND	ND	ND	NA	ND	33.16	33.97	33.57
12/2/2020	A	ND	ND	ND	NA	31.37	31.21	32.03	31.54
	B	ND	ND	ND	NA	32.93	33.51	33.01	33.15
	C	ND	ND	ND	NA	33.64	33.1	ND	33.37
	D	ND	ND	ND	NA	ND	ND	ND	NA
12/7/2020	A	ND	ND	ND	NA	32.69	33.15	32.96	32.93
	B	ND	ND	ND	NA	34.18	33.82	ND	34
	C	ND	ND	ND	NA	ND	33.18	35.03	34.12

- Not detected; NA- Not applicable

**Table S6.** Inhibiton control data for workflow 2 (ddPCR)

<b>Date</b>	<b>Inhibition efficiency (%)</b>			
	<b>Plant A</b>	<b>Plant B</b>	<b>Plant C</b>	<b>Plant D</b>
6/24/2020	101	101	100	-
7/1/2020	101	101	104	99
7/8/2020	96	104	98	101
7/15/2020	102	98	100	97
7/22/2020	99	98	103	99
7/29/2020	98	99	97	98
8/5/2020	94	98	NT	98
8/12/2020	100	98	100	97
8/19/2020	101	98	97	102
8/26/2020	97	100	99	97
9/2/2020	97	96	99	94
9/9/2020	95	92	-	98
9/15/2020	-	98	89	-
9/30/2020	98	97	94	91
10/7/2020	96	90	94	95
10/14/2020	98	95	94	97



10/28/2020	96	95	101	93
11/2/2020	98	105	100	-
11/4/2020	96	93	-	-
11/11/2020	93	94	101	96
11/16/2020	93	100	94	-
11/18/2020	95	93	97	-
11/23/2020	97	97	96	-
11/25/2020	94	93	94	96
11/30/2020	96	101	97	-
12/7/2020	101	96	94	-
12/9/2020	90	95	95	91
12/14/2020	-	95	91	-

**Table S7a.** Detailed reagents and methods utilized in workflow 2

<b>Kingfisher Magnetic bead extraction</b>			
<b>Quantity</b>	<b>Description</b>	<b>Vendor</b>	<b>Catalog Number</b>
1	Tip comb	Fisher Sci.	22-387-029
7	96 well, DWP	Fisher Sci.	22-387-031
800 $\mu$ L	Wash Buffer 1	Biomerieux	280130
1000 $\mu$ L	Wash Buffer 2	Biomerieux	280131
500 $\mu$ L	Wash Buffer 3	Biomerieux	280132
50 $\mu$ L	Magnetic silica	Biomerieux	280134
100 $\mu$ L	Buffer AE	Life Technologies	19077

**Extraction protocol for workflow 2**

***Lysate Plate***

- Dispense up to 950  $\mu$ L of supernatant of lysed sample into a 96 Deep Well plate containing 50 mL of magnetic beads.
- Incubate 96 well Deep Well plate containing sample and Magnetic Beads for 10 minutes at room temperature.
- While incubating, set up other plates.

***Wash 1 Plate***

- Dispense 400  $\mu$ L of Wash Buffer 1 into each well

***Wash 2 Plate***

- Dispense 400  $\mu$ L of Wash Buffer 1 into each well

***Wash 3 Plate***

- Dispense 500  $\mu$ L of Wash Buffer 2 into each well

***Wash 4 Plate***

- Dispense 500  $\mu$ L of Wash Buffer 2 into each well

***Wash 5 Plate***

- Dispense 500  $\mu$ L of Wash Buffer 3 into each well

***Eluate Plate***

- Dispense 100  $\mu$ L of AE Elution Buffer into each well

***Tip Comb Plate***

- Place 96 DW tip comb into plate

After 10-minutes of incubation, plates were added to the Kingfisher Instrument as the program directed.

**Table S7b.** Sample preparation for Reverse Transcription

<b>Quantity</b>	<b>Description</b>	<b>Vendor</b>	<b>Catalog Number</b>
50 $\mu$ L	Superscript VILO IV MM	Fisher Sci.	11756050
1	Foil seal	BioRad	1814040
1	96 well PCR plate	BioRad	12001925
1 $\mu$ L	Mouse Lung RNA	BioChain	R1334152-50

<b>RT MasterMIX</b>	<b><math>\mu</math>L</b>
Volume per rxn	70
Superscript IV VILO 5X	14
RNA template	35
Mouse Lung RNA	1
Water	20

## References

- (1) CDC, (2020). (<https://www.cdc.gov/healthywater/surveillance/wastewater-surveillance/testing-methods.html>)
- (2) Decaro, N.; Elia, G.; Campolo, M.; Desario, C.; Mari, V.; Radogna, A.; Colaianni, M. L.; Cirone, F.; Tempesta, M.; Buonavoglia, C. Detection of Bovine Coronavirus Using a TaqMan-Based Real-Time RT-PCR Assay. *J. Virol. Methods* 2008, *151* (2), 167–171. <https://doi.org/10.1016/j.jviromet.2008.05.016>.
- (3) Schlueter, V.; Schmolke, S.; Stark, K.; Hess, G.; Ofenloch-Haehnle, B.; Engel, A. M. Reverse Transcription-PCR Detection of Hepatitis G Virus. *J. Clin. Microbiol.* 1996, *34* (11), 2660–2664. <https://doi.org/10.1128/jcm.34.11.2660-2664.1996>.