

Supplementary Table - Flow cytometry panels			
Panel 1: Phenotyping of MNPs from PBMCs			
Fluorochrome	Marker	Company	Clone
FITC	CD83	Biolegend	HB15e
PerCp Cy5.5	CD123	BD	7G3
PE Cy7	CD1c	Miltenyi	AD5-8E7
PE Cy5	CD11c	BD	B-Ly6
PE TR	HLA-DR	Life Technologies	TU36
PE	CD141	Miltenyi	AD5-14H12
APC-Cy7	CD3	BD	SK7
	CD19	BioLegend	HIB19
	CD20	BD	L27
	CD56	BD	HCD56
	CD66abce	Miltenyi	TET2
AF700	CD16	BioLegend	3GE
APC	CD45	BD	HI30
BV786	CD62L	BD	SK11
BV650	CD86	BD	2331
BV605	CCR2	BioLegend	K036C2
BV510	CD14	BD	M5E2
V450	CCR7	BD	150503
DAPI	Live/Dead Blue	Thermo Fisher	Cat no: L34962

Panel 2: Phenotyping of PBMCs following TLR stimulation- TNF-release assay			
Fluorochrome	Marker	Company	Clone
FITC	CD83	Biolegend	HB15e
PerCp Cy5.5	CD123	BD	7G3
PE Cy7	CD1c	Miltenyi	AD5-8E7
PE Cy5	CD11c	BD	B-Ly6
PE TR	HLA-DR	Life Technologies	TU36
PE	TNF-α Detection Antibody	Miltenyi	(Reagent from 130-091-268)
APC-Cy7	CD3	BD	SK7
	CD19	BioLegend	HIB19
	CD20	BD	L27
	CD56	BD	HCD56
	CD66abce	Miltenyi	TET2
AF700	CD16	BioLegend	3GE
APC	CD141	Miltenyi	AD5-14H12
BV786	CD62L	BD	SK11
BV650	CD86	BD	2331
BV605	CCR2	BioLegend	K036C2
BV510	CD14	BD	M5E2
V450	CCR7	BD	150503
DAPI	Live/Dead Blue	Thermo Fisher	Cat no: L34962