

Supplementary Table 6 Multi-FACS flow panel antibodies and reagents

Steps/mastermix	Marker/reagent	Fluorochoime	Clone
Viability stain	Live/dead blue	Live/dead blue	N/A
Block Fc receptor	Human TruStain FcX (Fc Receptor Blocking Solution)	N/A	N/A
Mastermix 1	Brilliant stain buffer plus	N/A	N/A
	CD183/CXCR3	PE-CY5	1C6/CXCR3
	CD117	BUV805	104D2
	CD294/CRTH2	BUV615	BM16
	CD161	eFluor450	HP-3G10
Mastermix 2	Brilliant stain buffer plus	N/A	N/A
	CD185/CXCR5	Brilliant Violet 750	RF8B2
	ICOS	BB515	DX29
	CCR7	Brilliant Violet 785	G043H7
	CD196/CCR6	Brilliant Violet 711	G034E3
Mastermix 3	Brilliant stain buffer plus	N/A	N/A
	CD16	PerCP-Cy5.5	3G8
	CD20	Pacific Orange	HI47
	CD95	Brilliant Violet 650	DX2
	CD27	APC	O323
	CD38	APC-eFluor780	HIT2
	CD24	PE/Dazzle™ 594	ML5
	TIGIT	PE-Cy7	A15153G
	CD25	PE-Af700	CD25-3G10
	TCRgd	PerCP-eFluor710	B1.1
	CD127	APC-R700	HIL-7R-M21
	CD45	PerCP	H130
	CD123	Super Bright 436	6H6
	KLRG1	Alexa Fluor® 647	SA231A2
	CD4	CF568	C4/206
	CD19	Spark NIR 685	HIB19
	CD8	Qdot 800	3B5
	CD57	FITC	HNK-1
	CD14	Spark blue 550	63D3
	CD39	BUV496	TU66
	CD28	BUV563	CD28.2
	PD-1	Brilliant Violet 421	EH12.2H7
	HLA-DR	Brilliant Violet 570	L243
	CD3	Brilliant Violet 510	OKT3
	CD56	BUV737	NCAM16.2
	CD31	Brilliant Violet 605	WM59
	CD11c	BUV661	B-ly6
	IgD	Brilliant Violet 480	IA6-2
	CD45RA	BUV395	5H9