

Cell Population	Subpopulation	Markers
B cells	Memory	MS4A1, COCH, AIM2, BANK1, CD79A
	Resting	IL4R, CXCR4, BTG1, TCL1A, YBX3
	Plasmablasts	MZB1, TNFRSF17, CPNE5, POU2AF1
	Activated	FOSB, CD196 (ADT)
T cells	T Proliferating	MKI67, TOP2A, PCLAF, TYMS, RRM2
	CD4 T CM	CD4, TCF7, IL7R
	CD4 T Regulatory	CD4, CTLA4, FOXP3, IL2RA
	CD4 Activated	CD4, CCL5, FYB1, GZMK, IL32, TRAC, KLRB, GZMA, ZEB2
	CD8 Naïve	CCR7, CD8A, CD8B
	CD8 EM	CD8A, CD28, CD127, CD45RO, CCL5, NKG7
T innate cells	MAIT	TRAV1-2, SLC4A10, KLRB1, GZMK, IL7R, SLC4A10, CXCR6, CD45RO (ADT)
	NK	GNLY, TYROBP, NKG7, FGFBP2, KLRF1
	$\gamma\delta$ T	TRDC, TRGC2, KLRC1, NKG7, TRDV2, TRGV9, KLRG1
	NK Proliferating	GNLY, NKG7, MKI67, TOP2A, STMN1, TYMS
Myeloid cells	pDC	CLEC4C, IL3RA
	Classical Monocytes	CD14, LYZ, CD36
	NKG7 ⁺ Monocytes	LYZ, NKG7, GNLY
	Non-classical monocytes	FCGR3A, MS4A7, CDKN1C
	NEAT1 ^{hi} Monocytes	CD14, NEAT1

Supplemental table 2. Markers used to delineate immune cell subsets.