

Table S1. Population characteristics of the dataset used to train the 3DCNN model.

Patient characteristics	Value*
Total number of patients	5351
Number of patients in HFA24-2	5227
Number of patients in HFA10-2	2047
Total number of eyes	9564
Number of eyes in HFA24-2	9281
Number of eyes in HFA10-2	3668
Total number of sets of paired data	47653
Number of sets of paired data in HFA24-2	45050
Number of sets of paired data in HFA10-2	22230
Number of HFA24-2 tests per eye	5.99 ± 6.72
Number of HFA10-2 tests per eye	5.40 ± 4.80
Interval between HFA24-2 tests (days)	243 ± 347
Interval between HFA10-2 tests (days)	218 ± 226
OCT to HFA24-2 interval (days)	27.3 ± 43.7
OCT to HFA10-2 interval (days)	40.1 ± 57.2
Number of OCT tests per eye	4.98 ± 5.17
Mean deviation of the HFA24-2 (dB)	-8.88 ± 8.80
Mean deviation of the HFA10-2 (dB)	-11.2 ± 9.65
Age (years)	69.0 ± 14.8
SSI	8.58 ± 1.06
Focus (D)	-1.73 ± 3.10

SD, standard deviation; HFA, Humphrey Field Analyzer; SSI, signal strength index of OCT image; OCT, optical coherence tomography.

*Values are presented as N or mean ± standard deviation as appropriate.