

Figure S1




Figure S1 - Boxplots of control miRNA used to normalise data across experiments. The values are processed sequencing reads, as described in the Methods section..

Table 1 – Summary of statistical analysis of all variables in the full dataset. Tests between groups were conducted by t-test or non-parametric Kruskal–Wallis test for continuous variables. Data are presented as mean (standard deviation, STD) or median (interquartile range) for continuous variables. Significant values are given in red.

Type	Variable	Time	Normal distribution	AKI	non-AKI	p.value
AKI biomarkers	SCr	pre	No	85 (73-104.75)	88 (73-103)	0.848
AKI biomarkers	SCr	icu	No	96 (82.25-112.75)	92 (76-108)	0.172
AKI biomarkers	SCr	6-12hr	No	113 (99-130)	96 (82-115)	0.001
AKI biomarkers	SCr	24hr	No	122 (99.75-153.5)	96 (76-113)	<0.001
AKI biomarkers	SCr	48hr	No	132 (110-160)	91 (77.5-111)	<0.001
AKI biomarkers	SCr	72hr	No	116 (96-142.25)	90.5 (78.75-109.5)	<0.001
AKI biomarkers	SCr	96hr	No	111 (96-142.5)	85.5 (76-97.25)	<0.001
EV analysis	Annexin V	pre	No	10.348 (5.516-25.045)	8.198 (4.853-14.263)	0.263
EV analysis	Annexin V	icu	No	11.521 (7.352-17.58)	11.508 (6.367-19.875)	0.771
EV analysis	Annexin V	6-12hr	No	10.765 (4.69-16.372)	10.343 (5.889-18.661)	0.504
EV analysis	Annexin V	24hr	No	12.052 (8.602-22.61)	13.144 (4.198-21.769)	0.355
EV analysis	Annexin V	48hr	No	10.809 (4.114-18.435)	12.088 (6.875-18.915)	0.489
EV analysis	CD14 (monocytes)	pre	No	2.375 (0.89-5.362)	2.05 (0.785-5.94)	0.769
EV analysis	CD14 (monocytes)	icu	No	1.805 (0.697-5.075)	3.14 (1.24-5.03)	0.242
EV analysis	CD14 (monocytes)	6-12hr	No	1.66 (0.54-5.18)	1.78 (0.968-4.718)	0.297
EV analysis	CD14 (monocytes)	24hr	No	2.11 (0.767-4.29)	3.245 (0.965-5.707)	0.242
EV analysis	CD14 (monocytes)	48hr	No	2.21 (0.812-4.865)	1.73 (0.815-4.11)	0.698
EV analysis	CD142 (tissue factor)	pre	No	7.21 (2.458-15.31)	5.77 (2.685-25.245)	0.687
EV analysis	CD142 (tissue factor)	icu	No	5.905 (1.96-13.805)	8.45 (5.21-14.97)	0.127
EV analysis	CD142 (tissue factor)	6-12hr	No	5.705 (1.97-9.82)	7.875 (4.18-12.948)	0.1
EV analysis	CD142 (tissue factor)	24hr	No	4.72 (2.16-11.93)	10.08 (3.643-19.898)	0.194
EV analysis	CD142 (tissue factor)	48hr	No	5 (2.035-15.682)	6.32 (3.16-8.8)	0.724
EV analysis	CD144 (endothelial)	pre	No	3.48 (1.785-8.56)	2.865 (1.343-4.902)	0.301
EV analysis	CD144 (endothelial)	icu	No	4.64 (2.085-12.79)	4.86 (2.66-13.42)	0.809
EV analysis	CD144 (endothelial)	6-12hr	No	3.56 (2.04-7.95)	6.285 (2.77-10.485)	0.108
EV analysis	CD144 (endothelial)	24hr	No	5.27 (2.09-9.37)	5.71 (3.308-8.422)	0.723
EV analysis	CD144 (endothelial)	48hr	No	2.91 (1.21-6.293)	3.65 (2.105-6.87)	0.319
EV analysis	CD16 (granulocytes)	pre	No	3.64 (1.878-5.648)	2.49 (1.4-4.217)	0.186
EV analysis	CD16 (granulocytes)	icu	No	4.21 (2.57-7.232)	3.76 (2.32-5.4)	0.638
EV analysis	CD16 (granulocytes)	6-12hr	No	4.77 (2.23-6.6)	4.445 (2.475-7.033)	0.702
EV analysis	CD16 (granulocytes)	24hr	No	5.55 (2.62-7.54)	4.285 (2.378-5.97)	0.361
EV analysis	CD16 (granulocytes)	48hr	No	6.245 (3.6-8.025)	4.72 (3.26-6.35)	0.132
EV analysis	CD235a (red cells)	pre	No	3.47 (1.208-10.887)	2.645 (0.945-8.25)	0.418
EV analysis	CD235a (red cells)	icu	No	2.35 (0.915-8.76)	3.71 (1.92-6.1)	0.554
EV analysis	CD235a (red cells)	6-12hr	No	1.43 (0.61-4.61)	3.115 (1.092-7.067)	0.11
EV analysis	CD235a (red cells)	24hr	No	1.94 (0.82-4.47)	1.865 (0.627-5.352)	0.819
EV analysis	CD235a (red cells)	48hr	No	1.71 (0.597-5.092)	2.05 (1.18-6.66)	0.483
EV analysis	CD284 (myeloid cells)	pre	No	14.24 (8.02-24.81)	13.525 (7.902-29.665)	0.856
EV analysis	CD284 (myeloid cells)	icu	No	15.88 (4.47-25.787)	14.51 (7.39-28.38)	0.733
EV analysis	CD284 (myeloid cells)	6-12hr	No	11.14 (7.47-17.92)	15.87 (7.01-26.64)	0.371
EV analysis	CD284 (myeloid cells)	24hr	No	13.4 (5.1-31.51)	12.865 (5.295-21.247)	0.623
EV analysis	CD284 (myeloid cells)	48hr	No	10.795 (4.192-22.282)	8.77 (4.46-15.56)	0.5
EV analysis	CD3 (T-cells)	pre	No	5.93 (3.143-13.482)	7.135 (4.33-16.29)	0.259
EV analysis	CD3 (T-cells)	icu	No	6.46 (3.267-12.688)	7.88 (6.19-12.63)	0.236
EV analysis	CD3 (T-cells)	6-12hr	No	6.3 (3.74-10.19)	7.8 (5.355-14.19)	0.055
EV analysis	CD3 (T-cells)	24hr	No	9.07 (5.32-14.613)	8.48 (5.96-11.453)	0.423
EV analysis	CD3 (T-cells)	48hr	No	9.66 (4.093-22.778)	6.88 (5.225-11.11)	0.148
EV analysis	CD41 (platelets)	pre	No	0.255 (0.078-0.635)	0.13 (0.038-0.252)	0.05
EV analysis	CD41 (platelets)	icu	No	0.155 (0.062-0.483)	0.08 (0.04-0.2)	0.128
EV analysis	CD41 (platelets)	6-12hr	No	0.09 (0.03-0.45)	0.06 (0.032-0.305)	0.769
EV analysis	CD41 (platelets)	24hr	No	0.095 (0.05-0.185)	0.05 (0.02-0.18)	0.028
EV analysis	CD41 (platelets)	48hr	No	0.065 (0.022-0.242)	0.05 (0.025-0.295)	0.986
EV analysis	CD62E (endothelial)	pre	No	2.46 (1.088-6.087)	2.425 (1.55-4.64)	0.41
EV analysis	CD62E (endothelial)	icu	No	2.635 (1.032-5.715)	3.27 (1.77-5.12)	0.617
EV analysis	CD62E (endothelial)	6-12hr	No	2.7 (0.82-4.77)	3.15 (1.885-6.738)	0.176
EV analysis	CD62E (endothelial)	24hr	No	4.085 (1.302-6.535)	3.045 (1.24-4.452)	0.182
EV analysis	CD62E (endothelial)	48hr	No	3.51 (1.765-5.535)	3.79 (1.56-4.665)	0.453
EV analysis	MV concentration	pre	No	9.197 (4.887-22.753)	14.455 (5.036-20.823)	0.916
EV analysis	MV concentration	icu	No	11.915 (6.937-15.312)	8.585 (6.281-15.496)	0.225
EV analysis	MV concentration	6-12hr	No	8.262 (2.955-17.048)	5.952 (2.714-13.603)	0.414
EV analysis	MV concentration	24hr	No	5.458 (2.999-10.675)	5.697 (2.784-10.822)	0.806
EV analysis	MV concentration	48hr	No	7.137 (3.987-15.998)	6.323 (3.405-10.257)	0.179
miRNA	hsa-miR-1262	pre	No	8.281 (0-1986.495)	1 (0-6714.537)	0.566
miRNA	hsa-miR-1262	icu	No	28.677 (0.241-26293.808)	0 (0-4.311)	0.005
miRNA	hsa-miR-1262	24hr	No	8.006 (0-2846.239)	0 (0-88.788)	0.106
miRNA	hsa-miR-133b	pre	No	0.694 (0.021-22.034)	31.481 (0.126-61660.576)	0.01
miRNA	hsa-miR-133b	icu	No	26.371 (0.72-3860.3)	351.216 (2.716-2709.37)	0.541
miRNA	hsa-miR-133b	24hr	No	165.899 (2.676-3099.212)	113.148 (0.846-2516.746)	0.767
miRNA	hsa-miR-920	pre	No	0 (0-0)	0 (0-0.002)	0.056
miRNA	hsa-miR-920	icu	No	0 (0-0)	0 (0-0)	0.463
miRNA	hsa-miR-920	24hr	No	0 (0-0)	0 (0-0)	0.813
miRNA	hsa-miR-92a-1	pre	No	0.009 (0-0.758)	0.004 (0-1)	0.733
miRNA	hsa-miR-92a-1	icu	No	0.836 (0.008-9.821)	0.007 (0-5.022)	0.201
miRNA	hsa-miR-92a-1	24hr	No	0.785 (0-8.602)	0 (0-1)	0.064
Multiplate analysis	ADP	pre	Yes	703.615 (279.448)	779.026 (262.799)	0.194
Multiplate analysis	ADP	icu	No	457 (277-696)	659.5 (385.75-1031.75)	0.064
Multiplate analysis	ADP	6-12hr	Yes	610.113 (301.135)	797.000 (292.303)	0.003
Multiplate analysis	ADP	24hr	No	505.5 (359.75-780)	559 (473-840)	0.213
Multiplate analysis	ADP	48hr	No	441 (327-574)	413 (337.5-668)	0.457
Multiplate analysis	ASPI	pre	No	493 (176-786)	311.5 (178.25-671.25)	0.312
Multiplate analysis	ASPI	icu	No	379 (145-688)	296.5 (173.75-599.25)	0.938
Multiplate analysis	ASPI	6-12hr	No	445 (200-729)	374 (275.5-704.25)	0.822
Multiplate analysis	ASPI	24hr	No	370.5 (233-597.5)	367 (249-604.5)	0.837
Multiplate analysis	ASPI	48hr	No	267 (170-443)	270 (192.5-456)	0.663
Multiplate analysis	TRAP	pre	Yes	1168.962 (354.421)	1227.921 (300.653)	0.397
Multiplate analysis	TRAP	icu	Yes	1052.849 (526.757)	1162.95 (433.772)	0.272
Multiplate analysis	TRAP	6-12hr	Yes	1211.396 (406.153)	1342.4 (316.585)	0.084

Multiplate analysis	TRAP	24hr	Yes	949.926 340.307732648069	947.816 (373.635)	0.978
Multiplate analysis	TRAP	48hr	No	750 (572-1005)	700 (505-900.5)	0.422
Flow cytometry analysis	Activated GPIIb/IIIa (PAC-1)	pre	Yes	53.766 (28.359)	57.142 (25.778)	0.605
Flow cytometry analysis	Activated GPIIb/IIIa (PAC-1)	icu	No	41.4 (22.72-68)	37.195 (24.587-62.09)	0.549
Flow cytometry analysis	Activated GPIIb/IIIa (PAC-1)	6-12hr	No	35.885 (18.65-60.427)	24.56 (13.65-44.435)	0.043
Flow cytometry analysis	Activated GPIIb/IIIa (PAC-1)	24hr	No	58.03 (32.03-80.318)	45.695 (36.888-64.553)	0.433
Flow cytometry analysis	Activated GPIIb/IIIa (PAC-1)	48hr	No	46.93 (32.525-84.398)	41.835 (31.363-66.2)	0.398
Flow cytometry analysis	CD11b	pre	No	3.83 (0.76-12.015)	2.28 (0.83-3.73)	0.201
Flow cytometry analysis	CD11b	icu	No	8.28 (1.415-27.985)	5.45 (0.472-21.317)	0.372
Flow cytometry analysis	CD11b	6-12hr	No	4.36 (0.82-18.69)	3.645 (0.855-9.88)	0.787
Flow cytometry analysis	CD11b	24hr	No	6.31 (1.84-14.69)	2.38 (0.89-10.367)	0.251
Flow cytometry analysis	CD11b	48hr	No	3.84 (1.19-21.173)	1.42 (0.66-3.94)	0.357
Flow cytometry analysis	CD14/CD41	pre	No	7.245 (5.14-8.285)	6.235 (4.784-7.755)	0.294
Flow cytometry analysis	CD14/CD41	icu	No	4.269 (2.783-5.969)	3.33 (1.035-5.895)	0.302
Flow cytometry analysis	CD14/CD41	6-12hr	No	7.875 (2.829-14.554)	6.67 (1.376-11.7)	0.301
Flow cytometry analysis	CD14/CD41	24hr	No	7.38 (3.74-13.081)	4.7 (2.455-9.385)	0.201
Flow cytometry analysis	CD14/CD41	48hr	No	7.52 (2.327-10.803)	6.013 (2.596-11.314)	0.971
Flow cytometry analysis	CD16/CD41	pre	No	37.41 (16.993-53.438)	13.304 (2.419-39.443)	0.009
Flow cytometry analysis	CD16/CD41	icu	No	19.682 (5.096-43.63)	8.182 (1.396-26.928)	0.133
Flow cytometry analysis	CD16/CD41	6-12hr	No	17.158 (6.192-46.396)	5.823 (1.204-14.831)	0.006
Flow cytometry analysis	CD16/CD41	24hr	No	29.968 (13.244-48.674)	11.357 (1.816-26.631)	0.019
Flow cytometry analysis	CD16/CD41	48hr	No	30.452 (16.14-48.169)	17.124 (3.477-38.01)	0.075
Flow cytometry analysis	CD41/CD62P	pre	No	11.14 (7.87-20.16)	11.77 (9.06-18.94)	0.583
Flow cytometry analysis	CD41/CD62P	icu	No	10.86 (8.14-18.13)	10.12 (5.51-14.675)	0.393
Flow cytometry analysis	CD41/CD62P	6-12hr	No	10.07 (7.978-14.555)	9.935 (6.522-13.608)	0.572
Flow cytometry analysis	CD41/CD62P	24hr	No	11.73 (7.387-15.41)	9.82 (7.237-13.675)	0.316
Flow cytometry analysis	CD41/CD62P	48hr	No	12.34 (8.7-21.12)	11.185 (5-16.415)	0.259
Flow cytometry analysis	CD64/CD163	pre	No	7.476 (5.805-10.05)	6.79 (5.543-8.937)	0.378
Flow cytometry analysis	CD64/CD163	icu	No	8.75 (3.891-16.539)	7.509 (1.975-12.405)	0.308
Flow cytometry analysis	CD64/CD163	6-12hr	No	14.929 (6.672-28.555)	11.061 (5.195-27.611)	0.288
Flow cytometry analysis	CD64/CD163	24hr	No	16.21 (7.499-28.39)	7.509 (5.081-16.515)	0.05
Flow cytometry analysis	CD64/CD163	48hr	No	16.335 (7.939-30.59)	8.913 (6.177-24.41)	0.266
Circulating biomarkers	CXCL1	pre	No	2.768 (2.318-10.4)	3.199 (2.318-6.65)	0.341
Circulating biomarkers	CXCL1	icu	No	12.119 (2.462-30.367)	6.679 (2.682-14.154)	0.167
Circulating biomarkers	CXCL1	6-12hr	No	6.085 (2.318-16.255)	3.989 (2.318-9.676)	0.338
Circulating biomarkers	CXCL1	24hr	No	3.898 (2.318-13.169)	2.461 (2.318-6.826)	0.32
Circulating biomarkers	CXCL1	48hr	No	3.136 (2.318-14.852)	2.61 (2.318-7.86)	0.661
Circulating biomarkers	ICAM1	pre	No	1599162.035 (888565.143-2159956.672)	1050797.044 (503857.618-1563558.9)	0.046
Circulating biomarkers	ICAM1	icu	No	1513902.02 (696906.761-1998536.287)	1141869.135 (934889.756-1717702.125)	0.261
Circulating biomarkers	ICAM1	6-12hr	No	1491300.445 (911558.2-2225791.226)	1309041.26 (911398.265-2018495.947)	0.524
Circulating biomarkers	ICAM1	24hr	No	1679358.43 (939582.854-2518320.279)	1385257.45 (787516.012-1977158.473)	0.157
Circulating biomarkers	ICAM1	48hr	No	1948995.848 (1162919.005-3104333.25)	1012597.923 (837179.132-1688613.032)	0.002
Circulating biomarkers	IL10	pre	No	0.355 (0.06-1.07)	0.208 (0.06-0.937)	0.561
Circulating biomarkers	IL10	icu	No	23.377 (8.31-36.178)	26.424 (10.154-59.832)	0.222
Circulating biomarkers	IL10	6-12hr	No	8.146 (2.975-21.964)	6.219 (2.32-16.935)	0.528
Circulating biomarkers	IL10	24hr	No	3.085 (1.259-7.27)	1.721 (0.3-6.29)	0.172
Circulating biomarkers	IL10	48hr	No	1.311 (0.193-3.342)	1.345 (0.06-3.7)	0.745
Circulating biomarkers	IL6	pre	No	0.52 (0.52-0.524)	0.52 (0.52-0.52)	0.581
Circulating biomarkers	IL6	icu	No	53.786 (23.085-129.063)	49.056 (25.734-79.706)	0.535
Circulating biomarkers	IL6	6-12hr	No	45.864 (20.607-82.747)	46.023 (23.639-60.262)	0.742
Circulating biomarkers	IL6	24hr	No	34.397 (11.97-59.674)	25.535 (13.926-45.318)	0.542
Circulating biomarkers	IL6	48hr	No	28.154 (11.699-62.629)	30.19 (14.471-74.294)	0.939
Circulating biomarkers	IL8	pre	No	0.003 (0.003-0.826)	0.003 (0.003-0.957)	0.744
Circulating biomarkers	IL8	icu	No	3.138 (1.529-7.707)	4.662 (1.67-6.969)	0.692
Circulating biomarkers	IL8	6-12hr	No	3.492 (1.421-5.645)	3.409 (1.794-5.447)	0.939
Circulating biomarkers	IL8	24hr	No	2.082 (0.655-3.617)	1.71 (0.425-3.757)	0.709
Circulating biomarkers	IL8	48hr	No	1.848 (0.22-3.56)	1.241 (0.238-2.709)	0.438

eTable 2 - Significantly different miRs - logFC - log fold change; AveExpr - average expression; t - t-statistics; P.Value - unadjusted p-value; adj.P.Val - Benjamini-Hochberd-adjusted p-value, B - B-statistics

Time	miRNA	logFC	AveExpr	t	P.Value	adj.P.Val	B
pre	has-miR-668	-6.923	-30.550	-5.852	<0.001	<0.001	6.292
ICU	hsa-miR-518a-3p	4.541	-31.097	210.544	<0.001	<0.001	114.872
ICU	hsa-miR-133b	17.189	-26.748	5.342	<0.001	0.001	4.461
ICU	hsa-miR-92a	7.172	-30.102	4.933	<0.001	0.002	3.005
ICU	Hsa-miR-1262	9.919	-8.115	4.487	<0.001	0.007	1.583
ICU	Hsa-miR-920	6.052	-30.279	3.942	<0.001	0.034	-0.210

eTable3

eTable 3 – Summary of statistical analysis of all variables in patients receiving aspirin only. Tests between groups were conducted by t-test or non-parametric Kruskal-Wallis test for continuous variables. Data are presented as mean (standard deviation, STD) or median (interquartile range) for continuous variables. Significant values are given in red.

Type	Variable	Time	Normal	AKI	non-AKI	p-value
AKI biomarkers	Serum creatinine	pre	Yes	86.188 (19.42)	86.793 (22.823)	0.912
AKI biomarkers	Serum creatinine	icu	No	93 (80 - 103.75)	88 (76 - 101)	0.374
AKI biomarkers	Serum creatinine	6hr	No	106 (96 - 123.5)	94.5 (82 - 113.75)	0.016
AKI biomarkers	Serum creatinine	24hr	No	114 (97 - 141)	100.5 (77.5 - 113)	0.003
AKI biomarkers	Serum creatinine	48hr	No	128.5 (106.25 - 142.5)	91 (83 - 111)	0.001
AKI biomarkers	Serum creatinine	72hr	No	114 (91.5 - 133.5)	91 (78 - 109)	0.018
AKI biomarkers	Serum creatinine	96hr	No	105.5 (89.5 - 131)	85.5 (77.25 - 97.75)	0.011
EV analysis	Annexin V	pre	No	11.828 (6.252 - 30.099)	9.666 (5.183 - 14.799)	0.201
EV analysis	Annexin V	icu	No	11.651 (8.997 - 19.112)	13.585 (6.044 - 22.033)	0.993
EV analysis	Annexin V	6hr	No	11.135 (5.824 - 16.501)	11.255 (6.344 - 19.241)	0.636
EV analysis	Annexin V	24hr	No	15.241 (9.366 - 24.579)	11.471 (4.198 - 26.488)	0.216
EV analysis	Annexin V	48hr	No	9.477 (3.473 - 18.407)	12.088 (9.572 - 18.727)	0.279
EV analysis	CD14 (monocytes)	pre	No	3.98 (1.418 - 6.2)	2.34 (0.785 - 7.98)	0.946
EV analysis	CD14 (monocytes)	icu	No	1.92 (0.922 - 5.432)	3.06 (1.255 - 6.365)	0.506
EV analysis	CD14 (monocytes)	6hr	No	1.76 (0.66 - 5.275)	2.215 (1.102 - 5.228)	0.329
EV analysis	CD14 (monocytes)	24hr	No	2.335 (0.845 - 5.44)	4.135 (1.465 - 5.905)	0.325
EV analysis	CD14 (monocytes)	48hr	No	2.65 (1.225 - 4.46)	2.16 (1.29 - 5.79)	0.941
EV analysis	CD142 (tissue factor)	pre	No	10.49 (3.896 - 19.472)	12.27 (2.685 - 27.69)	0.759
EV analysis	CD142 (tissue factor)	icu	No	8 (2.287 - 17.94)	7.86 (4.828 - 14.15)	0.755
EV analysis	CD142 (tissue factor)	6hr	No	8.24 (2.185 - 10.78)	7.875 (4.3 - 12.583)	0.337
EV analysis	CD142 (tissue factor)	24hr	No	6.28 (2.435 - 16.005)	10.09 (4.715 - 20.123)	0.429
EV analysis	CD142 (tissue factor)	48hr	No	7.06 (2.93 - 19.935)	6.58 (3.9 - 8.95)	0.612
EV analysis	CD144 (endothelial)	pre	No	6.61 (2.645 - 11.08)	2.975 (1.343 - 5.803)	0.148
EV analysis	CD144 (endothelial)	icu	No	5.62 (2.783 - 13.555)	4.245 (2.135 - 14.55)	0.736
EV analysis	CD144 (endothelial)	6hr	No	3.79 (2.325 - 9.755)	6.94 (3.565 - 10.622)	0.255
EV analysis	CD144 (endothelial)	24hr	No	5.26 (2.208 - 10.375)	5.98 (3.433 - 8.745)	0.561
EV analysis	CD144 (endothelial)	48hr	No	2.77 (1.585 - 6.18)	3.7 (2.22 - 7.22)	0.226
EV analysis	CD16 (granulocytes)	pre	No	3.64 (2.056 - 5.663)	2.08 (1.298 - 3.848)	0.108
EV analysis	CD16 (granulocytes)	icu	No	4.41 (2.587 - 6.207)	3.005 (2.245 - 4.852)	0.073
EV analysis	CD16 (granulocytes)	6hr	No	4.92 (2.265 - 6.625)	4.275 (2.725 - 6.365)	0.894
EV analysis	CD16 (granulocytes)	24hr	No	5.785 (3.788 - 7.883)	4.375 (2.378 - 6.62)	0.329
EV analysis	CD16 (granulocytes)	48hr	No	6.75 (4.165 - 9.355)	4.77 (2.45 - 6.64)	0.06
EV analysis	CD235a (red cells)	pre	No	5.295 (1.388 - 15.085)	2.8 (1.078 - 9.93)	0.159
EV analysis	CD235a (red cells)	icu	No	3.99 (1.85 - 9.41)	3.785 (1.622 - 5.897)	0.742
EV analysis	CD235a (red cells)	6hr	No	2.11 (0.74 - 5.355)	2.79 (1.117 - 6.595)	0.499
EV analysis	CD235a (red cells)	24hr	No	2.345 (0.893 - 5.212)	2.375 (0.977 - 5.457)	0.994
EV analysis	CD235a (red cells)	48hr	No	1.38 (0.355 - 4.065)	2.68 (1.43 - 6.95)	0.232
EV analysis	CD284 (myeloid cells)	pre	No	19.35 (10.46 - 25.325)	15.05 (9.305 - 29.665)	0.643
EV analysis	CD284 (myeloid cells)	icu	No	17.32 (8.015 - 33.469)	15.37 (8.165 - 24.322)	0.751
EV analysis	CD284 (myeloid cells)	6hr	No	12.15 (5.035 - 19.92)	16.135 (7.94 - 26.913)	0.432
EV analysis	CD284 (myeloid cells)	24hr	No	15.81 (4.845 - 35.127)	11.005 (5.055 - 18.722)	0.263
EV analysis	CD284 (myeloid cells)	48hr	No	10.7 (1.81 - 19.53)	8.77 (5.05 - 15.43)	0.856
EV analysis	CD3 (T-cells)	pre	No	6.05 (3.235 - 13.698)	9.505 (4.883 - 16.29)	0.213
EV analysis	CD3 (T-cells)	icu	No	6.785 (3.845 - 11.373)	7.86 (6.568 - 10.957)	0.509
EV analysis	CD3 (T-cells)	6hr	No	7.29 (4.065 - 10.495)	8.745 (5.738 - 13.37)	0.216
EV analysis	CD3 (T-cells)	24hr	No	8.93 (5.92 - 14.33)	8.57 (6.215 - 11.265)	0.614
EV analysis	CD3 (T-cells)	48hr	No	8.88 (4.9 - 15.445)	7.59 (5.46 - 10.73)	0.342
EV analysis	CD41 (platelets)	pre	No	0.355 (0.1 - 1.235)	0.135 (0.045 - 0.37)	0.023
EV analysis	CD41 (platelets)	icu	No	0.165 (0.07 - 0.723)	0.06 (0.04 - 0.19)	0.096
EV analysis	CD41 (platelets)	6hr	No	0.11 (0.025 - 0.5)	0.06 (0.04 - 0.328)	0.643
EV analysis	CD41 (platelets)	24hr	No	0.11 (0.062 - 0.238)	0.05 (0.03 - 0.115)	0.032
EV analysis	CD41 (platelets)	48hr	No	0.06 (0.025 - 0.185)	0.05 (0.03 - 0.28)	0.791
EV analysis	CD62E (endothelial)	pre	No	2.22 (1.04 - 6.305)	2.72 (1.617 - 5.315)	0.278
EV analysis	CD62E (endothelial)	icu	No	3.13 (1.515 - 6.555)	4.27 (2.077 - 5.3)	0.813
EV analysis	CD62E (endothelial)	6hr	No	3.27 (1.44 - 6.64)	4.26 (2.277 - 7.778)	0.402
EV analysis	CD62E (endothelial)	24hr	No	4 (1.45 - 5.96)	4.15 (2.13 - 4.985)	0.743
EV analysis	CD62E (endothelial)	48hr	No	3.65 (2.065 - 4.61)	4.05 (2.52 - 4.84)	0.76
EV analysis	EV concentration	pre	No	9.707 (4.696 - 21.838)	15.02 (4.863 - 24.41)	0.866
EV analysis	EV concentration	icu	No	10.533 (7.043 - 12.528)	8.937 (6.355 - 15.465)	0.719
EV analysis	EV concentration	6hr	No	7.875 (2.694 - 16.305)	6.1 (2.732 - 13.237)	0.711
EV analysis	EV concentration	24hr	No	5.677 (3.363 - 11.437)	5.11 (3.367 - 10.08)	0.486
EV analysis	EV concentration	48hr	No	6.647 (3.498 - 15.153)	6.17 (3.04 - 8.16)	0.269
miRNA	hsa-miR-1262	pre	No	27.328 (0-1700.292)	0.509 (0-4883.161)	0.677
miRNA	hsa-miR-1262	icu	No	199.713 (14.878-167301.166)	0.001 (0-4.311)	0.016
miRNA	hsa-miR-1262	24hr	No	4.294 (0.001-3327.722)	0 (0-299.063)	0.159
miRNA	hsa-miR-133b	pre	No	0.694 (0.005-19.664)	16.693 (0.063-25155.195)	0.107
miRNA	hsa-miR-133b	icu	No	34.88 (1.629-26185.531)	423.669 (9.207-3857.218)	0.647
miRNA	hsa-miR-133b	24hr	No	171.48 (6.089-2755.387)	114.52 (2.97-1387.016)	0.567
miRNA	hsa-miR-920	pre	No	0 (0-0)	0 (0-0.001)	0.233
miRNA	hsa-miR-920	icu	No	0 (0-0)	0 (0-0)	0.204
miRNA	hsa-miR-920	24hr	No	0 (0-0)	0 (0-0)	0.524
miRNA	hsa-miR-92a-1	pre	No	0 (0-0.222)	0.004 (0-0.052)	0.596
miRNA	hsa-miR-92a-1	icu	No	0.792 (0.024-22.092)	0.01 (0-149.256)	0.667
miRNA	hsa-miR-92a-1	24hr	No	1.136 (0-20.842)	0 (0-11.398)	0.173
Multiplate analysis	ADP	pre	Yes	665.25 (237.21)	753.222 (283.415)	0.207

eTable3

Multiplate analysis	ADP	icu	No	427 (223.5 - 658)	548 (361 - 784)	0.125
Multiplate analysis	ADP	6hr	No	552 (475.5 - 867.5)	817 (561 - 1002)	0.025
Multiplate analysis	ADP	24hr	Yes	549.829 (222.828)	590.346 (265.239)	0.531
Multiplate analysis	ADP	48hr	No	362.5 (301.25 - 516)	399 (335.25 - 602.5)	0.336
Multiplate analysis	ASPI	pre	No	323 (165.5 - 772.5)	286 (114.5 - 465.5)	0.249
Multiplate analysis	ASPI	icu	No	251 (149 - 619.5)	232 (137 - 493)	0.599
Multiplate analysis	ASPI	6hr	No	445 (198.5 - 623.5)	365 (265 - 702)	0.826
Multiplate analysis	ASPI	24hr	No	310 (234 - 627)	342 (265.25 - 455.75)	0.983
Multiplate analysis	ASPI	48hr	No	223.5 (146.75 - 348.5)	225 (185.25 - 381.5)	0.475
Multiplate analysis	TRAP	pre	Yes	1158.656 (350.609)	1235.407 (324.391)	0.387
Multiplate analysis	TRAP	icu	Yes	1060.548 (535.746)	1115.552 (423.282)	0.66
Multiplate analysis	TRAP	6hr	Yes	1247.692 (392.467)	1331.424 (302.441)	0.311
Multiplate analysis	TRAP	24hr	Yes	952.829 (273.392)	842.2 (362.364)	0.205
Multiplate analysis	TRAP	48hr	No	673.5 (562.5 - 792.25)	640 (504 - 846.25)	0.628
Flow cytometry analysis	Activated GPIIb/IIIa (PAC-1)	pre	Yes	52.078 (29.695)	56.848 (26.03)	0.565
Flow cytometry analysis	Activated GPIIb/IIIa (PAC-1)	icu	No	45.84 (19.698 - 73.232)	33.74 (16.72 - 51.78)	0.293
Flow cytometry analysis	Activated GPIIb/IIIa (PAC-1)	6hr	No	48.04 (18.43 - 60.562)	26.78 (17.433 - 44.772)	0.104
Flow cytometry analysis	Activated GPIIb/IIIa (PAC-1)	24hr	No	50.28 (29.513 - 80.073)	47.01 (34.18 - 56.56)	0.71
Flow cytometry analysis	Activated GPIIb/IIIa (PAC-1)	48hr	No	40.525 (25.315 - 85.828)	43.965 (31.363 - 66.2)	1
Flow cytometry analysis	CD11b	pre	No	3.715 (0.603 - 17.413)	2.39 (0.6 - 4.685)	0.582
Flow cytometry analysis	CD11b	icu	No	11.485 (0.34 - 29.6)	9.46 (0.908 - 22.192)	0.774
Flow cytometry analysis	CD11b	6hr	No	1.76 (0.63 - 19.055)	4.775 (0.727 - 13.085)	0.864
Flow cytometry analysis	CD11b	24hr	No	6.31 (1.31 - 9.735)	2.51 (1.145 - 12.055)	0.808
Flow cytometry analysis	CD11b	48hr	No	3.8 (0.865 - 13.372)	2.335 (0.672 - 3.99)	0.799
Flow cytometry analysis	CD14/CD41	pre	No	7.245 (6.036 - 9.147)	6.35 (4.968 - 7.68)	0.18
Flow cytometry analysis	CD14/CD41	icu	No	4.269 (2.593 - 5.883)	3.35 (1.15 - 5.853)	0.434
Flow cytometry analysis	CD14/CD41	6hr	No	9.56 (2.703 - 15.674)	5.5 (1.362 - 11.865)	0.294
Flow cytometry analysis	CD14/CD41	24hr	No	7.764 (4.47 - 14.972)	5.717 (2.191 - 12.03)	0.32
Flow cytometry analysis	CD14/CD41	48hr	No	7.852 (2.27 - 9.811)	6.265 (2.995 - 15.774)	0.727
Flow cytometry analysis	CD16/CD41	pre	No	42.801 (20.605 - 52.947)	11.936 (2.44 - 39.315)	0.003
Flow cytometry analysis	CD16/CD41	icu	No	21.366 (4.58 - 39.51)	7.872 (1.689 - 26.794)	0.227
Flow cytometry analysis	CD16/CD41	6hr	No	19.569 (8.271 - 48.52)	4.781 (1.158 - 13.887)	0.014
Flow cytometry analysis	CD16/CD41	24hr	No	27.988 (10.389 - 38.298)	12.24 (2.051 - 30.151)	0.155
Flow cytometry analysis	CD16/CD41	48hr	No	30.452 (16.12 - 47.194)	23.408 (3.642 - 46.02)	0.357
Flow cytometry analysis	CD41/CD62P	pre	No	11.77 (8.123 - 24.302)	13.365 (9.012 - 18.95)	0.801
Flow cytometry analysis	CD41/CD62P	icu	No	11.46 (8.515 - 18.23)	10.02 (4.615 - 13.89)	0.17
Flow cytometry analysis	CD41/CD62P	6hr	No	10.925 (8.512 - 16.135)	10.62 (7.01 - 13.355)	0.504
Flow cytometry analysis	CD41/CD62P	24hr	No	12.34 (8.54 - 14.73)	10.11 (7.38 - 14.27)	0.741
Flow cytometry analysis	CD41/CD62P	48hr	No	12.34 (8.765 - 18.7)	10.23 (4.69 - 16.22)	0.316
Flow cytometry analysis	CD64/CD163	pre	No	7.238 (5.795 - 10.135)	6.77 (5.684 - 8.748)	0.493
Flow cytometry analysis	CD64/CD163	icu	No	5.359 (3.463 - 15.035)	7.614 (1.892 - 12.404)	0.745
Flow cytometry analysis	CD64/CD163	6hr	No	14.89 (7.224 - 34.879)	10.085 (4.702 - 20.658)	0.218
Flow cytometry analysis	CD64/CD163	24hr	No	15.659 (6.395 - 24.184)	8.527 (4.57 - 16.269)	0.238
Flow cytometry analysis	CD64/CD163	48hr	No	17.861 (7.676 - 29.341)	9.416 (6.618 - 21.383)	0.53
Circulating biomarkers	CXCL1	pre	No	3.744 (2.318 - 10.585)	3.199 (2.318 - 6.647)	0.269
Circulating biomarkers	CXCL1	icu	No	11.112 (2.615 - 34.176)	6.712 (3.332 - 13.968)	0.244
Circulating biomarkers	CXCL1	6hr	No	6.567 (2.318 - 15.784)	4.911 (2.318 - 9.637)	0.313
Circulating biomarkers	CXCL1	24hr	No	5.727 (2.318 - 15.051)	2.318 (2.318 - 6.65)	0.266
Circulating biomarkers	CXCL1	48hr	No	2.768 (2.318 - 17.557)	2.349 (2.318 - 6.65)	0.475
Circulating biomarkers	ICAM1	pre	No	1308557.275 (839981.001 - 1732249.1461197292.905)	543823.482 - 1532129.496	0.269
Circulating biomarkers	ICAM1	icu	No	1497879.046 (647305.841 - 2084625.5611095508.272)	806047.016 - 1643080.197	0.425
Circulating biomarkers	ICAM1	6hr	No	1298714.983 (829995.511 - 2210903.8461304970.189)	813370.504 - 2149794.037	0.921
Circulating biomarkers	ICAM1	24hr	No	1654346.776 (779856.163 - 2549254.1861465262.722)	852929.49 - 2208904.562	0.845
Circulating biomarkers	ICAM1	48hr	No	1750900.171 (1052195.998 - 3057463.001115578.167)	784847.679 - 1668053.479	0.035
Circulating biomarkers	IL10	pre	No	0.06 (0.06 - 0.922)	0.481 (0.06 - 1.074)	0.553
Circulating biomarkers	IL10	icu	No	23.207 (9.614 - 35.575)	36.22 (10.347 - 58.5)	0.157
Circulating biomarkers	IL10	6hr	No	8.049 (3.017 - 21.394)	6.421 (2.395 - 17.51)	0.803
Circulating biomarkers	IL10	24hr	No	5.028 (1.53 - 7.69)	2.365 (0.472 - 6.31)	0.197
Circulating biomarkers	IL10	48hr	No	1.548 (0.528 - 3.195)	1.026 (0.06 - 2.627)	0.319
Circulating biomarkers	IL6	pre	No	0.52 (0.52 - 0.521)	0.52 (0.52 - 0.52)	0.712
Circulating biomarkers	IL6	icu	No	52.094 (21.971 - 129.317)	43.429 (26.151 - 85.071)	0.681
Circulating biomarkers	IL6	6hr	No	46.908 (19.698 - 85.864)	52.963 (21.065 - 60.958)	0.825
Circulating biomarkers	IL6	24hr	No	46.319 (12.592 - 74.567)	27.611 (18.354 - 45.902)	0.38
Circulating biomarkers	IL6	48hr	No	29.903 (16.168 - 64.577)	30.337 (14.471 - 74.294)	0.779
Circulating biomarkers	IL8	pre	No	0.003 (0.003 - 0.386)	0.003 (0.003 - 0.807)	0.662
Circulating biomarkers	IL8	icu	No	2.839 (1.252 - 8.194)	5.291 (1.19 - 7.24)	0.656
Circulating biomarkers	IL8	6hr	No	3.552 (1.397 - 5.67)	3.79 (2.095 - 6.285)	0.693
Circulating biomarkers	IL8	24hr	No	2.003 (0.442 - 3.646)	1.83 (0.416 - 3.449)	0.797
Circulating biomarkers	IL8	48hr	No	1.326 (0.112 - 3.035)	1.326 (0.386 - 2.265)	0.628