

Extended Data Table 1 | Internal cohort description

		Total	Training	Validation	Test
Patients, <i>n</i>		26,067	19,868	1,611	4,588
Studies, <i>n</i>		33,927	26,655	2,142	5,130
Videos, <i>n</i>		1,230,490	1,033,146	83,064	114,280
Age (years), <i>mean (sd)</i>		67.1 (16.0)	67.5 (15.9)	67.2 (15.4)	65.4 (16.7)
Sex, <i>n (%)</i>	<i>Male</i>	16,826 (52.2)	13,201 (52.5)	1,047 (52.3)	2,578 (50.4)
	<i>Female</i>	15,438 (47.8)	11,949 (47.5)	955 (47.7)	2,534 (49.6)
Race, <i>n (%)</i>	<i>White</i>	24,603 (79.7)	19,204 (79.6)	1,533 (80.3)	3,866 (80.1)
	<i>Black</i>	4,384 (14.2)	3,453 (14.3)	253 (13.3)	678 (14.0)
	<i>Asian</i>	541 (1.8)	399 (1.7)	33 (1.7)	109 (2.3)
	<i>Other</i>	1,330 (4.3)	1065 (4.4)	89 (4.7)	176 (3.6)
Aortic Root Dimension (cm), <i>mean (sd)</i>		3.33 (0.47)	3.32 (0.47)	3.32 (0.48)	3.35 (0.48)
AV Peak Velocity (m/s), <i>mean (sd)</i>		1.83 (0.89)	1.87 (0.92)	1.82 (0.89)	1.65 (0.67)
AV Regurgitation, <i>n (%)</i>	<i>None-Trace</i>	22,014 (75.9)	16975 (75.5)	1363 (75.3)	3676 (77.9)
	<i>Mild</i>	5,082 (17.5)	3980 (17.7)	329 (18.2)	773 (16.4)
	<i>Moderate or greater</i>	1,917 (6.6)	1526 (6.8)	119 (6.6)	272 (5.8)
AV Stenosis, <i>n (%)</i>	<i>None</i>	22,754 (81.3)	17118 (79.5)	1422 (81.8)	4214 (89.6)
	<i>Mild-Moderate</i>	3,209 (11.5)	2616 (12.1)	192 (11.0)	401 (8.5)
	<i>Severe</i>	2,020 (7.2)	1809 (8.4)	125 (7.2)	86 (1.8)
AV Structure, <i>n (%)</i>	<i>Bicuspid</i>	264 (1.0)	207 (1.0)	5 (0.3)	52 (1.1)
E e' Ratio, <i>mean (sd)</i>		11.42 (5.05)	11.70 (5.17)	11.38 (5.04)	10.19 (4.31)
Ejection Fraction (%), <i>mean (sd)</i>		58.98 (11.72)	58.95 (11.95)	59.23 (11.63)	59.03 (10.62)
Global Longitudinal Strain (%), <i>mean (sd)</i>		-17.26 (4.15)	-16.79 (4.30)	-16.73 (4.06)	-18.91 (3.15)
Increased LA Size, <i>n (%)</i>	<i>None</i>	20,373 (63.1)	15595 (62.1)	1251 (62.4)	3527 (68.8)
	<i>Mild</i>	5,702 (17.7)	4543 (18.1)	354 (17.7)	805 (15.7)
	<i>Moderate or greater</i>	6,190 (19.2)	4992 (19.9)	400 (20.0)	798 (15.6)
Increased LV Size, <i>n (%)</i>	<i>None</i>	21,238 (86.3)	16249 (86.2)	1329 (88.4)	3660 (85.8)
	<i>Mild</i>	1,827 (7.4)	1396 (7.4)	93 (6.2)	338 (7.9)
	<i>Moderate or greater</i>	1,553 (6.3)	1201 (6.4)	82 (5.5)	270 (6.3)
Increased LV Wall Thickness,	<i>Any</i>	15,470 (52.4)	12509 (54.8)	990 (54.0)	1971 (40.2)

<i>n (%)</i>	<i>Moderate or greater</i>	3,575 (12.7)	3021 (13.9)	239 (13.7)	315 (6.5)
Increased RA Size, <i>n (%)</i>	<i>Any</i>	3,726 (29.7)	2988 (29.8)	231 (28.2)	507 (29.3)
Increased RV Size, <i>n (%)</i>	<i>None</i>	24,868 (84.2)	19229 (83.8)	1551 (85.3)	4088 (85.9)
	<i>Mild</i>	2,940 (10.0)	2313 (10.1)	170 (9.3)	457 (9.6)
	<i>Moderate or greater</i>	1,730 (5.9)	1417 (6.2)	98 (5.4)	215 (4.5)
IVSd (cm), <i>mean (sd)</i>		1.11 (0.30)	1.12 (0.30)	1.13 (0.34)	1.02 (0.25)
LA Volume (cm³), <i>mean (sd)</i>		59.41 (19.25)	59.87 (19.37)	59.73 (19.22)	57.29 (18.58)
LAIDs (cm), <i>mean (sd)</i>		4.05 (0.82)	4.06 (0.82)	4.01 (0.81)	4.02 (0.81)
LV Diastolic Dysfunction, <i>n (%)</i>	<i>None</i>	11,809 (50.9)	8583 (48.4)	702 (50.1)	2524 (62.0)
	<i>Mild</i>	8,086 (34.9)	6376 (36.0)	496 (35.4)	1214 (29.8)
	<i>Moderate or greater</i>	3,298 (14.2)	2762 (15.6)	203 (14.5)	333 (8.2)
LV ED Volume (cm³), <i>mean (sd)</i>		107.07 (47.76)	106.07 (47.86)	104.96 (46.06)	113.51 (47.45)
LV ES Volume (cm³), <i>mean (sd)</i>		45.77 (33.25)	45.42 (33.39)	44.22 (32.06)	48.42 (32.85)
LV Stroke Volume (cm³), <i>mean (sd)</i>		61.21 (23.55)	60.53 (23.54)	60.69 (23.23)	65.21 (23.34)
LV Systolic Dysfunction, <i>n (%)</i>	<i>None-Hyperdynamic</i>	25,702 (85.1)	19820 (84.5)	1607 (85.9)	4275 (87.6)
	<i>Mild</i>	2,225 (7.4)	1782 (7.6)	124 (6.6)	319 (6.5)
	<i>Moderate or greater</i>	2,273 (7.5)	1849 (7.9)	140 (7.5)	284 (5.8)
LV Wall Motion Abnormalities, <i>n (%)</i>		2195 (20.9)	1704 (20.6)	143 (21.4)	348 (22.6)
LVIDd (cm), <i>mean (sd)</i>		4.58 (0.72)	4.58 (0.73)	4.55 (0.73)	4.64 (0.69)
LVIDs (cm), <i>mean (sd)</i>		3.07 (0.75)	3.07 (0.76)	3.04 (0.76)	3.08 (0.70)
LVOT Diameter (cm), <i>mean (sd)</i>		2.04 (0.21)	2.04 (0.21)	2.05 (0.22)	2.05 (0.21)
Elevated LVOT PG, <i>n (%)</i>	<i>≥20 mmHg</i>	298 (1.1)	249 (1.2)	25 (1.5)	24 (0.6)
LVPWd (cm), <i>mean (sd)</i>		1.05 (0.22)	1.06 (0.22)	1.06 (0.22)	1.00 (0.20)
MV Regurgitation, <i>n (%)</i>	<i>None-Trace</i>	16013 (52.3)	12041 (50.7)	996 (52.3)	2976 (59.4)
	<i>Mild</i>	10013 (32.7)	7968 (33.6)	626 (32.9)	1419 (28.3)
	<i>Moderate or greater</i>	4612 (15.1)	3718 (15.7)	282 (14.8)	612 (12.2)
MV Stenosis, <i>n (%)</i>	<i>Any</i>	901 (3.8)	758 (4.1)	54 (3.7)	89 (2.1)
Pericardial Effusion, <i>n (%)</i>	<i>Any</i>	1027 (3.9)	935 (4.0)	71 (3.8)	21 (3.4)
Elevated RA Pressure, <i>n (%)</i>	<i>≥8 mmHg</i>	5662 (20.0)	722 (15.5)	4596 (21.0)	344 (19.5)
RA Transverse Dimension (cm),		3.83 (0.78)	3.83 (0.80)	3.83 (0.77)	3.80 (0.75)

<i>mean (sd)</i>					
RV S' (cm/s), mean (sd)		12.68 (3.29)	12.60 (3.30)	12.85 (3.35)	12.93 (3.21)
RV Systolic Dysfunction, n (%)		3156 (10.7)	2589 (11.3)	196 (10.8)	371 (7.8)
RV Systolic Pressure, mean (sd)		31.73 (13.00)	32.16 (13.30)	31.10 (13.12)	29.79 (11.05)
RVIDd (cm), mean (sd)		3.15 (0.64)	3.16 (0.64)	3.13 (0.63)	3.11 (0.61)
TAPSE (cm), mean (sd)		2.14 (0.54)	2.13 (0.54)	2.16 (0.53)	2.16 (0.52)
TV Peak Gradient (mmHg), mean (sd)		26.96 (11.42)	27.28 (11.67)	26.38 (11.63)	25.48 (9.76)
TV Regurgitation, n (%)	<i>None-Trace</i>	16191 (51.9)	12301 (50.9)	1004 (52.1)	2886 (56.9)
	<i>Mild</i>	10409 (33.4)	8199 (33.9)	630 (32.7)	1580 (31.1)
	<i>Moderate or greater</i>	4570 (14.7)	3666 (15.2)	294 (15.2)	610 (12.0)

Demographics and summary statistics of the internal dataset used for model development and evaluation. Data consists of echocardiogram videos acquired at YNHHS hospitals from 2016-2022. All summary statistics for demographics and labels are computed after excluding missing values and provided at the study level since age and labels are unique to each echocardiography study. AV = aortic valve; ED = end-diastolic; ES = end-systolic; IVSd = intraventricular septum thickness at diastole; LA = left atrium; LV = left ventricle; LAIDs = left atrial internal diameter at systole; LVIDd = left ventricular internal diameter at diastole; LVIDs = left ventricular internal diameter at systole; LVOT = left ventricular outflow tract; LVPWd = left ventricular posterior wall thickness at diastole; MAE = mean absolute error; PG = pressure gradient; RA = right atrium; RV = right ventricle; RVIDd = right ventricular internal diameter at diastole; RV S' = right ventricular systolic excursion velocity; sd = standard deviation; TV = tricuspid valve; YNHHS = Yale-New Haven Health System.

Extended Data Table 2 | Internal model performance on classification tasks

		AUC	AP
Pericardial Effusion	<i>Any</i>	0.914 (0.834, 0.974)	0.992 (0.980, 0.999)
Increased LV Size	<i>None</i>	0.907 (0.897, 0.918)	0.981 (0.978, 0.984)
	<i>Mild</i>	0.833 (0.816, 0.850)	0.298 (0.262, 0.341)
	<i>Moderate or greater</i>	0.945 (0.931, 0.953)	0.685 (0.604, 0.696)
Increased LV Wall Thickness	<i>Any</i>	0.847 (0.838, 0.856)	0.890 (0.881, 0.898)
	<i>Moderate-Severe</i>	0.912 (0.898, 0.926)	0.992 (0.991, 0.994)
LV Systolic Dysfunction	<i>None-Hyperdynamic</i>	0.964 (0.961, 0.972)	0.994 (0.993, 0.996)
	<i>Mild</i>	0.922 (0.910, 0.932)	0.407 (0.365, 0.454)
	<i>Moderate or greater</i>	0.984 (0.978, 0.987)	0.858 (0.824, 0.877)
LV Wall Motion Abnormalities	<i>Any</i>	0.877 (0.858, 0.895)	0.745 (0.709, 0.781)
LV Diastolic Dysfunction	<i>None</i>	0.846 (0.832, 0.853)	0.909 (0.894, 0.910)
	<i>Mild-Indeterminate</i>	0.802 (0.790, 0.813)	0.630 (0.605, 0.655)
	<i>Moderate or greater</i>	0.926 (0.910, 0.930)	0.532 (0.479, 0.573)
Increased RV Size	<i>None</i>	0.841 (0.821, 0.848)	0.967 (0.962, 0.970)
	<i>Mild</i>	0.776 (0.759, 0.793)	0.242 (0.217, 0.274)
	<i>Moderate or greater</i>	0.873 (0.864, 0.899)	0.336 (0.285, 0.397)
RV Systolic Dysfunction	<i>Any</i>	0.926 (0.915, 0.937)	0.993 (0.991, 0.994)
Increased LA Size	<i>None</i>	0.841 (0.828, 0.847)	0.921 (0.909, 0.923)
	<i>Mild</i>	0.753 (0.739, 0.767)	0.357 (0.328, 0.385)
	<i>Moderate or greater</i>	0.889 (0.879, 0.898)	0.593 (0.562, 0.627)
Increased RA Size	<i>Any</i>	0.772 (0.752, 0.792)	0.879 (0.862, 0.895)
AV Structure	<i>Bicuspid</i>	0.806 (0.749, 0.859)	0.997 (0.995, 0.998)
AV Stenosis	<i>None</i>	0.955 (0.944, 0.959)	0.995 (0.992, 0.995)
	<i>Mild-Moderate</i>	0.933 (0.924, 0.943)	0.540 (0.499, 0.588)
	<i>Severe</i>	0.986 (0.978, 0.989)	0.644 (0.484, 0.665)
AV Regurgitation	<i>None-Trace</i>	0.863 (0.848, 0.869)	0.953 (0.947, 0.957)
	<i>Mild</i>	0.790 (0.777, 0.803)	0.396 (0.370, 0.426)
	<i>Moderate or greater</i>	0.925 (0.908, 0.933)	0.527 (0.470, 0.569)
Elevated LVOT PG	≥ 20 mmHg	0.937 (0.878, 0.990)	0.492 (0.328, 0.690)
MV Stenosis	<i>Any</i>	0.963 (0.944, 0.978)	0.999 (0.998, 0.999)
MV Regurgitation	<i>None-Trace</i>	0.854 (0.835, 0.854)	0.886 (0.871, 0.890)
	<i>Mild</i>	0.733 (0.721, 0.745)	0.489 (0.466, 0.513)
	<i>Moderate or greater</i>	0.905 (0.896, 0.915)	0.604 (0.567, 0.635)
TV Regurgitation	<i>None-Trace</i>	0.828 (0.819, 0.838)	0.854 (0.842, 0.864)
	<i>Mild</i>	0.723 (0.709, 0.734)	0.503 (0.480, 0.524)
	<i>Moderate or greater</i>	0.890 (0.887, 0.907)	0.613 (0.581, 0.647)
Elevated RA Pressure	≥ 8 mmHg	0.814 (0.801, 0.829)	0.526 (0.496, 0.560)

Multi-task classification performance on the internal YNHHS test set. Values in parentheses represent bootstrapped 95% confidence intervals. AP = average precision; AUC = area under the receiver operating characteristic curve; AV = aortic valve; LA = left atrium; LV = left ventricle; LVOT PG = left ventricular outflow tract pressure gradient; RA = right atrium; RV = right ventricle; TV = tricuspid valve; YNHHS = Yale-New Haven Health System.

Extended Data Table 3 | Internal model performance on regression tasks

	MAE	R²
Ejection Fraction (%)	4.440 (4.336, 4.551)	0.654 (0.636, 0.672)
Global Longitudinal Strain (%)	1.893 (1.780, 2.012)	0.399 (0.312, 0.473)
LV ED Volume (cm³)	18.969 (18.475, 19.461)	0.700 (0.687, 0.713)
LV ES Volume (cm³)	10.408 (10.063, 10.761)	0.761 (0.745, 0.777)
LV Stroke Volume (cm³)	11.696 (11.429, 11.981)	0.569 (0.552, 0.585)
IVSd (cm)	0.130 (0.127, 0.133)	0.478 (0.437, 0.511)
LVPWd (cm)	0.115 (0.113, 0.117)	0.449 (0.433, 0.465)
LVIDs (cm)	0.338 (0.332, 0.345)	0.596 (0.579, 0.611)
LVIDd (cm)	0.378 (0.370, 0.385)	0.506 (0.492, 0.520)
LVOT Diameter (cm)	0.135 (0.132, 0.138)	0.322 (0.307, 0.337)
E e' Ratio	1.966 (1.922, 2.009)	0.628 (0.613, 0.644)
RV Systolic Pressure (mmHg)	5.985 (5.843, 6.144)	0.461 (0.438, 0.483)
RVIDd (cm)	0.397 (0.387, 0.406)	0.300 (0.280, 0.318)
TAPSE (cm)	0.336 (0.329, 0.343)	0.333 (0.315, 0.350)
RV S' (cm/s)	1.911 (1.859, 1.960)	0.343 (0.319, 0.365)
LAIIDs (cm)	0.398 (0.390, 0.405)	0.605 (0.592, 0.616)
LA Volume (cm³)	9.412 (9.218, 9.598)	0.593 (0.577, 0.606)
RA Transverse Dimension (cm)	0.467 (0.457, 0.478)	0.368 (0.344, 0.391)
AV Peak Velocity (m/s)	0.310 (0.303, 0.318)	0.593 (0.578, 0.607)
TV Peak Gradient (mmHg)	5.572 (5.439, 5.720)	0.401 (0.378, 0.422)
Aortic Root Dimension (cm)	0.282 (0.277, 0.287)	0.428 (0.416, 0.440)

Multi-task regression performance on the internal YNHHS test set. Values in parentheses represent bootstrapped 95% confidence intervals. AV = aortic valve; ED = end-diastolic; ES = end-systolic; IVSd = intraventricular septum thickness at diastole; LA = left atrium; LV = left ventricle; LAIDs = left atrial internal diameter at systole; LVIDd = left ventricular internal diameter at diastole; LVIDs = left ventricular internal diameter at systole; LVOT = left ventricular outflow tract; LVPWd = left ventricular posterior wall thickness at diastole; MAE = mean absolute error; RA = right atrium; RV = right ventricle; RVIDd = right ventricular internal diameter at diastole; RV S' = right ventricular systolic excursion velocity; TV = tricuspid valve; YNHHS = Yale-New Haven Health System.

Extended Data Table 4 | External model performance on EchoNet-LVH

<i>Classification Tasks</i>			
		AUC	AP
Increased LV Size	<i>None</i>	0.892 (0.883, 0.894)	0.964 (0.960, 0.965)
	<i>Mild</i>	0.845 (0.838, 0.852)	0.513 (0.495, 0.533)
	<i>Moderate-Severe</i>	0.982 (0.979, 0.985)	0.744 (0.693, 0.764)
Increased LV Wall Thickness	<i>Any</i>	0.842 (0.835, 0.848)	0.939 (0.935, 0.942)
	<i>Moderate-Severe</i>	0.887 (0.879, 0.896)	0.992 (0.991, 0.992)
<i>Regression Tasks</i>			
		MAE	R²
IVSd (cm)		0.134 (0.132, 0.136)	0.450 (0.437, 0.462)
LVPWd (cm)		0.130 (0.128, 0.131)	0.349 (0.335, 0.364)
LVIDs (cm)		0.358 (0.353, 0.363)	0.704 (0.696, 0.712)
LVIDd (cm)		0.381 (0.376, 0.386)	0.634 (0.627, 0.642)

Multi-task performance on the external EchoNet-LVH dataset. Values in parentheses represent bootstrapped 95% confidence intervals. AP = average precision; AUC = area under the receiver operator characteristic curve; IVSd = intraventricular septum thickness at diastole; LV = left ventricle; LVIDd = left ventricular internal diameter at diastole; LVIDs = left ventricular internal diameter at systole; LVPWd = left ventricular posterior wall thickness at diastole; MAE = mean absolute error.

Extended Data Table 5 | External model performance on EchoNet-Dynamic

<i>Classification Tasks</i>			
		AUC	AP
LV Systolic Dysfunction	<i>None-Hyperdynamic</i>	0.886 (0.884, 0.897)	0.946 (0.942, 0.951)
	<i>Mild</i>	0.775 (0.766, 0.785)	0.347 (0.330, 0.366)
	<i>Moderate or greater</i>	0.941 (0.939, 0.949)	0.764 (0.760, 0.793)
<i>Regression Tasks</i>			
		MAE	R²
Ejection Fraction (%)		5.526 (5.447, 5.620)	0.613 (0.599, 0.626)
LV ED Volume (cm³)		25.701 (25.256, 26.169)	0.310 (0.295, 0.327)
LV ES Volume (cm³)		14.629 (14.281, 14.980)	0.493 (0.475, 0.510)

Multi-task performance on the external EchoNet-Dynamic dataset. Values in parentheses represent bootstrapped 95% confidence intervals. AP = average precision; AUC = area under the receiver operator characteristic curve; ED = end-diastolic; ES = end-systolic; LV = left ventricle; MAE = mean absolute error.

Extended Data Table 6 | Transfer learning performance on EchoNet-Dynamic

Model	Backbone	Initialization	MAE	R² (%)	Epochs to converge
3D	3DResNet18	Kinetics-400	5.63 (5.41, 5.87)	61.9 (57.8, 65.3)	5
2+1D	ConvNeXt-B	EchoCLIP	5.40 (5.17, 5.61)	64.8 (61.2, 68.2)	17
2+1D	ConvNeXt-T	Random	5.62 (5.37, 5.85)	59.8 (56.2, 63.5)	15
2+1D	ConvNeXt-T	ImageNet	5.37 (5.13, 5.61)	64.3 (60.4, 67.6)	9
2+1D	<i>ConvNeXt-T</i>	<i>PanEcho (ours)</i>	<i>4.69 (4.49, 4.86)</i>	<i>73.8 (70.9, 76.5)</i>	<i>2</i>

EF regression results when fine-tuning deep learning models with different architectures and weight initializations. Models were fine-tuned on the official training set and evaluated on the official test set of EchoNet-Dynamic. “Epochs to converge” refers to the number of epochs after which the best validation performance was reached before training was terminated via early stopping. Values in parentheses represent bootstrapped 95% confidence intervals. EF = ejection fraction; MAE = mean absolute error.

Extended Data Table 7 | Transfer learning performance on EchoNet-Pediatric

Model	Backbone	Initialization	MAE	R² (%)	Epochs to converge
3D	3DResNet18	Kinetics-400	4.82 (4.71, 4.94)	55.7 (52.8, 58.3)	13.7 ± 5.6
2+1D	ConvNeXt-B	EchoCLIP	5.16 (5.02, 5.30)	41.0 (37.4, 44.3)	12.7 ± 6.0
2+1D	ConvNeXt-T	Random	4.91 (4.78, 5.03)	51.1 (48.2, 54.1)	9.6 ± 3.1
2+1D	ConvNeXt-T	ImageNet	4.47 (4.36, 4.58)	61.2 (58.6, 63.6)	10.6 ± 3.4
2+1D	<i>ConvNeXt-T</i>	<i>PanEcho (ours)</i>	<i>3.87 (3.79, 3.95)</i>	<i>73.3 (71.1, 75.2)</i>	<i>5.5 ± 1.9</i>


Pediatric EF regression results when fine-tuning deep learning models with different architectures and weight initializations. Models were fine-tuned via 10-fold cross-validation given the official splits provided by EchoNet-Pediatric, then evaluated on the union of all held-out predictions. “Epochs to converge” refers to the mean ± standard deviation number of epochs (across the 10 cross-validation runs) after which the best validation performance was reached before training was terminated via early stopping. Values in parentheses represent bootstrapped 95% confidence intervals. EF = ejection fraction; MAE = mean absolute error.

Extended Data Table 8 | Description of multi-task labels


Task Name	Task Type	Classes
Aortic Root Dimension	Regression	
AV Peak Velocity	Regression	
AV Regurgitation	Multi-class classification	None-Trace, Mild, Moderate or greater
AV Stenosis	Multi-class classification	None, Mild-Moderate, Severe
AV Structure	Binary classification	Bicuspid
E e' Ratio	Regression	
Ejection Fraction	Regression	
Global Longitudinal Strain	Regression	
Increased LA Size	Multi-class classification	None, Mild, Moderate or greater
Increased LV Size	Multi-class classification	None, Mild, Moderate or greater
Increased LV Wall Thickness	Multi-label classification	Any, Moderate or greater
Increased RA Size	Binary classification	Any
Increased RV Size	Multi-class classification	None, Mild, Moderate or greater
IVSd	Regression	
LA Volume	Regression	
LAIDs	Regression	
LV Diastolic Dysfunction	Multi-class classification	None, Mild-Indeterminate, Moderate or greater
LV ED Volume	Regression	
LV ES Volume	Regression	
LV Stroke Volume	Regression	
LV Systolic Dysfunction	Multi-class classification	None-Hyperdynamic, Mild, Moderate or greater
LV Wall Motion Abnormalities	Binary classification	Any
LVIDd	Regression	
LVIDs	Regression	
LVOT Diameter	Regression	
LVOT PG	Binary classification	≥ 20 mmHg
LVPWd	Regression	
MV Regurgitation	Multi-class classification	None-Trace, Mild, Moderate or greater
MV Stenosis	Binary classification	Any
Pericardial Effusion	Binary classification	Any
RA Pressure	Binary classification	≥ 8 mmHg
RA Transverse Dimension	Regression	
RV S'	Regression	
RV Systolic Dysfunction	Binary classification	Any
RV Systolic Pressure	Regression	
RVIDd	Regression	
TAPSE	Regression	
TV Peak Gradient	Regression	
TV Regurgitation	Multi-class classification	None-Trace, Mild, Moderate or greater

Description of task types and subclasses. While multi-class classification involves discriminating mutually exclusive classes, multi-label classification involves discerning potentially overlapping classes. AV = aortic valve; ED = end-diastolic; ES = end-systolic; IVSd = intraventricular septum thickness at diastole; LA = left atrium; LV = left ventricle; LAIDs = left atrial internal diameter at systole; LVIDd = left ventricular internal diameter at diastole; LVIDs = left ventricular internal diameter at systole; LVOT = left ventricular outflow tract; LVPWd = left ventricular posterior wall thickness at diastole; MAE = mean absolute error; PG =

pressure gradient; RA = right atrium; RV = right ventricle; RVIDd = right ventricular internal diameter at diastole; RV S' = right ventricular systolic excursion velocity; TV = tricuspid valve.



Extended Data Fig. 1 | PanEcho performance under a simplified imaging protocol. Multi-task evaluation of *PanEcho* on the internal YNHHS test set leveraging at most one video from each of the following key views per study: PLAX, PSAX at the papillary muscle, A4C, A5C, and A2C. Error bars and values in parentheses represent bootstrapped 95% confidence intervals. **(a)** For classification tasks, AUC values are presented. The grey dashed line represents the performance of random guessing. * = moderate or greater; ^ = severe. **(b)** For regression tasks, normalized MAE (using the mean of ground truth measurements) is visually presented to account for varying scales and units. Raw MAE is also presented in text beside each task. AUC = area under the receiver operating characteristic curve; AV = aortic valve; A2C = apical 2-chamber; A4C = apical 4-chamber; A5C = apical 5-chamber; ED = end-diastolic; ES = end-systolic; IVSd = intraventricular septum thickness at diastole; LA = left atrium; LV = left ventricle; LAIDs = left atrial internal diameter at systole; LVIDd = left ventricular internal diameter at diastole; LVIDs = left ventricular internal diameter at systole; LVOT = left ventricular outflow tract; LVPWd = left ventricular posterior wall thickness at diastole; MAE = mean absolute error; PG = pressure gradient; PLAX = parasternal long axis; PSAX = parasternal short axis; RA = right atrium; RV = right ventricle; RVIDd = right ventricular internal diameter at diastole; RV S' = right ventricular systolic excursion velocity; TV = tricuspid valve; YNHHS = Yale-New Haven Health System.



Extended Data Fig. 2 | Task-specific view relevance. View relevance scores depicting the relative importance of each echocardiographic view (column) for a given task (row). For each task, a score of 1 indicates the most relevant view, and values in parentheses represent the rank importance of a view. **(a)** View relevance scores for all 39 tasks and all major 2D echocardiographic views. **(b)** View relevance scores for tasks known to be aided by Doppler imaging when including color Doppler videos (with potentially unknown view) as a “view”. AUC = area under the receiver operating characteristic curve; AV = aortic valve; A2C = apical 2-chamber; A3C = apical 3-chamber; A4C = apical 4-chamber; A5C = apical 5-chamber; ED = end-diastolic; ES = end-systolic; IVSd = intraventricular septum thickness at diastole; LA = left atrium; LV = left ventricle; LAIDs = left atrial internal diameter at systole; LVIDd = left ventricular internal diameter at diastole; LVIDs = left ventricular internal diameter at systole; LVOT = left ventricular outflow tract; LVPWd = left ventricular posterior wall thickness at diastole; MAE = mean absolute error; PLAX = parasternal long axis; PSAX = parasternal short axis; PG = pressure gradient; RA = right atrium; RV = right ventricle; RVIDd = right

ventricular internal diameter at diastole; RV S' = right ventricular systolic excursion velocity; TV = tricuspid valve.