

**Supplemental Table 1: Accuracy for Negative Samples. Real Time PCR Compared to Reference Method(s)**

Sample #	Specimen	NPM1 Type A Mutant RT-PCR Result	NPM1 Type A Mutant RT-PCR Value (% NCN, NPM1 Type A Mutant/ ABL) [Copies NPM1 MutA, Copies ABL1; Copies NPM1 MutA, Copies ABL1 ]	NGS Reference Result	ddPCR Reference Result (% NCN, NPM1 Type A Mutant/ ABL1 )	Pathol. Diagnosis
N15	BM	Negative	Undetectable [0,19000; 0,22300]	Negative	Negative: 0.005	AML
N17	BM	Negative	Undetectable [0,26400; 0,26800]	Negative	Negative: 0.003	AML
N20	PB	Negative	0.00 [0.009,41000 ; 0,44700]	Negative	Negative: 0.001	AML
N21	BM	Negative	Undetectable [0,66500; 0,66000]	Negative	Negative: 0.000	AML
N23	BM	Negative	Undetectable [0,31400 ; 0,33400]	Negative	Negative: 0.003	AML
N25	PB	Negative	Undetectable [0, 74500; 0,79800]	Negative	Negative: 0.000	AML
N32	BM	Negative	Undetectable [0, 9360; 0,8920]	Negative	Negative: 0.000	AML
N33	BM	Negative	Undetectable [0, 40600; 0,45000]	Negative	Negative: 0.005	AML
N34	PB	Negative	Undetectable [0, 54800; 0,50500]	Negative	Negative: 0.004	AML
N35	BM	Negative	Undetectable [0, 11700; 0,12000]	Negative	Negative: 0.000	AML

NOTES: Copies of NPM MUTA and ABL are shown for both duplicates in the RT-PCR assay.

%NCN = %Normalized Copy Number = (NPM1 Mutant Type A Copies/ ABL Copies)\*100


**Supplemental Table 2: Intra-Run Reproducibility**

Sample #	%, NCN:	Copy Number: [NPM1 MutA, ABL]	%, NCN:	Copy Number: [NPM1 MutA, ABL]	%, NCN	Copy Number: [NPM1 MutA, ABL]	%NCN Mean	%NCN SD	<u>Max Log Fold Change</u>
N16	ND	[0,54200 ]	ND	[ 0,53000 ]	ND	[ 0,39200 ]	0	0	/
N21	ND	[0,55400 ]	ND	[ 0,59000 ]	ND	[ 0,56500 ]	0	0	/
N22	ND	[0,53400 ]	ND	[ 0,55500 ]	ND	[ 0,43200 ]	0	0	/
P1, High+	496.51	[139000,27900 ]	566.75	[ 151000,26700 ]	745.57	[ 157000;21100 ]	602.94	128.41	0.18
P11, High+	659.38	[520000,78900 ]	776.07	[588000,75800]	871.57	[ 562000;64500 ]	769.01	106.27	0.12
P13, High+	691.32	[36800,5320 ]	617.79	[ 36900,5970]	695.95	[ 36800;5290 ]	668.35	43.85	0.05
P1, Low+	0.47	[26,5880]	0.33	[18,5740]	0.14	[8,5610]	0.31	0.17	0.53
P11, Low+	0.098	[38,38800]	0.13	[49,38000]	0.095	[34,36300]	0.11	0.02	0.14
P16, Low+	0.15	[100,67700]	0.23	[137,59000]	0.22	[156,70200]	0.20	0.04	0.19
Pos Cell Line, Low+	0.15	[225,150000]	0.16	[203,123000]	0.17	[206,122000]	0.16	0.01	0.05

NOTES: %NCN = %Normalized Copy Number = (NPM1 Mutant Type A Copies/ ABL Copies)\*100

<b>Supplemental Table 3: Analytic Sensitivity. Cell Lines.</b>					
<b>Sample; RNA Dilution (% Pos Control RNA)</b>	<b>NPM Type A Mutant, Number of Copies</b>	<b>ABL, Number of Copies</b>	<b>% NCN, NPM1 Type A Mutant/ABL</b>	<b>Mean, % NCN</b>	<b>SD</b>
<b>Cell Lines Dilution (% Pos Control RNA)</b>					
OCI-AML3_100%	449078.875	74524.29	602.594	603.331	1.042
OCI-AML3_100%	453850.406	75132.38	604.068		
OCI-AML3+MV411_10%	31229.891	48039.60	65.009	57.588	10.495
OCI-AML3+MV411_10%	26652.297	53127.45	50.167		
OCI-AML3+MV411_1%	3235.177	46716.25	6.925	6.738	0.264
OCI-AML3+MV411_1%	3376.028	51527.98	6.552		
OCI-AML3+MV411_10-1%	330.779	47363.50	0.698	0.671	0.039
OCI-AML3+MV411_10-1%	322.527	50119.39	0.644		
OCI-AML3+MV411_10-2%	32.126	56495.57	0.057	0.071	0.020
OCI-AML3+MV411_10-2%	43.579	51327.15	0.085		
OCI-AML3+MV411_10-3%	<b>3.231</b>	<b>54812.52</b>	<b>0.006</b>	<b>0.006</b>	<b>0.000</b>
OCI-AML3+MV411_10-3%	<b>3.091</b>	<b>55290.85</b>	<b>0.006</b>		
OCI-AML3+MV411_10-4%	1.817	68993.86	0.003	0.002	0.002
OCI-AML3+MV411_10-4%	0.366	73434.43	0.000		
OCI-AML3+MV411_10-5%	0.009	79998.87	0.000	0.000	0.000
OCI-AML3+MV411_10-5%	0.000	74978.66	0.000		
MV-4-11	0.001	70750.01	0.000	0.000	0.000
MV-4-11	0.000	69508.93	0.000		

# NPM1 MutA Standard Curves, n=19



# ABL Standard Curves, n=27

