

# Supplemental Appendix

## Table of Contents

1. AMA-1 hemi-nested PCR protocol for amplicon deep-sequencing.....	pages 3, 4
2. Figure 1: Bioinformatics workflow .....	5
3. Table 1: Detailed explanation of skip rule criteria.....	6
4. Table 2: Haplotype sequences and frequencies.....	7

## 1. AMA-1 hemi-nested PCR protocol for amplicon deep-sequencing

Primer	Sequence (5'→3')
Ama1OF	GCTGAAGTAGCTGGA ACTCAA
Ama1F+a*	acactctttcctacacgacgctcttccgatct- CCATCAGGGAAATGTCCAGT
Ama1R+a	gtgactggagttcagacgtgtgctcttccgatct- TTTCCTGCATGTCTTGAACA

\*a = adapter

### AMA1 Nested PCR Protocol

25µl Reaction

Reagent	Initial Concentration	Volume Added (µl)	Final Concentration
FastStart High Fidelity Buffer	10X with 18mM MgCl <sub>2</sub>	2.5	1X (1.8mM MgCl <sub>2</sub> )
dNTPs	10mM	0.5	200uM each dNTP
Roche High Fidelity Taq	5U / µl	0.25	2.5 U
Pf AMA1OF primer	10uM	0.5	167nM
Pf AMA1R primer	10uM	0.5	167nM
Template DNA	-	5	-
Water	-	15.75	-

Use same master mix protocol for Round 1 and Round 2, substituting “AMA1OF Primer” with “AMA1F-a”, “AMA1R primer” with “AMA1R-a primer” and “Template DNA” with “Round 1 PCR Product” in Round 2.

### Cycling Conditions

1<sup>st</sup> round PCR

Cycling Conditions R1	Temperature	Time
Step 1	95 C	2 min
Step 2	95 C	30 sec
Step 3	55 C	30 sec
Step 4	72 C	1 min
Step 5	Repeat step 2-4 total 20X	
Step 6	72 C	10 min
Step 7	4 C or 12 C	forever

2<sup>nd</sup> round PCR

<b>Cycling Conditions R2</b>	<b>Temperature</b>	<b>Time</b>
Step 1	95 C	2 min
Step 2	95 C	30 sec
Step 3	55 C	30 sec
Step 4	72 C	1 min
Step 5	Repeat step 2-4 total 10X, 20X, or 25X*	
Step 6	72 C	10 min
Step 7	4 C or 12 C	forever

\*10X for samples with parasite density between 10,000 and 1,000 parasites/ $\mu$ L

20X for samples with parasite density between 1,000 and 10 parasites/ $\mu$ L

25X for samples with parasite density between 10 and 0.1 parasites/ $\mu$ L


**INDEXING PCR**

50 $\mu$ L reaction

<b>Reagent</b>	<b>Volume Added (ul)</b>
NEBNext MasterMix	25 $\mu$ L
Sample (2 <sup>nd</sup> round PCR product)	15 $\mu$ L
Barcode	10 $\mu$ L

<b>Temperature</b>	<b>Time</b>	<b>Cycles</b>
98 C	30 sec	
98 C	10 sec	8-22 cycles
65 C	75 sec	
65 C	5 min	
10 C	forever	

## 2. Figure 1: Bioinformatics workflow



### 3. Table 1. Detailed explanation of skip rule criteria

**By clone:**

- 1) "Skips" = dates on which haplotypes recovered from a genotyped sample did not include a clone seen previously in that individual, or where qPCR was negative
- 2) 4 skips in a row were required to call the infection with that clone cleared
- 3) If there were  $\leq 3$  skips and the same haplotype is detected again, infection with that clone is presumed to have continued throughout the dates of the skips
- 4) If there were 4 skips and the same haplotype was detected again, it was considered a new infection with that clone

**By infection event:**

- 1) "Skips" = dates on which haplotypes recovered from a genotyped sample did not include any of the haplotypes included in an infection event, or where qPCR was negative
- 2) 4 skips in a row were required to call the infection event cleared
- 3) If there were  $\leq 3$  skips and any haplotype from the same infection event was detected again, that infection event is presumed to have continued throughout the dates of the skips
- 4) If there were 4 skips and haplotypes from a previous infection event were detected again, it was considered a new infection event

